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- 1950** Spectral photometer
- 1955** Flame photometer
- 1958** Piston-stroke pipette; Patent in 1960
- 1962** Microliter centrifuge & first automatic analyzer
- 1978** Hand-held dispenser Multipipette with Combitips
- 1985** Cell injector / Micromanipulator
- 1986** Sterile consumables for molecular biology
- 1997** PCR Cycler with gradient technology
- 1999** Electronic pipette
- 2003** Liquid Handling Workstations
- 2004** Realtime PCR
- 2007** Fermentors, Shakers, CO₂-Incubators & Freezers
- 2012** Parallel bioreactors and software solutions
- 2013** Integrated single-use bioreactors

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Alimentos & Bebidas con Eppendorf

Food research

The amount of arable land is not growing with the world population. Thus, applied research has to show ways to increase nutrient levels in crops, optimize crop rotation or make crops more drought tolerant. The seasonal crop samples needed for this research are invaluable.

Food production

Downstream processing of food has to be optimized to minimize nutrient loss. Nutrient levels have to be tested before and after new processing steps. Heating, cooling, packaging—everything can have an influence on the quality of food and beverages.

Food analysis

Food analysis and food quality control require reproducible workflows—in every detail. Only the best lab products give you the confidence in the data you need. Your results can make a big difference.

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	PP	PS	Relevant areas of application
Transparency	Medium	High	Transmission measurements
Temperature stability	ca. 120 °C	up to ca. 60 °C	Incubation, storage, autoclaving
Resistance to organic solvents	High	Low	Nucleic acid purification, protein analytic, compounds, assays
Mechanical strength	High	Low	Centrifugation, automation
Binding of biomolecules	Low	Low to high depending on the type	Applications using nucleic acids and proteins

Comparison of properties of the materials PP and PS

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
News

Plastics hamper DNA assays

Chemicals leaching from lab plastic throw off results.

Alla Katsnelson

Is it DNA or is it just chemicals leaching out of the tube?

dra_schwartz/iStockphoto

Biologists using standard plastic test tubes to gauge the concentration of DNA and proteins in their samples may be getting wildly incorrect readings because chemicals are leaching out of the containers.

An established way of assessing the concentrations and some key properties of DNA and proteins is to measure the

Published online 9 December 2008 | Nature | doi:10.1038/news.2008.1291

News

More biologists report plastic contamination

Chemicals from lab equipment are ruining experiments worldwide.

Daniel Cressey

Andrew Holt describes them as "horror stories".

Ever since Holt published a paper describing how standard pieces of plastic laboratory equipment leached experiment-ruining chemicals into his assays, researchers across the world have been in touch about their own experiences. And manufacturers are starting to take steps that could provide a workaround for the problem.


Reports of chemicals in plastic labware ruining experiments have come in from all over the world.

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Eppendorf Safe-Lock Tubes; Typical values for trace metal

The values in the table indicate typical values of trace metal concentrations which are obtained after incubating Eppendorf Safe-Lock Tubes with conc. nitric acid for 1 hour (see: Materials and Methods backside).

As the indicated values were determined in a one-time measurement, they cannot be guaranteed for every lot of Eppendorf Safe-Lock Tubes. Rather, they give the customer an idea to what extent trace elements can be eluted from the Eppendorf Tubes.

Eppendorf Safe-Lock Tubes	Trace metal release (ng/µl)								
	Al	Cd	Cr	Cu	Hg	Mn	Ni	Pb	Zn
Safe-Lock Tube 0.5 ml	0.003	<0.00002	<0.00005	<0.0001	<0.001	<0.00005	<0.00005	<0.00005	<0.001
Safe-Lock Tube 1.5 ml	0.002	<0.00002	<0.00005	<0.0001	<0.001	<0.00005	<0.00005	<0.00005	<0.001
Safe-Lock Tube 2.0 ml	0.002	<0.00002	<0.00005	<0.0001	<0.001	<0.00005	<0.00005	<0.00005	<0.001
Safe-Lock Tube Protein LoBind 0.5 ml	0.006	<0.00002	<0.00005	<0.0001	<0.001	<0.00005	<0.00005	<0.00005	<0.001
Safe-Lock Tube Protein LoBind 1.5 ml	0.004	<0.00002	<0.00005	<0.0001	<0.001	<0.00005	<0.00005	<0.00005	<0.001
Safe-Lock Tube									

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
Fig. 2: 384-well Eppendorf twin.tec PCR plate filled in a chess board-like pattern



Fig. 3: Amplification of the human β -globin gene across one fourth of a 384-well PCR plate [2].

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
Figure 4b: High quality genomic DNA
On a 1 % agarose gel 7.5 μ L out of 75 μ L mouse tail gDNA eluates were separated electrophoretically. The obtained DNA displayed a high molecular weight as indicated by the absence of low molecular weight smear. DNA size standard: Lambda *HindIII* (Fermentas).

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
Figure 4a: Reproducibility of DNA purification
Genomic DNA was isolated from 24 mouse tail samples, 10 mg each. The average yield was 17.3 μ g gDNA with a CV value of 2.09 %, giving evidence for a high level of consistency over 24 replicates. Each symbol represents the DNA yield of one sample; grey bars denote the average yield \pm 1 SD.

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> Alto rendimiento


Figure 5: DNA yields achieved with different methods
 Genomic DNA was isolated with different methods from 10 mg mouse tail samples. For each method 6 samples were processed in parallel. Other supplier's methods a and b were based on silica column technology. DNA yields obtained by the magnetic bead based manual and automated methods were considerably higher than the yields from silica column based methods.

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
Figure 6: No cross contamination detectable
 5 mg tissue samples were processed in a checkerboard pattern along with PBS, re-using the tips. The final eluates were diluted 1:10 and used in a real-time PCR to amplify a 212 bp fragment. Specific amplification was only detectable for positions with tissue lysates (red); from PBS filled positions no specific amplification was obtained (blue). Real-time PCR was performed on a Roche LightCycler system.

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Paula Pérez perez.p@eppendorf.es