

# Molecular solutions in the Agro-alimentary lab


MRAMA  
November 20<sup>th</sup> 2018

**Elisa Ruiz**  
Application Specialist  
Technical Service Department




# Agenda

- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- Molecular biology workflow with food samples
  - DNA Extraction and purification: Critical step
  - Maxwell® extraction instrument
  - Pathogen detection by real time PCR
- Application: Analysis of pathogens
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions



# Agenda

- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- Molecular biology workflow with food samples
  - DNA Extraction and purification: Critical step
  - Maxwell<sup>®</sup> extraction instrument
  - Pathogen detection by real time PCR
- Application: Analysis of pathogens
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions





# Promega: Mission and Vision

Mission: To provide reliable and personalized solutions that allows the progress in biomedical research, agro-alimentary industry, forensics and molecular diagnostic by offering quality products and best possible customer service.

Vision: To facilitate molecular biology techniques and tools for agroalimentary labs. No agroalimentary lab without PCR

# Promega Corporation

Biotech Multinational, global manufacturer

Sales~ \$400M, >4000 products portfolio,

~1400 employees

Presence in R&D markets


-  Promega Headquarters
-  Promega Manufacturing Facilities
-  Branch Office
-  Distributor



CELEBRATING OUR FIRST  
**40** YEARS  
**Promega**  
1978 - 2018



# We support...



# Standards of Quality



ISO 13485:2003



ISO 9001:2015

Sistema de Gestión  
ISO 9001:2015

www.tuv.com  
ID 9105047804



ISO 18385:2016




REACH Compliant



WEEE Compliant

# Agenda


- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- Molecular biology workflow with food samples
  - DNA Extraction and purification: Critical step
  - Maxwell<sup>®</sup> extraction instrument
  - Pathogen detection by real time PCR
- Application: Analysis of pathogens
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions





# Food Analyses

I. We offer a molecular biology based solution to control different parameters in food samples (DNA extraction kits + qPCR detection kits).



Microorganisms



Species identification, GMOs, allergens


II. Optimized and validated methods. Easily to accredit with ENAC.

III. Consulting and continuous technical support  (open to organize trainings, developing new possible applications)



# Agenda

- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- **Molecular biology workflow with food samples**
  - DNA Extraction and purification: A critical step
  - Maxwell<sup>®</sup> extraction instrument
  - Pathogen detection by real time PCR
- Application: Analysis of pathogens
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions




# Molecular biology workflow with food samples




- ✓ Fast method
- ✓ Minimal steps and less hands-on time
- ✓ Efficient with food inhibitors
- ✓ DNA ready to be used directly in a variety of downstream applications (real-time PCR, microarrays, or sequencing: NGS and Sanger)
- ✓ Adaptable Capabilities (16, 48, 96 sample format)

- ✓ Fast
- ✓ High specificity
- ✓ High sensibility
- ✓ Decrease of false negatives








# DNA Extraction and purification: A critical step

## Poor DNA quality:

- Unable to analyze sample
- Incorrect target assessment (pathogens, GMOs, allergens..)
- False +/-
- Misinterpretation of results
- Irreproducible results
  - Variation in event data
  - Inconsistency
  - Delays
  - Additional costs



*“Quality Data Begins with Quality Analytes”*

# Automated extraction system




*Maxwell*

- Extraction system based in **Paramagnetic Particles** which bind to the DNA/RNA
- Pre-filled cartridges with the reagents
- Extraction time: 20-60 minutes

Maxwell® HT option: Promega offers same chemistry adapted to high throughput formats for use with liquid handling platforms (Eg. Tecan, Hamilton, etc...)



Maxwell® RSC Instrument



Maxwell® HT in Liquid Handlers

# Pathogen detection by real time PCR

- **Food Pathogens**

- *Salmonella spp.*
- *Listeria monocytogenes*
- STEC
- *Campylobacter*
- *Vibrio*
- *Yersinia*
- Virus: HAV, HEV, Norovirus, ASFV, etc..

- Environmental pathogens


- *Legionella*

- Animal health pathogens

- Plant health pathogens

- Food industry


- Spoilage bacteria: beer, wine, juice, meats, etc..





# Agenda

- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- Molecular biology workflow with food samples
  - DNA Extraction and purification: A critical step
  - Maxwell<sup>®</sup> extraction instrument
  - Pathogen detection by real time PCR
- **Application: Analysis of pathogens**
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions



# Application: Analysis of *Listeria* + *Salmonella* in food samples

Efficiency and productivity

2 parameters

BPW  
(*Salmonella*)

HF  
(*Listeria*)

1 DNA extraction



1 PCR reaction

Bactoplex  
*Listeria* + *Salmonella*



2 results



Enrichment: 18h to  
26h



Extraction: 30 min



PCR: 1 h 25 min




# Application: Analysis of Legionella in water samples

UNE 100030:2017: Recognition of the PCR as alternative method to the culture


## Step 1. Sample prep

Sample Concentration:



Incubation-Lysis

Add the concentrate sample to the Maxwell® cartridge



## Step 2. Maxwell® Extraction System




Select the protocol and put the rack into the Maxwell®



~ 1h

## Step 3. Real-Time PCR Detection

Legiofast  
SPECIES




~1h 30min



# Agenda

- Promega Presentation
  - Mission and Vision
  - Promega Corp.
  - Our Standard of Quality
- Food analysis
- Molecular biology workflow with food samples
  - DNA Extraction and purification: A critical step
  - Maxwell<sup>®</sup> extraction instrument
  - Pathogen detection by real time PCR
- Application: Analysis of pathogens
  - Analysis of *Listeria* + *Salmonella* in food samples
  - Analysis of *Legionella* in water samples
- Application: Analysis of GMOs, food fraud and allergens
- Conclusions



# Analysis of GMOs, food fraud and allergens

## GMOs Analysis

- Food and seed imported to the EU from countries where there are genetically modified crops, they are suspect to have GMO traces.
- The EU requires this products must assessed and correctly labeled.
- When the content of GMO is >0.9%, normative requires label it as GMO product. If product has a content of 0.9% or lower is defined as non GMO

## Food Fraud

- Ambiguity in the practices of labelling
- Label fraud on products (adulterations, substitutions of ingredients)
- Concerns by consumers (price implications, diet restrictions, halal products)

## Allergens


- Increasing numbers of allergies and intolerances
- Labeling normative more stringent



# Application: Analysis of GMOs, food fraud and allergens

## Step 1. Sample prep

- Add CTAB, Rnase y Proteinase K.
- Incubate 30 min
- Centrifuge 10 min.
- Add 300µl of lysate and 300µl of Lysis Buffer into well #1 of the Maxwell® RSC cartridge



## Step 2. Maxwell® Extraction System



Select method and put the rack into the Maxwell®




## Paso 3. Detection kit (GMO, allergens, ID species).



Sample prep: 40 min

Extraction: 40 min


PCR: ~1h 30min



# Conclusions

- Molecular Biology increase capabilities of an agroalimentary lab
- Broad use of real time PCR for analyses detection.  
In next 5 years all agroalimentary labs will have a thermocycler in their facilities.
- Quality of DNA is crucial to get reliable and consistent results.

**Promega offers the technology necessary for the analysis of pathogens, GMO and allergens in food samples along with consulting and continuous technical support.**







**Promega**

**Thank you**  
**Questions, any comments?**

**[esp\\_techserv@promega.com](mailto:esp_techserv@promega.com)**

