



Comprehensive microbiology solutions







Illumina next-generation sequencing (NGS) offerings for microbiology






Data analysis & Bioinformatics solutions

BaseSpace™ Sequence Hub is a powerful cloud-based computing environment where NGS data can be easily stored, analyzed, and shared. Sequencing data from Illumina instruments are streamed in real-time over the internet to BaseSpace. The data is automatically converted to standard file formats (FASTQ) for analysis by applications (apps). Users launch apps in BaseSpace with just a few mouse clicks to analyze and visualize their data.




Whole-genome sequencing

-  **Bacterial Analysis Pipeline:** Identifies species, multilocus sequence type, plasmids, virulence, and antimicrobial resistance genes in bacteria.
-  **SPAdes Genome Assembler:** An open source tool for *de novo* sequencing. This app is designed to assemble small genomes.
-  **SRST2:** Reports sequence types from a MLST database and/or reference genes from a database of sequences for virulence genes, resistance genes, and plasmid replicons.
-  **Velvet *de novo* Assembly:** Use the Velvet assembler for *de novo* assembly of bacterial genomes.



16S rRNA sequencing

-  **16S Metagenomics:** Performs taxonomic classification of 16S rRNA targeted amplicon reads using an Illumina-curated version of the GreenGenes taxonomic database.
-  **Kraken Metagenomics:** Assigns taxonomic labels to short DNA sequences with high sensitivity and speed using exact alignments of k-mers and a novel classification algorithm.
-  **QIIME:** An open-source bioinformatics pipeline for performing microbiome analysis from raw DNA sequencing data.

Shotgun metagenomics

-  **GENIUS Metagenomics:** Know Now: CosmolID's curated genome database provides rapid, accurate, and actionable bacterial identification at the species, subspecies, and/or strain level.
-  **Kraken Metagenomics:** Assigns taxonomic labels to short DNA sequences with high sensitivity and speed using exact alignments of k-mers and a novel classification algorithm.
-  **MetaPhlAn (metagenomic phylogenetic analysis):** A computational tool for profiling the composition of microbial communities from metagenomic shotgun sequencing data.

Virology

-  **DeepChek® (HBV, HCV, HIV):** Performs deep sequencing analyses and reports on subtyping, genotyping, and inferred levels of resistance.
-  **Kraken Metagenomics:** Assigns taxonomic labels to short DNA sequences with high sensitivity and speed using exact alignments of k-mers and a novel classification algorithm.

Getting started is easy

When you partner with Illumina, you become part of a community with more than 10,000* publications in microbiology and virology. Our 'starter bundle' packages provide training, library prep, sequencing instrument, and reagent kits that support a wide range of sample volumes. We also offer a program that allows you to trade in your sequencer for an Illumina system.

Explore the exciting discoveries our NGS solutions are enabling. Visit www.illumina.com/MicroStories and www.illumina.com/MicroWebinars.



*Illumina maintains an up-to-date database of all published scientific articles that use Illumina technology.

Illumina • 1.800.809.4566 toll-free (US) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2018 Illumina, Inc. All rights reserved. Pub. No. 1270-2016-003 QB 6390

illumina®

Microbial Whole-Genome Sequencing with the iSeq™ 100 Sequencing System

Fast and efficient sequencing that provides uniform coverage and genome assembly for microbial species.

Highlights

- **Streamlined Workflow**
Access a comprehensive workflow from DNA to data
- **Optimized Library Prep**
Obtain robust, consistent results over a wide range of DNA input, even at low DNA input amounts (1 ng)
- **Comprehensive Coverage**
Produce sequencing data with uniform coverage for viruses, bacteria, and other microbes

Introduction

Next-generation sequencing (NGS) has been established as an important tool in microbiology research for analysis of small genomes (≤ 5 Mb), including bacteria, viruses, and other microbes. Microbial NGS, including whole-genome sequencing (WGS) and targeted resequencing, enables mapping and *de novo* assembly of novel organisms, completing genomes of known organisms, and comparing genomes across samples.

The development of Nextera™ chemistry shortened and simplified library preparation by consolidating DNA fragmentation and adapter tagging steps into a single reaction (termed tagmentation) and eliminating the need for library quantitation before pooling and sequencing.¹ The Nextera DNA Flex Library Preparation Kit represents the next step in the evolution of Illumina library prep. In addition to speed and efficiency gains in the workflow, Nextera DNA Flex offers exceptional flexibility for sample input type and amount and robust, consistent preparation of sequencing-ready libraries.



Figure 1: The iSeq 100 System—The iSeq 100 System harnesses the power of NGS in the most affordable, compact benchtop sequencing system in the Illumina portfolio.

The latest innovation in NGS is poised to advance microbiology genomics research. The compact iSeq 100 System (Figure 1) combines complementary metal-oxide–semiconductor (CMOS) technology with proven Illumina sequencing by synthesis (SBS) chemistry to deliver high-accuracy data with fast time to results. The iSeq 100 System is part of a streamlined NGS workflow for targeted and whole-genome microbial sequencing (Figure 2).




Figure 2: Microbial Sequencing Workflow—Microbial sequencing on the iSeq 100 System is part of a streamlined, comprehensive NGS workflow that includes Nextera DNA Flex library preparation, sequencing, and data analysis.

Simple, Integrated Workflow

Microbial sequencing on the iSeq 100 System is part of an integrated NGS workflow that includes library preparation with the Nextera™ DNA Flex Library Preparation Kit, proven Illumina sequencing, and push-button data analysis in BaseSpace™ Sequence Hub (Figure 2). The entire workflow proceeds from DNA to data in less than 24 hours.

Optimized Library Prep

A major advance in Illumina library prep chemistry and key feature of the Nextera DNA Flex Library Preparation Kit is On-Bead Tagmentation, which uses bead-linked transposomes (BLTs) to mediate simultaneous DNA fragmentation and the tagging of Illumina sequencing primers (Figure 3).




Figure 3: Nextera On-Bead Tagmentation Chemistry — (A) BLTs mediate tagmentation. (B) Reduced-cycle PCR amplifies sequencing ready DNA fragments and adds indexes and adapters. (C) Sequencing-ready fragments are washed and pooled.



To learn more about Nextera DNA Flex and On-Bead Tagmentation, read the Nextera DNA Flex Library Preparation Kit Data Sheet at www.illumina.com/nextera-dna-flex

Sequencing on the iSeq 100 System

After preparation, libraries are loaded into a prefilled reagent cartridge for sequencing on the iSeq 100 System. Starting a run on the iSeq 100 System is as easy as load and go with less than five minutes of setup. The iSeq 100 System integrates clonal amplification, sequencing, and data analysis into a single instrument. The intuitive user interface provides guidance through every step of the run setup and run initiation processes, allowing researchers to perform various sequencing applications with minimal user training and minimal set up time.



The iSeq 100 System harnesses proven Illumina SBS chemistry, used to generate more than 90% of the world's sequencing data.² Illumina SBS chemistry is used in all Illumina sequencing systems, enabling researchers to compare data across systems and scale their studies to higher throughput systems.

Easy, Flexible Data Analysis

The iSeq 100 System offers several data analysis options, including onboard and cloud-based data analysis. The Local Run Manager software, an onboard analysis software, features modular architecture to support current and future assays. Local Run Manager software supports planning sequencing runs, tracking libraries and runs with audit trails, and integration with onboard data analysis modules.

Alternatively, sequence data can be instantly transferred, analyzed, and stored securely in BaseSpace Sequence Hub, the Illumina genomics computing environment. BaseSpace Sequence Hub features a rich ecosystem of commercial and open-source apps for downstream data analysis, including the Integrative Genomics Viewer and BWA Aligner apps (Table 1).

Table 1: BaseSpace Apps for Microbial Sequencing Data Analysis

BaseSpace App	Description
 Integrative Genomics Viewer	The Integrative Genomics Viewer (IGV) app displays alignments and variants from multiple samples for performing complex variant analysis.
 BWA Aligner	The BWA Aligner app aligns samples (FASTQ files) to a reference genome using the Burrows-Wheeler Aligner maximal exact match (BWA-MEM) algorithm.

Comprehensive Coverage

To demonstrate the comparable performance of the iSeq 100 System to other sequencing systems in the Illumina portfolio in the genome assembly of microbial organisms, input genomic DNA from three different bacterial species with varying GC content (Table 2) were prepared with the Nextera DNA Flex Library Preparation Kit. Libraries were sequenced using paired-end 2 × 151 bp reads on the iSeq 100 System, MiniSeq™ System, and MiSeq™ System.

The iSeq 100 System delivers similar uniformity of coverage across different bacterial species, as compared to the MiniSeq and MiSeq Systems (Figure 4). These results support the exceptional performance of the iSeq 100 System for targeted and whole-genome microbial sequencing.

Table 2: GC Content of Sequenced Microbial Genomes

	<i>B. cereus</i>	<i>E. coli</i>	<i>R. sphaeroides</i>
Genome Size	~ 5.4 Mb	~ 4.6 Mb	~ 4.1 Mb
% GC Content	~ 35%	~ 51%	~ 69%




Figure 4: Consistent Uniformity of Coverage —The iSeq 100 System delivers similar uniformity of coverage across different three bacterial species, as compared to the MiniSeq and MiSeq Systems.

Summary

The iSeq 100 System is part of a fully supported solution for targeted and whole-genome microbial sequencing that includes simplified library preparation with the Nextera DNA Flex Library Preparation Kit, sequencing, and user-friendly data analysis. The iSeq 100 System delivers the same data quality as larger benchtop sequencers in a smaller footprint with faster run times, making it an ideal, cost-effective solution for small-scale microbiology NGS applications.

Ordering Information

Library Prep	Catalog No.
Nextera DNA Flex Library Prep Kit (24 samples)	20018704
Nextera DNA Flex Library Prep Kit (96 samples)	20018705
Nextera DNA CD Indexes (24 indexes, 24 samples)	20018707
Nextera DNA CD Indexes (96 indexes, 96 samples)	20018708
Sequencing System	Catalog No.
iSeq 100 System	20021532
Sequencing Reagents	Catalog No.
iSeq 100 i1 Reagents (300 cycles single kit)	20021533
iSeq 100 i1 Reagents (300 cycles quad kit)	20021534

Learn More

To learn more about the iSeq 100 System, visit www.illumina.com/systems/sequencing-platforms/iseq.html

To learn more about microbial whole-genome sequencing, visit www.illumina.com/microbiology.html

References

1. Illumina. Nextera XT DNA Library Preparation Kit Data Sheet. www.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet_nextera_xt_dna_sample_prep.pdf. Accessed September 2017.
2. Data calculations on file. Illumina, Inc., 2017.

16S metagenomics sequencing with the iSeq™ 100 System

Fast and efficient microbial sequencing on the most affordable Illumina sequencing system

Highlights

- Culture-free, NGS-based microbial analysis
Identify and compare bacterial populations from diverse microbiomes
- Cost-efficient microbial metagenomics
Study bacterial populations quickly and affordably
- Simple one-button data analysis
Analyze sequencing data easily with the 16S Metagenomics BaseSpace™ App

Introduction

Metagenomic studies are commonly performed by analyzing the prokaryotic 16S ribosomal RNA (16S rRNA) gene, which is approximately 1500 bp long and contains nine variable regions interspersed between conserved regions. Variable regions of 16S rRNA are frequently used for phylogenetic classification of genus or species in diverse microbial populations.¹⁻⁴

The choice and number of 16S rRNA regions to sequence are areas of debate, and the region of interest might vary depending on requirements such as experimental objectives, design, and sample type. This application note describes a comprehensive workflow that combines the Illumina demonstrated protocol for 16S metagenomics sequencing (Part # 15044223) with the iSeq 100 System (Figure 1) and secondary analysis using BaseSpace Sequence Hub.



Figure 1: The iSeq 100 System—The iSeq 100 System harnesses the power of NGS in the most affordable, compact benchtop sequencing system in the Illumina portfolio.

Simple, integrated workflow

16S metagenomics sequencing on the iSeq 100 System is part of an integrated next-generation sequencing (NGS) workflow that includes library preparation of the 16S V3 and V4 amplicon, proven high-quality Illumina sequencing, and push-button data analysis in BaseSpace Sequence Hub (Figure 2). The entire workflow proceeds from DNA to data in less than 30 hours.




Figure 2: 16S metagenomics sequencing workflow—16S metagenomics sequencing on the iSeq 100 System is part of a streamlined, comprehensive NGS workflow that includes library preparation, sequencing, and data analysis.

Table 1: Primer sequences for 16S metagenomics sequencing

Name	Sequence ^a
16S amplicon PCR forward primer ^b	5'- TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG - CCTACGGGNGGCWGCAG -3'
16S amplicon PCR reverse primer ^b	5'- GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG - GACTACHVGGGTATCTAATCC -3'

a. International Union of Pure and Applied Chemistry (IUPAC) nucleotide nomenclature: N = any base; W = A or T; H = A or C or T; V = A or C or G
 b. Primer sequence before the hyphen is Illumina overhang adapter sequence. Primer sequence after the hyphen corresponds to locus-specific sequence.

Library preparation

The 16S metagenomics sequencing workflow begins with PCR amplification of the V3 and V4 regions of the 16S rRNA gene using a bacterial primer pair selected from the scientific literature (Table 1).⁵ Illumina sequencing adapters and dual-index barcodes are then added to the generated amplicons using the Nextera™ XT DNA Index Kit. Libraries are normalized and pooled, and are ready for sequencing (Figure 3).




Figure 3: 16S V3 and V4 amplicon chemistry—Locus-specific primers with Illumina overhang adapters are used to amplify the V3 and V4 region of the 16S rRNA gene from genomic DNA. Sequencing adapters and dual-index barcodes are added, and libraries are normalized and pooled before sequencing.

Sequencing on the iSeq 100 System


After preparation, libraries are loaded into a prefilled reagent cartridge for sequencing on the iSeq 100 System. Starting a run on the iSeq 100 System is as easy as load and go with less than five minutes of setup. The iSeq 100 System integrates clonal amplification on a single instrument. The intuitive user interface provides guidance through every step of the run setup and run initiation processes, allowing researchers to perform various sequencing applications with minimal user training and set up time. The iSeq 100 System harnesses proven Illumina SBS chemistry, used to generate more than 90% of the world's sequencing data.⁶

Illumina SBS chemistry is used in all Illumina sequencing systems, enabling researchers to compare data across systems and scale their studies to higher throughput systems.

Easy, flexible data analysis

Sequence data can be instantly transferred, analyzed, and stored securely in BaseSpace Sequence Hub, the Illumina genomics computing environment. BaseSpace Sequence Hub features a rich ecosystem of commercial and open-source apps for downstream data analysis. The 16S Metagenomics App performs taxonomic classification of 16S rRNA targeted amplicon reads using a version of the GreenGenes taxonomic database curated by Illumina (Table 2).

Table 2: 16S Metagenomics BaseSpace App

BaseSpace App	Description
 16S Metagenomics App	Sample comparisons can be performed using the 16S Metagenomics App which enables analysis of 16S rRNA amplicon sequencing data and provides interactive visualization of taxonomic classification and relative abundance.

Experimental methods and results

To demonstrate the exceptional performance of the iSeq 100 System as part of a demonstrated protocol for 16S rRNA amplicon sequencing, data generated on the iSeq 100 System was compared against data generated on the MiSeq™ System for bacterial classification and relative abundance.

Methods

Samples and library preparation

Microbial genomic DNA samples were obtained from two sources for library prep, sequencing, and analysis. One source was the American Type Culture Collection (ATCC) 20 Strain Staggered Mix Genomic Material (ATCC MSA-1003). This mock microbial community comprises a staggered distribution of genomic DNA prepared from bacterial strains that were selected based on relevant attributes such as Gram stain, GC content values, and sporulation attributes. Real-world environmental samples were also obtained as part of a collaboration with academic researchers at the University of California, San Diego. Libraries were prepared following the 16S metagenomic sequencing library preparation workflow.⁷ Prepared libraries were normalized and pooled before sequencing.

Table 3: Comparison of multiplexing capacity by sequencing system

Sequencing system	Multiplexing capacity ^a			Run data quality	
	PF paired reads ^b	15K reads per sample	100K reads per sample	% PhiX	% Reads ≥ Q30 ^c
iSeq 100 System ^d 2 × 150 bp	4M	267	40	5	94.1
MiSeq System ^e 2 × 300 bp	25M	1667	250	10-25	74.9

- a. Based on recommended 15K-100K reads per sample for analysis with 16S Metagenomics BaseSpace App.
- b. Based on published instrument specifications.
- c. Average of Read 1 and Read 2 data.
- d. iSeq 100 System: v1 > 80% bases higher than Q30 at 2 × 150 bp.
- e. MiSeq System: v3 > 70% bases higher than Q30 at 2 × 300 bp.

Sequencing and data analysis

Prepared and pooled libraries were run at varying read lengths on the iSeq 100 and MiSeq Systems. Sequencing results were analyzed using the 16S Metagenomics App in BaseSpace Sequence Hub.

Comparison of multiplexing capacity on the iSeq 100 System

The multiplexing capacity of the iSeq 100 and MiSeq Systems shows a high sample multiplexing ability across all instruments based on the need for 15K-100K reads per sample for the 16S Metagenomics App in BaseSpace Sequence Hub. The iSeq 100 System is able to take advantage of the 384 indexes available with Nextera XT DNA Index kits (Table 3).

Comparable Q30 scores with the iSeq 100 System

Sequencing low-diversity libraries, such as those used for 16S rRNA sequencing, is challenging due to unbalanced base composition, causing a large percentage of the clusters to show the same base during each cycle. The high signals caused by the imbalance result in low Q-scores even though the base calling accuracy is not necessarily poor. Therefore, a 5% PhiX spike-in enables error rate calculations that allow verification of base calling accuracy over the course of the run, for all PhiX clusters, which can be extrapolated to the samples. Comparing the Q30 scores for the iSeq 100 System to the MiSeq System shows robust performance across all systems and run types, with the iSeq 100 System having higher Q30 scores while using less PhiX input (Table 3).



Quality score (Q-score): A metric in NGS that predicts or estimates the probability of an error in base calling. A Quality score (Q-score) serves as a compact way to communicate very small error probabilities. A high Q-score implies that a base call is more reliable and less likely to be incorrect.

Q30: A Q-score predicting that one in 1000 base calls will be incorrect. Q30 is widely considered a benchmark for high-quality data. A successful run will produce between 75-95% bases with Q30 scores or higher depending on the sequencing system, read length, and sequencing library quality.

Characterization of microbial composition across sequencing systems

To demonstrate the exceptional performance of the iSeq 100 System as part of a demonstrated workflow for 16S metagenomics sequencing, mock community and real-world samples were interrogated on the iSeq 100 and MiSeq Systems.

Characterization of ATCC Microbiome Standards

In order to compare performance across systems, the 20 Strain Staggered Mix Genomic Material (ATCC MSA-1003) was sequenced across multiple systems and run types. Analysis of the sequencing data with the 16S Metagenomics App on the iSeq 100 System identified all members of the bacterial community and showed comparable performance to the MiSeq System with fewer reads, shorter reads, and lower PhiX spike-in (Figure 4). The use of 2 × 300 bp read length for the MiSeq System showed equivalent identification of bacteria using the 16S Metagenomics App, demonstrating the robustness of the analysis tools with lower read lengths ideal for the iSeq 100 System (Figure 4).




Figure 4: Comparative analysis of microbial composition of ATCC samples across systems—Analysis of microbial composition of ATCC samples with the iSeq 100 System results in excellent genera coverage as compared to the MiSeq System.

Characterization of real-world samples

Microbial composition of real world soil samples (Figure 5 A and B) and fecal samples (Figure 5 C and D) were compared using the 16S metagenomics sequencing workflow with the iSeq and MiSeq Systems. The community profiles of all samples tested were highly concordant between the iSeq 100 and MiSeq Systems (Figure 5). These results further reinforce the use of shorter read lengths on the iSeq 100 System for 16S metagenomics applications with real-world samples.

To further interrogate the real-world samples, the 10 highest represented genera were compared between Fecal Sample 1 and Soil Sample 1 to demonstrate the difference in bacterial identity as well as the difference in distribution for the highest genera found in these two sample types (Figure 6). No overlap is seen in the top represented genera between the two samples and the fecal samples show the bacterial community is more heavily dominated by a smaller number of bacterial genera.




Figure 5: Comparative analysis of microbial composition of real world samples across sequencing systems—Analysis of fecal and soil samples for bacterial representation was highly concordant between the iSeq 100 and MiSeq Systems. Each axis is the fractional representation of each genera in each sample plotted against each other.

Summary

Using the 16S metagenomics workflow with the iSeq 100 System, microbiologists can achieve genus-level sensitivity for metagenomic surveys of bacterial populations. In this study, the Illumina workflow was used to study microbial populations in ATCC Microbiome Standards comprised of mock communities and in real-world samples. 16S metagenomic studies comprise one of many applications empowered by the iSeq 100 System. Illumina solutions support researchers during every step of the process, from DNA isolation through data analysis, enabling a range of applications for microbial genomics.

Illumina, Inc. • 1.800.809.4566 toll-free (US) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

© 2018 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html. Pub. No. 770-2018-009-B QB# 6012

10 highest represented genera




Figure 6: Analysis of bacterial populations in distinct microbiomes—Analysis of bacterial composition and distribution in fecal and soil samples showed no overlap in the top represented genera between the two samples. This confirms the samples are from two distinct microbiomes. **Note:** *Clostridium* was found in both fecal and soil samples tested but was not found to be one of the 10 highest represented genera in these samples.

Learn More

To learn more about the iSeq 100 System, visit www.illumina.com/systems/sequencing-platforms/iseq.html

To learn more about 16S metagenomics sequencing, visit www.illumina.com/areas-of-interest/microbiology/microbial-sequencing-methods/16s-rrna-sequencing.html

References

1. Woese CR, Fox GE. Phylogenetic structure of the prokaryotic domain: The primary kingdoms. *PNAS*. 1977;74(11):5088-5090.
2. Gray MW, Sankoff D, Cedergren JR. On the evolutionary descent of organisms and organelles: a global phylogeny based on a highly conserved structural core in small subunit ribosomal RNA. *Nucleic Acids Res*. 1984;12(14):5837-5852.
3. Woese CR, Kandler O, Wheelis ML. Towards a natural system of organisms: proposal for the domains Archaea, Bacteria, and Eucarya. *PNAS*. 1990;87(12):4576-4579.
4. Weisburg WG, Barns SM, Pelletier DA, Lane DJ. 16S ribosomal DNA amplification for phylogenetic study. *J Bacteriol*. 1991;173(2):697-703.
5. Klindworth A, Pruesse E, Schweer T, et al. Evaluation of general 16S ribosomal RNA gene PCR primers for classical and next-generation sequencing-based diversity studies. *Nucleic Acids Res*. 2013;41(1):e1.
6. Data calculations on file. Illumina, Inc., 2017.
7. Illumina. 16S Metagenomic Sequencing Library Preparation User Guide. support.illumina.com/documents/documentation/chemistry_documentation/16s/16s-metagenomic-library-prep-guide-15044223-b.pdf. Accessed April 2018.

Table 2: MiSeq System Configurations

Flow Cell	No. of Reads	Read Length	Output	No. of 16S Samples Per Run
600-cycle V3 standard flow cell	25 M	2 × 300 bp	15 Gb	Hundreds of 16S samples
500-cycle V2 standard flow cell	15 M	2 × 250 bp	8 Gb	
300-cycle V2 micro flow cell	4 M	2 × 150 bp	1.2 Gb	Tens of 16S samples
500-cycle V2 nano flow cell	1 M	2 × 250 bp	0.5 Gb	

Library Preparation

The Illumina 16S Metagenomic Sequencing Library Preparation Guide is an easy-to-follow protocol for preparing DNA libraries. It is optimized to target the V3 and V4 regions of the 16S rRNA gene, although it can be adapted to target other variable regions. The 16S Metagenomic Sequencing Library Preparation Guide leads users through each step of library preparation, from genomic DNA to sequencing-ready libraries. All necessary reagents are listed, including the required primer sequences that target the V3 and V4 regions of the 16S rRNA gene. These primers can also be modified to target different regions of the 16S gene, or altered for custom applications. The 27 samples from the reservoir were prepared using the 16S library preparation protocol and the Nextera® XT DNA Index Kit⁶ for cost-effective sample multiplexing.

Sequencing


The MiSeq System can deliver 2 × 300 bp reads and up to 50 million paired-end reads, generating up to 15 Gb of data. The flexible system enables microbiologists to scale studies from one to hundreds of samples. Micro and nano flow cell options and accompanying reagents are available to support lower-throughput experiments by optimizing sample volume and coverage needs (Table 2).

Samples from the reservoir were loaded onto a MiSeq reagent cartridge and then onto the instrument. Automated cluster generation and a 2 × 300 bp paired-end sequencing run were performed. The resulting sequence reads were equally distributed across the samples, demonstrating uniform coverage.

Data Analysis


Illumina has removed much of the complexity from sequencing data analysis. Following the Illumina workflow, researchers can analyze sequencing data generated on the MiSeq System either on the instrument or in BaseSpace. MiSeq Reporter software is able to analyze data on the sequencer or on a standalone computer. Alternatively, data can be transferred, analyzed, stored, and shared with collaborators in BaseSpace. BaseSpace can deliver analyzed sequences in as little as 12 hours following the 16S workflow, and BaseSpace applications (apps) provide access to a growing collection of analysis tools.

Figure 2: 16S Metagenomics App



The 16S Metagenomics App features an intuitive user interface to simplify analysis of environmental samples.

Figure 3: Relative Abundance by Taxonomic Level



The interactive sunburst chart in the BaseSpace 16S Metagenomics App shows the relative abundance of bacterial species. Users can select any classification to magnify a taxonomic level of interest. This chart shows the relative abundance of species in a single sample from the Eagle Creek reservoir.

The reservoir samples were analyzed using the BaseSpace 16S Metagenomics App (Figure 2). The app delivers all phylogenetic data—including coverage statistics and detected species—in intuitive, easy-to-analyze reports. Sequencing reads are classified against the Greengenes⁷ database, achieving up to species-level sensitivity.

Results

The 16S Metagenomics App delivers highly interactive visualizations for exploring taxonomic classifications. The sunburst classification chart provides a detailed view of the relative abundance of bacterial species within each taxonomic level. Researchers can select a category to magnify a particular level of interest and explore the diversity of any sample (Figure 3).

Exploring the Microbial Communities Within and Around Us

Next-generation sequencing with the MiSeq® System enables researchers to study the microbiota of humans, model organisms, and clinical environments.

Introduction

We are never really alone. Each of us is a host to flourishing populations of microorganisms arranged in communities referred to collectively as microbiota. Researchers are finding that regardless of whether they're located in the gut, skin, or airways, these communities possess great diversity that can change as we age, in response to certain diseases, changes in diet, or the ingestion of therapeutic drugs. While some pathogenic microorganisms can lead to disease or even death, many are essential to human health and well-being.

New microbial profiling approaches, such as 16S ribosomal RNA (rRNA) sequencing on the MiSeq System, have led to a greater understanding of our microbial communities and their interactions with us. Christopher Taylor, PhD, is part of the Louisiana State University Health Science Center (LSUHSC) Microbial Genomics Resource Group, an organization that supports microbial genomics with scientific expertise and research services. As an Associate Professor in the Department of Microbiology, Immunology, and Parasitology, Dr. Taylor uses rRNA and DNA sequencing approaches to investigate microbes of importance to human health.

iCommunity spoke with Dr. Taylor about his microbiome projects and how the MiSeq System has enabled his studies.

Q: How did you become involved in metagenomic studies?

Christopher Taylor (CT): I have a computer science and mathematics background and got involved with computational biology when I was in graduate school, where I was part of the US National Human Genome Research Institute's Encyclopedia of DNA Elements (ENCODE) project¹⁻³. I began my career as a faculty member at the University of New Orleans by focusing on applying high-throughput DNA and RNA sequencing in biological studies. One of my early projects was in collaboration with Dr. Erik Flemington, a virologist at Tulane University. We became interested in the RNA sequence reads from human cancer cell lines that did not map back to the human genome. In many labs at the time, the typical workflow was to map as many sequencing reads as possible back to the host genome, and then discard the remaining 15–20% of the reads. We wanted to look more closely at the nonmapping reads to see if we could find any viral, bacterial, or other recognizable sequences.⁴⁻⁶ This is still an active collaboration, and our most recent paper shows that there is a lot of microbial contamination in existing RNA sequencing data sets.⁷

Now that I'm at the LSUHSC School of Medicine, there's more of a health care focus to my work, and I've become immersed in research on microbial communities. My primary focus over the last 4 years has been using 16S rRNA sequencing to study the different microbial

communities that populate model organisms, humans, and the environments in which they live.

Q: What microbiomes are you studying?

CT: We have various ongoing studies looking at gut, vaginal, airway, and environmental microbiota. In a recent collaborative research study with Drs. Michael Ferris and Duna Penn at Children's Hospital of New Orleans, we used sequencing to look at the gut microbiota of infants in the neonatal intensive care unit, particularly premature infants suffering from necrotizing enterocolitis.⁸⁻⁹ Using 16S rRNA sequencing, we found that these infants have altered fecal microbiota characterized by a very low diversity in gut microbial communities, which might make them more susceptible to developing necrotizing enterocolitis.

Q: Have any of your studies looked at how diet impacts the gut microbiome?

CT: We've performed several studies where we've used sequencing to identify diet-associated variations in the gut microbiomes of mice.¹⁰⁻¹¹ In a recent collaborative study with Drs. Hans-Rudolf Berthoud, Annadora Bruce-Keller, Michael Salbaum, and David Welsh, we performed an antibiotic knockdown of the microbial gut community in a group of mice that had been on a standard mouse chow diet. By oral gavage, we then transplanted in the microbiota from mice that had been fed either a high-fat diet (HFD) or a standard mouse chow diet. Sequencing-based phylogenetic analysis using the MiSeq System confirmed the presence of a very distinctive difference in microbiota between the groups. The mice given HFD microbiota also showed



Christopher Taylor, PhD is an Associate Professor at the Louisiana State University Health Sciences Center.

For Research Use Only. Not for use in diagnostic procedures.

AGAATGATAACAGTAAACACACTTCTGTGTTAACCTTAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTTGCCAACGACGAAAGAAATGATAACAGTAACACACTTCTGTGTTAACCTT
TCAACGTACCGTAAACGAAACGTTATCAATTGAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTTGCCAACGACGAAAGAAATGATAACAGTAACACACTTCTGTGTTAACCTT
CGACGAAAGAAATGATAACAGTAAACACACTTCTGTGTTAACCTTAAGATTACTTGATCCACTGATTCAAACGTACCGTAAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTT
AACGTTACCATTAAAGAGCTACCGTTGCCAACGAAACGTTATCAATTGAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTTGCCAACGACGAAAGAAATGATAAC
AGATTGATAACAGTAAACACACTTCTGTGTTAACCTTAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAACG
GATTACTTGATCCACTGATTCAAACGTAAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAGCAACGACCG
CGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGATTACTTGATCCACTGATTCAAACGTACCGTAAACGAAACGTTATCAATTGAGACTAAATATAACGTTACCATTAAAGAGCTACCGTTGCCAACGAAAGAAATGATAACAGT

significant and selective differences in laboratory measures of mouse behavior, such as fear conditioning, maze solving, and anxiety display. This data reinforced the link between gut dysbiosis and obesity-associated changes in neurocognitive behavior in mice.

Sequencing-based phylogenetic analysis using the MiSeq System confirmed the presence of a very distinctive difference in microbiota between the groups.

Q: What other microbiome studies have you performed that enabled you to see the impact that microbiomes have on disease?

CT: We have conducted several genitourinary tract studies in collaboration with Drs. David Martin and Michael Ferris at LSUHSC. In one study, we looked at a disease called bacterial vaginosis and the role of *Gardnerella vaginalis*, which was originally thought to be the etiologic agent for that disease. A recent paper had shown that there are 2 fundamental *G. vaginalis* genotypes, a biofilm-forming pathogenic variant that was found in all patients with BV, and a commensal variant that does not form biofilms.¹²

In our study, we obtained samples from 53 women and their male sexual partners and used 16S rRNA sequencing to investigate the 2 *G. vaginalis* genotypes and their impact on disease.¹³ During that investigation, a graduate student in my lab, Murat Eren, developed a new method called oligotyping that allowed us to differentiate between various *G. vaginalis* strains based only on 16S sequencing data.

Dr. Martin had observed a phenomenon in the clinic where women diagnosed with bacterial vaginosis would finish a course of antibiotics, and be cured, but would then return to the clinic several weeks later with the same symptoms. Our idea was that the bacteria were harbored by the male sexual partner and transmitted back to the woman during sexual activity following the completion of antibiotic treatment. When we looked at the bacterial 16S oligotyping profiles of *G. vaginalis* from the women and men in the study, we found that we could determine which participants were sexual partners, without a priori knowledge and with a high degree of certainty. The ability to look this deeply into strain level variation using only 16S rRNA sequencing was unprecedented and has since been extended with a method called Minimum Entropy Decomposition by Murat and his colleagues and applied to several interesting microbial community studies.¹⁴

Q: What is the value of 16S sequencing?

CT: 16S sequencing enables us to identify bacterial species present in a community and to separate similar strains using techniques like oligotyping. We're currently working on an airway study in collaboration with Dr. David Welsh where we're using 16S sequencing data to look at relatively subtle differences between microbiota sampled from the oral cavity with microbiota from tissue samples taken further down the airway with different sampling brushes. Compared to gut microbiota, the bacterial burden of airway microbiota is a lot lower. The further down the airway you go, the fewer bacteria you will find. We want to determine how these microbiota transit

between these different environments and what this might tell us about the composition of microbiota living deeper in the airways.

Q: What types of environmental microbiome studies are you performing?

CT: We are planning a 16S sequencing study using the MiSeq System to analyze the metagenomic landscape of the Intensive Care Units (ICUs) in the new University Medical Center that is scheduled to open in mid-2015. We're very interested in looking at environmental samples from the new trauma ICU and medical ICU before patients are moved in and as these ICUs are put into service. Based on a trial run in another operating medical facility using qPCR to quantify 16S ribosomal RNA, we found the largest microbial counts came from the floor and from the lever on the hand sanitizer. In the new medical center, we're proposing to sample patient rooms and common areas where health care workers will be walking and using equipment, such as computer terminals, hand sanitizers, and scanners for tracking drug administration.

We also want to look at how antibiotic-resistant organisms end up entering and moving around in these new ICUs. One of our collaborators on the study is proposing a simultaneous culture-based investigation to look for particular antibiotic-resistant organisms. We're hoping to associate our 16S sequencing data with antibiotic-resistant and nonresistant versions of the cultured microorganisms and better understand how to control and prevent the spread of antibiotic-resistant microorganisms.

Q: What types of sequencing approaches did you use before obtaining a MiSeq System?

CT: In some of my initial research, we used genome tiling microarrays, but by the mid-2000s they had become obsolete with the advent of high-throughput sequencing. When I arrived at LSU, the School of Medicine had a pyrosequencing system that was being used to obtain long read lengths for 16S rRNA studies. This pyrosequencing system was older technology and was fraught with problems due to the large number of complex steps required for sample preparation. On many occasions, we struggled with this technology and had to perform multiple repeated sequencing runs attempting to resolve questionable results. It was also a long process, when you consider the time required for sample preparation, sequencing, and the weeks-long, intensive computation that was required to de-noise and prepare the data for analysis.

Q: How has the MiSeq System improved your workflow?

CT: We obtained a MiSeq System in August 2013, and the difference was like night and day. The cartridge-based MiSeq System eliminated many sources of variation where things could go wrong in sequencing preparation. Our sequencing throughput has been reduced from weeks to days. Due to the consistent and high-quality data generated by the MiSeq System, we were able to move forward with several of our projects. We didn't have to go back and repeat the same sequencing run over and over, trying to get usable data.

For our first MiSeq sequencing run, we used samples from the mouse gut repopulation study I mentioned earlier. When we ran these samples on the MiSeq System, the results were some of the best data we had ever seen. The beta diversity plot groupings obtained for the HFD mice versus the control diet mice were incredibly distinct and provided a better separation than we could have expected. We were

AATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAAAGTATCAATGAGACTAAATATTAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT
ACGTACCCTAAACGAAAGTATCAATGAGATTACTTGATCCACTGATTCAACGTACCCTAAACGAAAGTATCAATGAGACTAAATATTAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAG
ACGAAAGATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT
GTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT
AATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT
TACTTGATCCACTGATTCAACGTACCCTAAAGATTACTTGATCCACTGATTCAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT
TATCAATTGAGACTAAATATTAACGTACCCTAAAGAGTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGT

NGS Enables Efficient, High-Quality Microbiome Research

None of this work would have been possible without NGS systems, such as the MiSeq System, and bioinformatics programs that help researchers make sense of the data. “The irony of the situation is that collecting the samples is now the biggest problem I have in performing these studies,” Dr. Mills said. “Processing the samples and running them on a sequencer is easier than going to a facility and spending 8 hours obtaining 300 swabs, all documented so we know exactly where the sample was taken. In one of the papers we’re submitting soon, the winery had to collect about 2000 fermentation samples for us. That’s a large number of samples and they all have to be labeled correctly and stored somewhere. It’s a lot of work.”

For Dr. Mills’ research, the MiSeq System provides an optimum combination of speed and cost. “I like the MiSeq System because it’s a scaled down version of a high-throughput sequencer,” Dr. Mills stated. “It turns around samples quickly and I know that the quality is good. While I’d love to see a larger read, I’m happy with the reads that I get.”

His work has implications beyond identifying the microscopic creatures that help to make some of our favorite foods and wine. It might also help us understand how microbes move through our world.

“Our goal is to develop a microbial map of a whole facility,” Dr. Mills added. “One of the beauties of using NGS tools is that we can process a large number of samples. Our MiSeq runs are often 500 samples or more. In the old days, plating 500 samples would not have been feasible. We could never have performed the kind of high-throughput studies that we can today. With NGS data, we are opening up a completely new window, enabling us to view the microbial landscape of food facilities.”

References

1. Bokulich NA and Mills DA. Facility-specific “house” microbiome drives microbial landscapes of artisan cheesemaking plants. *Appl Environ Microbiol.* 2013;79(17): 5214–5123.
2. Bokulich NA, Ohta M, Richardson PM, Mills DA. Monitoring seasonal changes in winery-resident microbiota. *PLoS One.* 2013;8(6):e66437.
3. Bokulich NA, Thorngate JH, Richardson PM, Mills DA. Microbial biogeography of wine grapes is conditioned by cultivar, vintage, and climate. *Proc Natl Acad Sci USA.* 2014;111(1):139–148.

Learn more about the Illumina system mentioned in this article and the art of cheesemaking:

- MiSeq System, www.illumina.com/systems/miseq.html
- Illumina SciMon Video: Cheesemaking: Ancient Art or the Simplicity of Science, www.youtube.com/watch?v=MOQUzPcVp4&feature=youtu.be&list=PLKRu7cmBQLahAeoVCcb-pM1oVCKFkggFh



Follow the Illumina Agrigenomics Community

Sequencing

All 11 libraries were pooled together for cluster generation and sequencing. Libraries were loaded onto a reagent cartridge and clustered on the NextSeq 500 System. A paired-end, 2 × 150 bp sequencing run was performed using the NextSeq 500 High-Output Kit. A single sequencing run generated 400 million reads in 29 hours, corresponding to an average of 40 million reads per sample after quality filtering (Table 1). Base calls generated by the NextSeq 500 System were converted to FASTQ files for metagenomic analysis.

Data Analysis

Data analysis is often the most challenging part of a metagenomic sequencing study. Several programs and pipelines are available to perform metagenomic analysis. These programs are optimized for different study objectives, such as taxonomic profiling, assessing microbial composition, or identifying functional genes and pathways.

The analysis presented in this study used the publicly available MG-RAST program, an automated analysis platform that provides various tools for data visualization enabling users to assess species composition or functional abundance. FASTQ files generated by the NextSeq 500 System and sample metadata (Table 2) were uploaded to the MG-RAST server. Overlapping paired-end reads were joined and then processed in MG-RAST using the default trimming and filtering settings to allow comparison to other data sets in MG-RAST.

Table 1: Sequencing Data

MG-RAST ID	Metagenome Name	Base Count ^a	Sequence Count ^b
4554374.3	1305-530	5,061,868,863	33,395,202
4554375.3	1305-531	5,566,909,825	37,789,429
4554376.3	1305-532	8,409,674,885	57,196,837
4554377.3	1305-533	4,341,474,376	31,003,711
4554378.3	1307-521	5,184,242,473	32,673,383
4554379.3	1307-522	14,173,422,835	85,796,130
4554380.3	1307-523	8,369,469,023	52,993,984
4554381.3	1307-524	2,864,009,160	19,306,387
4554382.3	1310-522	8,688,763,557	55,503,302
4554383.3	1310-523	10,943,344,357	69,814,786
4554384.3	1310-524	7,776,089,825	46,838,779

a. All sequenced bases

b. Number of sequencing reads

Table 2: Environmental Metadata

MG-RAST ID	Sampling Depth	Sampling Date	Water Temperature (°C)	pH	Air Temperature (°C)	Location Coordinates ^a	Salinity (ppm)
4554374.3	Surface	23 May 2013	21.05	8.72	13.3	39.827, -86.303	0.21
4554375.3	3 meters	23 May 2013	21.00	8.73	13.3	39.827, -86.303	0.21
4554376.3	6 meters	23 May 2013	13.77	8.27	13.3	39.827, -86.303	0.22
4554377.3	Floor	23 May 2013	10.78	7.50	13.3	39.827, -86.303	0.20
4554378.3	Surface	25 July 2013	26.50	8.48	18.9	39.826, -86.304	0.21
4554379.3	3 meters	25 July 2013	26.02	8.40	18.9	39.826, -86.304	0.23
4554380.3	6 meters	25 July 2013	23.03	7.60	18.9	39.826, -86.304	0.22
4554381.3	Floor	25 July 2013	13.86	7.16	18.9	39.826, -86.304	0.23
4554382.3	3 meters	23 October 2013	15.09	7.80	4.4	39.826, -86.304	0.24
4554383.3	6 meters	23 October 2013	15.06	7.72	4.4	39.826, -86.304	0.24
4554384.3	Floor	23 October 2013	14.99	7.68	4.4	39.826, -86.304	0.24

a. Global Positioning System (GPS) coordinates

For Research Use Only. Not for use in diagnostic procedures.

AATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCCGTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTAA
 ACGTACCCGTAAACGAAACGATCAATTAAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCCGTTCTGTAAACCTTAAG
 ACGAAAGAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCCG
 GATACCAATTAAGAGCTACCCGTTCAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAAGCTACCCGTTCAACGAAAGAATGATAACA
 AATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAAGCTACCCGTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTAA
 TACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTACCAACGAA
 ATCAATTGAGACTAAATATTAACGTACCAATTAAGATTACTTGATCCACTGATTCAACGTACCCGTAACGAAACGATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCCGTTCAACGAAAGAATGATAACAGTAA




Figure 2: Comparative Genus Abundance Over Time—Comparisons of relative genus abundance (shown for the 15 most abundant genera) demonstrated a sharp decline in *Arthrobacter* populations on the reservoir floor in July 2013. These data were generated using MG-RAST.

Results

Analysis of species abundance with MG-RAST revealed drastic changes in microbial composition on the reservoir floor and an increase in the *Rhodococcus* genus (a soil bacteria), at the 3 meter depth in July 2013 (Figure 2). For the reservoir floor, the population declines might be correlated with an algaecide treatment that occurred on 31 May 2013. It is possible that the copper in the algaecide treatment caused the sharp decline in species from the *Arthrobacter* genus. This hypothesis is based on the observation that *Arthrobacter globiformis* often forms a symbiotic relationship with the *Anabaena* genus of nitrogen-fixing cyanobacteria and the aquatic ferns in the *Azolla* genus, both part of the algae family.⁶ After algaecide treatment, both *Arthrobacter globiformis* and *Azolla* experienced a decline; however, *Azolla* is not shown in Figure 2, because it is not 1 of the 15 most abundant genera in the samples. The hypothesis that the gram-positive *Anabaena* species might have played a role in the changes to community composition requires further analysis. Likewise for the 3 meter sampling depth, where analysis revealed an increase in the *Rhodococcus* genus, further analysis is required to assess root cause.

Conclusions

This study demonstrates how shotgun metagenomic sequencing with the NextSeq 500 System can reveal information about the microbial composition of a particular environment. The data presented show a correlation between the environmental metadata and the species composition in the Eagle Creek reservoir. This study provides new insights into the biological processes potentially associated with algal blooms, sampling depth, and seasonal differences in freshwater environments.

Learn More

To learn more about the NextSeq Series, visit www.illumina.com/nextseq.

For more information about the use of Illumina technology in microbial genomics, visit www.illumina.com/microbiology.

References

1. Meta-G-Nome- DNA Isolation Kit (www.epibio.com/applications/nucleic-acid-purification-extraction-kits/dna-purification-genomic/meta-g-nome-dna-isolation-kit) Accessed 21 April 2014.
2. Nextera XT DNA Library Prep Kit (www.illumina.com/products/nextera_xt_dna_sample_prep_kit.ilmn) Accessed 21 April 2014.
3. NextSeq 500 Desktop Sequencer (www.illumina.com/systems/nextseq-sequencer.html) Accessed 21 May 2015.
4. Meyer F, Paarmann D, D'Souza M, et al. The metagenomics RAST server - a public resource for the automatic phylogenetic and functional analysis of metagenomes. *BMC Bioinformatics*. 2008;9:386.
5. Citizens Energy Group (www.citizensenergygroup.com) Accessed 21 April 2014.
6. Martens T, Gram L, Grossart HP, et al. Bacteria of the *Roseobacter* clade show potential for secondary metabolite production. *Microb Ecol*. 2007;54:31-42.



AGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC
TCAACGTACCCTAACGAACGTATCATTAAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGAACGAAAAAAGAAATATAACAGTAACACACTTCTGTTAACCTT
CGACGAAAAGAAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTCAACGAGGAAAAGAAATGATA
AACGTACCATTAAAGAGCTACCGTGCAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGAGGAAAAGAAATGATA
AGAAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC
GATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGAGGAAAAGAAATGATAACAG
CGTATCAATTGAGACTAAATATAACGTACCATTAAAGATTACTTGATCCACTGATTCAACGTACCCTAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAACGAAAAAAGAAATATAACAG