

BWT AQU@SENSE MB

The Key to a Successful
Contamination Control Strategy
in Pharmaceutical Water

TIME TO TAKE NEW PATHS

The AQU@Sense MB allows you to meet the requirements for continuous germ load monitoring in pharmaceutical water fully automatically.

Flow Cytometry

Flow cytometry is not just a high-precision method for counting all living microorganisms; it can do much, much more.

- It detects individual microorganisms, microorganisms in larger conglomerates, or biofilm fragments and tells you the exact number of living cells.
- It creates a fingerprint of a microbiological population and can show any changes to it.

The Benefits

- Time saved for approval and troubleshooting
- Greater process reliability through trend analyses
- No "false positive" or "false negative" results

"Out of the Box" Measurements

- All the solutions required for up to 1000 measurements are placed in a hermetically sealed cartridge,
- which takes just a few minutes to replace.
- Up to 99% of the cartridge is recycled.

Online or Offline

Fixed installation or manual sampling

Specific

The staining method provides maximum specificity for living organisms. No disruption caused by particles.

Electronic Documentation

21 CFR Part 11-compliant electronic documentation

Fully Automatic


Up to 1000 measurements without any manual intervention

Proven Technology

The standard for drinking water and food for many years now

Primary Validation

In accordance with Ph.Eur. and USP requirements



FLOW CYTOMETRY


Quick, exact, reproducible and specific for intact cells

Explanatory
Video


Precise identification of intact microorganisms by evaluating scattered light (blue line) and fluorescence signals (green and red lines)


The "DOT plot" compares the two fluorescence signals. The area under the gate shows a signal for each intact organism. This DOT plot is a kind of fingerprint for a population.

→ The bottom line: flow cytometry counts every single living cell


You See Practically Nothing

But It Is There

Many living microorganisms do not grow on the plate, but they can be stained.

In Abundance

Most microorganisms occur as conglomerates. They would each form just a single colony.

CONTAMINATION CONTROL STRATEGY

Ph. Eur. 9.2 Chapter 5.1.6 and USP
Manufacturing practice (GMP) is a production...
Many countries have legisl...


Annex 1 sets it out plainly: every system operator needs to have control over all their processes at all times, and must be able to prove that they are functioning correctly. That means seamless, risk-based monitoring has to be ensured.

The Three Levels of a CCS

The basis for all the other steps on the path towards establishing a comprehensive CCS is accurate knowledge of the actual situation in the system.

In a data-driven approach, that means:

- Recording knowledge of the current situation for all process steps promptly
- Carrying out an assessment based on a risk analysis
- Taking suitable measures (CAPA) to correct or prevent errors


The Plate Count Is Not the Solution

The conventional plate count (HPC) has many disadvantages, particularly regarding the time required to obtain measurement results. The use of an alternative method is recommended in Ph. Eur., Annex 1, provided that the information can deliver a scientifically sound measurement for assessing microbiological quality. The AQU@Sense MB meets these requirements.

Alternatives Are Already Available Today

A variety of technologies for determining biological quality quickly and accurately are already available on the market.

The AQU@Sense MB uses flow cytometry technology, which is probably the most advanced technology currently available and can also be used for online bioburden monitoring. It does not just show the presence of microbiological contamination, but also provides an exact cell count. Furthermore, the measurement includes additional information about the population of the microorganisms.

A typical measurement result of flow cytometry.

TECHNOLOGY YOU CAN TRUST

Primary validation has shown that flow cytometry significantly outperforms the conventional plate count in some areas under comparable conditions.

Flow cytometry is an established microbial counting technique. It has been used in laboratories and to assess drinking water for many years now. The AQU@Sense MB is the first flow cytometer that enables online monitoring of water quality in pharmaceutical systems. The basis for the validation of this application is a primary validation in accordance with the applicable regulations (Ph.Eur. 9.2, chapter 5.1.6 and USP 41, chapter 1223). This was carried out by BWT in collaboration with bNovate, Eindhoven University of Technology (TU/e) and the University of Applied Sciences Northwestern Switzerland (FHNW).


Results

The AQU@Sense MB passed the primary validation in full!

Parameter	Results
Accuracy	✓
Precision	✓
Specificity	✓
Quantification limit	✓
Linearity	✓
Measurement range	✓
Robustness	✓
Ruggedness	✓

Your Next Steps...

Risk-benefit analysis
Validation for intended use

ONE INSTRUMENT, MULTIPLE APPLICATIONS

The AQU@Sense MB can be used online or offline as a standalone device.

IN PROCESS CONTROL

Requirements for the method

- Specific
- Reproducible
- Quick
- Cost-effective

Benefits of AQU@Sense MB


- Keep a global overview and control each step of production
- Monitoring of all process steps for biological contamination, from the incoming water to the end product
- Early detection (trending) of signs of contamination with microorganisms to take preventative measures
- Reduced costs for water samples and associated tests
- Time saved during investigations

Trend Monitoring

The AQU@Sense MB lets you monitor microbiological contamination changes promptly, long before critical values are reached.

Sanitation Monitoring

- Sanitization effect immediately demonstrable
- Faster approval following maintenance

QUALITY CONTROL UNIT

Requirements for the method

- Validated
- Specific
- Reproducible
- Correct


Benefits of AQU@Sense MB

- Immediate proof of quality
- Immediate intervention when problems occur
- Rapid analysis

Assess different results correctly

- Total water exposure (background signal)
- Detection and trend analysis of individual cells from conglomerates or biofilm fragments.


Trend Monitoring

Detecting changes in the system continuously and rapidly.

Sanitation Monitoring

Sanitization effect immediately demonstrable

SO TIME IS ON YOUR SIDE

Do you want to investigate yesterday's problem
or shape tomorrow's success?

The AQU@Sense MB gives you the data you need, quickly, reliably and conclusively. The AQU@Sense MB is not just a high-precision technology for counting all living microorganisms; it can do much, much more.




- It detects both individual organisms and organisms in larger conglomerates or biofilm fragments and tells you the exact number of living cells contained within them.
- It creates a fingerprint of a microbiological population in your system and can show any changes to it.

Automatic Digital Documentation

- No manual processes requiring time-consuming documentation.
- Analogue and digital interfaces to control systems
- Simple visualisation of results for further evaluations

Minimal Maintenance

- Menu-guided cartridge replacement in minutes following 1000 measurements
- No calibration

SERVICE AND SUPPORT

Minimal Maintenance

The AQUA@Sense MB is designed in such a way that scheduled maintenance is required just twice per year. It is performed by qualified BWT staff as part of the fully documented AQUA@Service.

Benefits

- No calibration required
- SST based on standard maintenance
- High availability as maintenance can be planned
- Transparency with maintenance and operating costs known in advance

Service Range

- Integration planning
- System integration
- Documentation
- Control system programming
- Automation adjustment
- Testing
- Project execution
- Change control
- Expertise/suggestions for improvement

TECHNICAL DATA

General


Measuring principle	Flow cytometry
Degree of protection	IP65
Ambient temperature	+5 °C to +35 °C
Humidity	10 to 90% RH
Dimensions (W × D × H)	350 × 240 × 373 mm
Weight	14 kg

Connections

Outputs	4 digital outputs and 2 × 4 to 20 mA analogue outputs, freely configurable
Inputs	4 digital inputs
Sampler input	¼" Swagelok tube fitting, male
Sampler output	¼" Swagelok tube fitting, male
Energy supply	DC 18 V, 1.4 A, ax. 20 W
Memory card	32 GB
Data export	USB and Ethernet

Measurements

Sample volume	90 µl
Measurement type	Continuous or manual
Measurement interval	30 min. to 6 hrs
Analysis duration	20 min.
Cartridge contents	Max. 1000 measurements
Sample flow	100 to 200 ml per min.
Max. sample temperature	40 °C
Sanitation capacity	Ozone up to 100 ppb Hot water up to 85 °C (without running measurement)

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The microbiological control of the quality of the water used in the pharmaceutical industry is analyzed in this article from the perspective of the application of rapid microbiological monitoring equipment (RMMs), considering the requirements for its implementation and its different possibilities of use, which They will open new avenues for monitoring and reassessing process flows, operations and work procedures.

KEYWORDS: GMP; RMM; OWBA; Microbiology; Sanitation

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TABLA 1.

(FUENTE: PH. EUR 5.1.6) TABLE 5.1.6.-1 - TASKS TO BE UNDERTAKEN DURING THE VALIDATION PROCESS

Activity	Normally carried out by	
	Supplier	User
Primary validation	+	-
URS (instrument, application)	-	+
Description of the technique	+	-
Risk benefit analysis	-	+
Design qualification (DQ)	-	+
Installation qualification (IQ)	-	+
Operational qualification (OQ)	-	+
Performance qualification (PQ):		
- verification of primary validation data given by the supplier;	-	+
- verification for the intended use (e.g. sterility testing, TAMC/TYMC, ...);	-	+
- method suitability test	-	+

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FIGURA 2. Resultados del citómetro de flujo antes y después de flamear los puntos de muestreo. Fuente: BWT AQUA AG
FIGURA 3. Resultados de un citómetro de flujo del experimento de parada y arranque. Fuente: BWT AQUA AG
Fuente: BWT AQUA AG