## OXOID BIOCHEMICAL IDENTIFICATION SYSTEM - MONO

*A rapid colorimetric test for the differentiation of Listeria monocytogenes from other Listeria species.*

### Intended use

The Oxoid Biochemical Identification System (O.B.I.S.) mono is a rapid colourimetric test for the determination of D-alanyl aminopeptidase (DALAase). It has been designed for the differentiation of presumptive *Listeria* species that have been isolated from selective media and plated on to a secondary medium for further biochemical testing.

### Principle of the test

The O.B.I.S. mono test offers a rapid screening method for differentiation of *Listeria monocytogenes* from other *Listeria* species. This reduces the need for full biochemical identification of all suspect colonies.

*Listeria* species, with the exception of *Listeria monocytogenes*, possess the enzyme D-alanyl aminopeptidase<sup>1,2,3</sup>.

Oxoid has developed a new system for aminopeptidase testing which uses a non-carcinogenic substrate. This is in response to health concerns associated with amino acid conjugates of  $\beta$ -naphthylamine<sup>4,5</sup> as these are potent carcinogens<sup>6</sup>.

D-alanyl-7-amido-4-methylcoumarin (DALA) is provided as a suspension. An acidic solution of dimethylaminocinnamaldehyde is used as a colour developer. If the substrate is hydrolysed by DALAase, free 7-amino-4-methylcoumarin (7AMC) combines with the developer to produce a purple Schiff's base<sup>7</sup>.

### Specimens

The test is designed for use from purity plates, not from primary isolation media, as the colonies on primary isolation media are too small to carry out an effective test.

Pick colonies which have typical *Listeria* morphology from selective *Listeria* isolation media such as Oxford Agar (CM0856), PALCAM Agar (CM0877) and chromogenic *Listeria* Media, and streak onto a purity plate.

O.B.I.S. mono tests can be performed from purity plates recommended by international standards, such as Tryptone Soya Agar (CM0131), Tryptone Soya Yeast Extract Agar (TSA-YE)<sup>8,9</sup> or a recognised chromogenic *Listeria* Medium<sup>8</sup>.

### Limitations of the test

O.B.I.S. mono is intended for the detection of DALAase in Gram-positive, catalase-positive, oxidase-negative short rod shaped bacteria, capable of growing on selective *Listeria* primary isolation media. It can be used as a screen to differentiate *Listeria monocytogenes* from other *Listeria*.

Occasionally aesculin positive *Bacillus* species may grow on *Listeria* isolation media. Bacilli are DALAase positive, but the colonies are different from *Listeria* and bacilli appear as large rods upon Gram staining.

Due to the small size of *Listeria* colonies on primary isolation media it is not possible to carry out the test directly from these plates. Use of multiple colonies from the primary isolation plate is not recommended as this may lead to a mixed culture and an incorrect result. International standards recommend sub-culturing presumptive *Listeria* species on to purity plates TSA (CM131), TSA-YE<sup>8,9</sup> or a recognised chromogenic *Listeria* Medium<sup>8</sup>.

Occasionally certain *Listeria* when grown on non-selective media form colonies that are difficult to remove with an inoculating loop. Ensure that there is sufficient colonial material on the inoculating loop before inoculating the reaction area. Failure to pick up sufficient material on the loop will result in a negative DALA test.

O.B.I.S. mono provides presumptive identification of *Listeria monocytogenes*, but does not replace full biochemical testing.

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2. Kämpfer, P. (1992) Differentiation of *Corynebacterium* spp., *Listeria* spp., and related organisms by fluorogenic substrates. *Journal of Clinical Microbiology* 30, 1067-1071.
3. Clark, A.G., McLauchlin, J. (1997) Simple color tests based on an alanine aminopeptidase reaction which differentiate *Listeria monocytogenes* from other *Listeria* species. *Journal of Clinical Microbiology* 35, 2155-2156.
4. Morgan, J. W. (1987) Evaluation of a rapid method for the determination of L-pyrrolidonyl- $\beta$ -naphthylamide hydrolysis. *Laboratory Medicine* 18, 682-683.
5. Bascomb, S. and Manafi, M. (1998) Use of enzyme tests in characterization and identification of aerobic and facultatively anaerobic Gram-positive cocci. *Clinical Microbiology Reviews* 11, 318-340.
6. Anonymous (2000) 2-naphthylamine. In Croner's *Hazardous substances: carcinogens guide*. pp. 2-291. Kingston-upon-Thames: Croner.CCH Group Limited.
7. Druggan, P., Roberts, P.B. and Swaine, D. (1999) A rapid chromogenic method for the differentiation of *Citrobacter* spp. and *Salmonella* spp. directly from enteric media. Abstract of the annual meeting of the American Society for Microbiology 1999, C444, p.71. Washington: ASM Press.
8. International Standards Organisation (1996) ISO 11290-1, *Microbiology of food and animal feeding stuffs – horizontal method for the detection and enumeration of Listeria monocytogenes - part 1: detection method*. London: British Standards Institute.
9. FDA (1998) Appendix 3: Trypticase Soy Agar with 0.6% Yeast Extract (TSAYE). In *Food and Drug Administration Bacteriological Analytical Manual, 8th edition, revision A* ed. Jackson, G. J. pp. 54-55. Gaithersburg: AOAC International.

## OXOID BIOCHEMICAL IDENTIFICATION SYSTEM - PYR

*The Oxoid Biochemical Identification System (O.B.I.S.) - PYR is a rapid colorimetric test for the determination of PYRase activity in streptococci and Citrobacter spp.*

### Principle of the test

PYRase activity distinguishes Group A *streptococci* and *enterococci* from other streptococcal groups including Group D *streptococci* (previously called faecal *streptococci*)<sup>1,2</sup>, Thus offering a rapid diagnostic alternative to time-consuming culture methods. In addition, the same enzymatic activity may be used to aid differentiation of *Salmonella* spp. from *Citrobacter* spp. and other *Enterobacteriaceae*<sup>3,4</sup>.

The PYRase test has traditionally been based on the use of a  $\beta$ -naphthylamide peptide. This is a potent carcinogen<sup>5</sup>. Oxoid has developed a new system using a non-carcinogenic substrate, in response to associated health concerns.

The O.B.I.S. - PYR test uses Test Cards impregnated with L-pyrroglutamic acid 7-amino-4-methyl-coumarin (7AMC) and dimethylamino-cinnamaldehyde for the detection of PYRase activity. The enzymatic hydrolysis of this substrate by enterococci, Group A *streptococci* and *Citrobacter* spp. produces a purple colour following the addition of the Developing Solution<sup>6</sup>.

### Specimens

For details of specimen collection and treatment a standard text book should be consulted<sup>7</sup>.

When identifying enterococci and Group A *streptococci*, fresh primary or secondary cultures grown overnight on non selective media such as blood agar give best results. Colonies tested must be Gram-positive cocci and catalase-negative. In case of insufficient growth, a subculture should be performed.

When identifying *Citrobacter* spp. from *Salmonella* spp. and *Enterobacteriaceae*, colonies from non-selective media such as XLD Medium, MLCB Agar, Desoxycholate Citrate Agar, Salmonella Shigella Agar, Brilliant Green Agar or Hektoen Enteric Agar may be tested. Colonies should be Gram-negative, oxidase-negative and urease-negative<sup>7</sup>.

## Expected results

Organism	Lancefield group	O.B.I.S. PYR
<i>Streptococcus pyogenes</i>	A	+
<i>Streptococcus agalactiae</i>	B	-
<i>Streptococcus</i> Group C	C	-
<i>Enterococcus faecalis</i>	D	+
<i>Enterococcus faecium</i>	D	+
<i>Streptococcus</i> Group F	F	-
<i>Streptococcus</i> Group G	G	-
<i>Citrobacter</i>	-	+
<i>Salmonellae</i>	-	-
<i>Escherichia coli</i>	-	-

### Limitations of the test

O.B.I.S.-PYR is intended for the detection of PYRase activity in Gram-positive, catalase negative cocci and *Citrobacter* spp. Less commonly encountered isolates of *lactococci* and *aerococci* may be PYRase positive. The confirmation of *enterococci* or Group A *streptococci* can be achieved by serological grouping with a suitable test, e.g. Oxoid Streptococcal Grouping Kit (DR 595), Oxoid Strep Plus (DR 575) or Oxoid Dryspot Streptococcal Grouping Kit (DR 400).

Some strains of *Enterobacter cloacae* are PYRase negative.

The reactions with O.B.I.S. - PYR are a marker for enzyme activity and atypical strains may occasionally occur.

A slight colour change may develop with negative reactions. This is restricted to the immediate site of inoculation. Refer to the positive and negative controls to aid interpretation.

Incubation beyond 20 seconds following addition of the Developing Solution may produce non-specific colour reactions. It is important therefore that the test is read as indicated.

### References

1. Ellener, P.D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of rapid colorimetric method for the presumptive identification of Group A *streptococci* and *enterococci*. *J. Clin. Microbiol.* **22**. 880-881.
2. Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D *streptococci*. *Appl. Microbiol.* **27**. 107-109
3. Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl- $\beta$ -naphthylamide hydrolysis test for differentiation of members of the families *Enterobacteriaceae* and *Vibrionaceae*. *J. Clin. Microbiol.* **31**. 1946-1948.
4. Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. *Letters in Applied Microbiology.* **28**. 175-178.
5. Bascomb, S. and Manafi, M. (1998). Use of enzyme tests in the characterization and identification of aerobic and facultatively anaerobic Gram-positive cocci. *Clinical Microbiological Reviews.* **11**. 318-340.
6. Data on file at Oxoid.
7. Faklam, R. R. and Washing, J. A. (1991). *Streptococci* and *aerococci* pp. 238-258, In Lunette, I. H. Balles, A., Hausler (Jr), W. J. and Shadomy, H. J. (eds) *Manual of Clinical Microbiology*, 5<sup>th</sup> edition. Amer. Soc. For Microb., Washington D.C.

## OBIS TEST FOR *Salmonella* DIFFERENTIATION FROM SELECTIVE MEDIA

OBIS Salmonella from Oxoid Limited is a simple, five-minute test used to differentiate *Salmonella* species from micro-organisms of similar colonial appearance grown on selective culture media.

This rapid chromogenic test quickly and reliably identifies false positive colonies such as *Citrobacter*, *Proteus* and *Morganella* species, providing a negative result within 5 minutes and considerably reducing the number of samples to require further lengthy testing for identification purposes. OBIS Salmonella thus eliminates potentially positive results more quickly, resulting in faster sample turnaround and more effective use of resources.

*Salmonella* presents a particular problem in food testing regimes since other non-pathogenic micro-organisms commonly found in food, such as *Citrobacter*, *Proteus* and *Morganella* species, may present with similar morphology on standard diagnostic culture media. Such presumptive positive colonies must then be further tested to eliminate or confirm the presence of *Salmonella*. Traditionally, this relies on expensive biochemical and serological tests and can take up to 24 hours to obtain a result.

In the meantime, until a negative result is obtained, the release of the product is delayed, production may be suspended and often the plant may be cleaned down as a precaution, resulting in considerable cost and concern for the manufacturer.

In the majority of cases it is sufficient to confirm that suspect colonies are NOT *Salmonella*, without identifying the species. It was in recognising this potential that Oxoid Limited developed OBIS Salmonella (OBIS: Oxoid Biochemical Identification System).




OBIS Salmonella incorporates two simple biochemical reactions on a single disposable test card. The first, PYR, detects PYRase activity in *Citrobacter* species and the second, NPA, demonstrates deaminase activity in *Proteus* and *Morganella* species (see illustration). Positive activity is demonstrated by a vivid colour reaction (PYR - purple, NPA - orange), indicating the presence of non-*Salmonella* species. Colonies that test negative for both reactions are confirmed as presumptive *Salmonella*.

By confirming that suspect colonies are not *Salmonella* in just 5 minutes, OBIS Salmonella allows batches to be released and production to continue up to 24 hours earlier than by traditional methods. Furthermore, fewer API tests (to identify *Salmonella*) are required, saving time and resources in the laboratory.

In trials involving 1000 samples grown on a variety of selective culture media, 470 produced presumptive positive *Salmonella* colonies. Of these OBIS Salmonella screened out 406 false positives, reducing the number of presumptive positive plates requiring further testing to 64. All positive *Salmonella* samples were detected (ref1).

### References:

1. Coleman, D.J., Barnwell, A.P. and Raggett, J. Birmingham PHLS West Midlands, Poster Presentation, Oxoid Limited.

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1. Hoosain, N. and Lastovica A.J. (2005). Evaluation of the Oxoid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. Sept 4-8, 2005, Gold Coast, Queensland, Australia. Griffith University.
2. Carlone, G.M., Valadez, M.J. and Pickett, M.J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. *Journal of Clinical Microbiology* 16 (6), 1157-1159.
3. Smith, C.M., Colborne, N.R., Stephens, P.J. and Druggan, P. (2006) A simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. Feb 8<sup>th</sup> 2006, Croke Park Conference Centre, Dublin, Ireland.

**APPROVED**

*By vzkmae at 4:57 pm, Apr 24, 2007*

ID0800M

# **OXOID BIOCHEMICAL IDENTIFICATION SYSTEM (O.B.I.S.) CAMPY**

## **Intended Use**

The Oxoid Biochemical Identification System (O.B.I.S.) campy is a simple rapid colorimetric test for the detection of L-alanyl aminopeptidase. It has been designed for the differentiation of *Campylobacter*, *Helicobacter* and *Arcobacter* species from other Gram-negative organisms and incorporates a Gram-lysis test which will demonstrate Gram status.

## **Principle of the Test**

The O.B.I.S. campy test will differentiate species of *Campylobacter*, *Helicobacter* and *Arcobacter* from all other Gram-negative organisms<sup>1</sup>. Unlike other Gram-negative organisms, *Campylobacteraceae* do not possess the enzyme L-alanyl aminopeptidase (L-ALA). The O.B.I.S. campy test incorporates a rapid test to detect this enzyme and a Gram-lysis reagent to rapidly determine the Gram status.

First, the Gram-lysis test (or a Gram stain) must be carried out. This test differentiates between Gram-positive and Gram-negative bacteria<sup>2</sup>. The test is carried out on a glass slide. Sodium Hydroxide (0.5M) is used to lyse the cell wall of Gram-negative organisms and release the DNA. The DNA forms a 'string' which can be seen when the loop is raised from the surface of the slide. This reaction does not occur with Gram-positive organisms.

Once the organism has been identified as Gram-negative, the L-ALA test can be carried out.

Each O.B.I.S. reaction card has been impregnated with the L-ALA substrate (L-alanyl-7-amino-4-methylcoumarin) in each of the six reaction zones. An acidic solution of dimethylaminocinnamaldehyde (DMAC) is used as a colour developer. If the substrate has been hydrolysed by the organism, the free 7-amino-4-methylcoumarin will combine with the developer to produce a purple Schiff's base.

## **Components of the O.B.I.S. campy Kit (ID0800M)**

Each kit contains the following reagents with enough material for 60 tests:

- ID0803M** O.B.I.S. campy Test Cards – one resealable pouch containing 10 cards, each with six reaction zones.
- ID0804M** O.B.I.S. campy Buffer – one white capped dropper bottle containing 7 ml of solution.
- ID0221M** O.B.I.S. DMAC Developer – one purple capped dropper bottle containing 7 ml of 0.5% w/v dimethylaminocinnamaldehyde in 1M hydrochloric acid.
- ID0802M** 0.5M Sodium Hydroxide – one flat white capped bottle containing 6 ml of 0.5M NaOH (sufficient for 600 Gram-lysis screening tests).
- ID0898** Paddle pastettes – pack of 60.  
Instruction leaflet.

## Materials required but not included

Clean glass slides.  
Sterile 10 µl plastic disposable inoculating loops.  
Positive and negative quality control organisms.  
Scissors

## Precautions

This product is for *in vitro* diagnostic use only.

Do not use the O.B.I.S. campy reagents beyond the stated expiry date.

The DMAC Developer contains a weak acid and will stain the skin. The 0.5M Sodium Hydroxide Solution is corrosive and may cause burns. Do not breathe fumes/vapour. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection.

Used Test Cards and inoculating loops should be disposed of as biohazardous waste. This should be incinerated or autoclaved at 121°C for at least 15 minutes.

*Campylobacter* are pathogens. A low infective dose can cause gastro-enteritis with potentially serious complications. Take appropriate precautions when handling potentially contaminated samples.

## Storage and Opening

The O.B.I.S. campy kit must be stored at 2°C to 8°C. Allow the pouches to reach room temperature before use to prevent condensation forming on the Test Cards. Open the pouches by cutting at the notch between the end seal and the clip-lock opening.

Remove the number of Test Cards required and reseal the pouch. Use within 60 minutes.

If fewer tests are required than the number on the Test Card, cut the card and return the unused portion to the pouch. Do not return used Test Cards to the pouch.

Precipitation will occur after long-term storage of the Sodium Hydroxide – this does not affect performance. A sterile loop should always be used. Discard if there is any sign of contamination.

## Quality Control Procedure

The following procedure should be performed each time the kit is used:

### Gram-lysis Test

1. **Positive Control** – Use a known Gram-negative (Gram-lysis positive) organism, such as *Campylobacter jejuni* ATCC®33291™. Follow the method given in the test procedure.
2. **Negative Control** – Use a known Gram-positive (Gram-lysis negative) organism such as *Bacillus cereus* ATCC®11778™. Follow the method given in the test procedure.

## L-ALA Test

1. **Positive Control** – use a known L-ALA positive organism, such as *Pseudomonas aeruginosa* ATCC®27853™. Follow the method given in the test procedure. Ensure that a purple colour forms around the colony material within 20 seconds.
2. **Negative Control** – Use a known L-ALA negative organism such as *Campylobacter jejuni* ATCC®33291™. Follow the method given in the test procedure. Ensure that no purple colour forms within 20 seconds.

## Specimens

The test is designed for use from purity plates such as Columbia Blood Agar. Primary isolation media should not be used as the colonies may be too small or too few to carry out an effective test. Pick colonies which have typical *Campylobacter* morphology from selective *Campylobacter* isolation media and streak onto Columbia Blood Agar. Incubate in a microaerobic atmosphere for 48 hours then conduct the O.B.I.S. test.

## Test Procedure and Interpretation of Results

### Gram-lysis Test

1. Take a 10 µl loopful of 0.5M Sodium Hydroxide solution (white flat-capped bottle) and place onto a clean, glass slide.
2. Using a sterile loop, take a small amount of material from a purity plate.
3. Mix the sample into the NaOH on the slide for up to one minute.
4. At intervals, carefully raise the loop from the mixture to check for the presence of a 'string' between the loop and the mixture.
5. Record the result:
  - A positive result is characterised by a viscous appearance to the mixture and the presence of the DNA string. This indicates that the organism is **Gram-negative**.
  - A negative result is characterised by the formation of a cell suspension with no string and indicates that the organism is **Gram-positive**.

This test may also be carried out using a 0.5M potassium hydroxide (KOH) solution. Alternatively, a traditional Gram stain may be used to ascertain the Gram status of the organism.

## L-ALA Test

Once the organism has been identified as Gram-negative, the L-ALA portion of the test can be carried out.

1. Remove one of the Test Cards from the pouch.
2. Using the paddle end of an unused plastic pastette, transfer the equivalent of 5x1 mm diameter colonies from a purity plate to the test area.
3. Spread the sample across the reaction zone (inside a circle) of a Test Card.
4. Dispense one drop of O.B.I.S. campy buffer (white capped dropper bottle) onto each of the inoculated reaction zones.

5. Wait for 30 seconds.
6. Dispense one drop of O.B.I.S. DMAC Developer (purple cap) onto each of the inoculated reaction zones.
7. The appearance of a purple colour around the original colony material within 20 seconds is a positive L-ALA reaction. A positive reaction indicates the organism is **not** a *Campylobacter*, *Helicobacter* or *Arcobacter* species.

If no colour develops around the original colony material within 20 seconds this is a negative reaction. A negative reaction indicates the organism is a presumptive *Campylobacter*, *Helicobacter* or *Arcobacter* species.


### Interpretation Chart

Typical reactions:

Organism	0.5M NaOH	L-ALA
<i>Campylobacteraceae</i>	+ (i.e. Gram neg)	-
Gram-positives	- (i.e. Gram pos)	-
Other Gram-negatives	+ (i.e. Gram neg)	+

### Venn Diagram

The Venn Diagram demonstrates how *Campylobacter* and *Campylobacter*-related organisms can be differentiated from other common bacteria using the Gram-lysis and L-ALA tests of O.B.I.S. campy.

## Limitations of the Test

O.B.I.S campy is intended for the detection of L-alanyl aminopeptidase in Gram-negative organisms. It can be used for the presumptive identification of *Campylobacter* and related organisms from pure culture.

The reaction with the O.B.I.S. campy L-ALA test is a marker for enzyme activity and atypical strains may occasionally occur. *Bacteroides ureolyticus* can give the same reactions as *Campylobacter* species. However, their colonial morphology and anaerobiosis aid differentiation.

Over time, the DMAC Developer may give a slightly pink reaction with the campy buffer even in negatives. However, this is easily differentiated from the clearly purple reaction seen around the colony material in a positive.

Do not use nichrome wire loops to inoculate cards as this material can interfere with the test.

## Performance

In an internal study, 46 *Campylobacter* species and 6 *Arcobacter* species were tested. All yielded a Gram-negative, L-ALA negative reaction. 252 other species (non-*Campylobacter*, *Arcobacter* or *Helicobacter*) were tested. Only one organism gave a result similar to *Campylobacter*, and was due to an atypical Gram-lysis reaction. This resulted in a sensitivity and specificity of 100% and 99.6% respectively<sup>3</sup>.

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## References

1. Hoosain, N. and Lastovica A. J. (2005). Evaluation of the Oxid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. Sept 4-8, 2005, Gold Coast, Queensland, Australia. Griffith University.
2. Carlone, G. M., Valadez, M. J. and Pickett, M. J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. *Journal of Clinical Microbiology* 16 (6), 1157-1159.
3. Smith, C. M., Colborne, N. R., Stephens, P .J. and Druggan, P. (2006) a simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. Feb 8th 2006, Croke Park Conference Centre, Dublin, Ireland.

OXOID, Wade Road, Basingstoke, Hampshire, RG24 8PW, England.  
Tel: +44 1256 841144. Fax: +44 1256 463388. Telex: 858793 OXOID G.