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1. ISO 6579:2002: Microbiology of food and animal feeding stuffs. Horizontal method for detection of *Salmonella* spp.
2. ISO 16140:2003: Microbiology of food and animal feeding stuffs. Protocol for the validation of alternative methods
3. Statistical Techniques in Bioassay; Z. Govindarajulu; 2001
4. A novel chromogenic ester agar medium for detection of Salmonellae. Cooke VM et al. *Appl Environ Microbiol.* 1999 Feb ; 65(2): 807-12

RAPID TWO DAY ISOLATION OF SALMONELLA USING SINGLE SELECTIVE ENRICHMENT AND BRILLIANCE SALMONELLA AGAR

Jeff J. Kozickowska, Dorn L. Clarka, Roy P. Radcliff¹, Christine M. Hilbert², and Alastair Thomas³

¹Marshfield Food Safety 1000 N. Oak Ave, Marshfield WI 54449 USA ²Thermo Fisher Scientific 12076 Santa Fe Drive, Lenexa KS 66215 USA and ³Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

ABSTRACT

Salmonella is a gram-negative, rod-shaped, motile bacterium with a widespread occurrence in animals, especially in poultry and swine. The most common food sources contaminated with *Salmonella* are poultry, beef, pork, eggs, and fresh produce¹. Detection is critical as *Salmonella* infections cause an estimated 1.4 million human illnesses and 400 deaths annually in the United States². It has been reported an infectious dose of as low as 17 cells may cause infection³. The data presented from this study describes the evaluation of Oxoid ONE Broth-Salmonella and Brilliance™ Salmonella Agar chromogenic medium as a detection method for *Salmonella* in ground beef, ground chicken, lettuce, shrimp, and shell eggs. The selected foods were inoculated with *Salmonella* serovars at levels ranging from low (1 colony forming unit (CFU)/25g) to high (1.1 CFU/g) and analyzed for growth.

The method was also evaluated for inclusivity using multiple *Salmonella* serovars (n=102) from various sources. Exclusivity testing was performed using closely related bacterial species (n=30) to determine the specificity of the described method. Overall, when compared to the standard USDA and FDA reference methods there was no significant difference in sensitivity between Brilliance Salmonella Agar and the standard reference method. However, when Oxoid ONE Broth-Salmonella was used as the enrichment step prior to plating on Brilliance Salmonella Agar, the time to detection was reduced from 3-4 days to 2 days. For the inclusivity study, 98 *Salmonella* serovars grew on Brilliance Salmonella Agar. For the exclusivity study, 30 non-*Salmonella* spp. were tested to determine the specificity of the described method and only one, *Cronobacter sakazakii*, demonstrated typical growth on the chromogenic medium, but only when not enriched in ONE Broth-Salmonella prior to plating.

INTRODUCTION

Brilliance Salmonella Agar is a novel chromogenic medium for the detection and identification of *Salmonella* spp. in food. When used in conjunction with Oxoid ONE Broth-Salmonella it reduces *Salmonella* detection time from 3-4 days to 2 days. This study evaluated and compared the performance of the Salmonella Rapid Culture method to the reference U.S. Department of Agriculture-Food Safety Inspection Service (USDA/FSIS) and U.S. Food and Drug Administration's Bacteriological Analytical Manual (FDA/BAM) methods for the detection of *Salmonella* spp. in raw ground beef, ground chicken, lettuce, shrimp, and shell eggs. Results demonstrate the Salmonella Rapid Culture method to be superior to the reference methods in incubation time and ease of identification.

METHOD COMPARISON

Ground beef, ground chicken, lettuce, shrimp, and shell eggs were inoculated with approximately 1 CFU/25g and allowed to equilibrate per AOAC instructions. Side-by-side samples were set up and run using Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar chromogenic medium and the appropriate FDA/USDA reference for each matrix. Three different lots of Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar were evaluated. Stability and lot-to-lot variation were evaluated simultaneously. Organisms were initially grown in BHI (brain-heart infusion broth) and incubated overnight. Inoculation levels of *Salmonella* serovars Typhimurium and Enteritidis were approximately 1-10 CFU into 225ml of ONE Broth. A non-*Salmonella* organism (*E. coli* O157:H7) was inoculated into 225ml of Oxoid ONE Broth-Salmonella at a concentration achieved with overnight growth in BHI.

INCLUSIVITY: 102 *Salmonella* (table 3) strains were cultured in Oxoid ONE Broth-Salmonella and plated on Brilliance Salmonella Agar according to manufacturer's instructions.

EXCLUSIVITY: All strains were cultured overnight in BHI. Strains were then plated on Brilliance Salmonella Agar and incubated per manufacturer's instructions.

LOT TO LOT AND STABILITY: Ruggedness : Organisms were initially grown in BHI and incubated overnight. Inoculation levels of *Salmonella* serovars Typhimurium and Enteritidis were approximately 1-10 CFU into 225ml of Oxoid ONE Broth-Salmonella. A non-*Salmonella* organism (*E. coli* O157:H7) was inoculated into 225ml of Oxoid ONE Broth-Salmonella at a concentration achieved with overnight growth in BHI. Variations were introduced as follows:

Ruggedness parameters tested:

- Incubation temperature of Oxoid ONE Broth-Salmonella: 40, 42, 44°C.
- Incubation temperature of Brilliance Salmonella Agar plates: 33, 35, 37, 39, 41°C
- Incubation time of Brilliance Salmonella Agar plates: 21, 24, 27 hrs.

RESULTS

The results obtained for ground beef, ground chicken, lettuce, shrimp, and shell eggs demonstrate that the Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar are equivalent to the reference method in the foods tested (table 1). When specificity was evaluated 98/102 *Salmonella* serovars were identified using this method (table 3) and 30/30 of the non-*Salmonella* showed no growth or atypical growth. The three different lots of Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar evaluated gave equivalent results (table 2). Variations in incubation time and temperature of the Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar did not effect the performance of the method.


Table 1: Method Comparison

INOCULATING ORGANISM	LEVEL	MPN/25g	NO. TEST PORTIONS	REFERENCE METHOD		TEST KIT	
				POSITIVE	PRESUMP. POSITIVE	POSITIVE	PRESUMP. POSITIVE
<i>Salmonella Muenchen</i>	Low	0.9	20	10	8	8	8
	Control	0	5	0	0	0	0
<i>Salmonella Montevideo</i>	Low	2.3	20	17	15	15	15
	Control	0	5	0	0	0	0
<i>Salmonella Heidelberg</i>	Low	5.75	20	15	14	14	14
	Control	0	5	0	0	0	0
<i>Salmonella Choleraesuis</i>	Low	<0.75	20	9	10	10	10
	Control	0	5	0	0	0	0
<i>Salmonella Typhimurium</i>	Low	0.9	20	8	11	11	11
	Control	0	5	0	0	0	0

Table 2: Ruggedness and Lot to Lot Stability

	Lot 1:08A14-1	Lot 2: 08B06-1	Lot 3 08B13-1
	EXPIRE: Broth:2-11-08 Plates:1-28-08	EXPIRE: Broth:3-5-08 Plates:2-27-08	EXPIRE: Broth:2-13-08 Plates:2-13-08
<i>Salmonella Typhimurium</i>	5/5 positive	5/5 positive	5/5 positive
<i>Salmonella Enteritidis</i>	5/5 positive	5/5 positive	5/5 positive
<i>E. coli</i> O157:H7	5/5 negative	5/5 negative	5/5 negative

References: 1. Foley, S.L., A.M. Lynne, and R. Nayak. (2008). *Salmonella* challenges: prevalence in swine and poultry and potential pathogenicity of such isolates. *J. Anim. Sci.* 86 (14 Suppl):E149-62. Epub (2007) Oct. 2. Voetsch A.C., T.J. Van Gilder, F.J. Angulo, et al. (2004). FoodNet estimate of the burden of illness caused by nontyphoidal *Salmonella* infections in the United States. *Clin. Infect. Dis.*; **38** (Suppl. 3):S127-34. 3. Blaser, M.J., and L.S. Newman. (1982). A Review of Human Salmonellosis: I. Infective Dose *Reviews of Infectious Diseases*, Vol. 4, No. 6 pp. 1096-1106.

Table 3: Inclusivity

GENUS	SEROVAR	RESULT	GENUS	SEROVAR	RESULT
<i>Salmonella</i>	4,27:undetermined	Typical Growth	<i>Salmonella</i>	jena	Typical Growth
<i>Salmonella</i>	4,5,12: i:-	Typical Growth	<i>Salmonella</i>	Kahla	Typical Growth
<i>Salmonella</i>	4,5,12:b:- variant L(+)	Typical Growth	<i>Salmonella</i>	Kentucky	Typical Growth
<i>Salmonella</i>	4,5,12:i:-	Typical Growth	<i>Salmonella</i>	Kitenge	Typical Growth
<i>Salmonella</i>	6,7:k:-	Typical Growth	<i>Salmonella</i>	Kunzendorf	Typical Growth
<i>Salmonella</i>	abaetetuba	Typical Growth	<i>Salmonella</i>	Litchfield	Typical Growth
<i>Salmonella</i>	Abony	Typical Growth	<i>Salmonella</i>	London	Typical Growth
<i>Salmonella</i>	abortusovis	No growth*	<i>Salmonella</i>	Maarsen	Typical Growth
<i>Salmonella</i>	adelaide	Typical Growth	<i>Salmonella</i>	Maartensdijk	Typical Growth
<i>Salmonella</i>	Agona	Typical Growth	<i>Salmonella</i>	Maastricht	Typical Growth
<i>Salmonella</i>	Albany	Typical Growth	<i>Salmonella</i>	Mbandaka	Typical Growth
<i>Salmonella</i>	Anatum	Typical Growth	<i>Salmonella</i>	Meleagridis	Typical Growth
<i>Salmonella</i>	Arizonae	Typical Growth	<i>Salmonella</i>	Menden	Typical Growth
<i>Salmonella</i>	bareilly	Typical Growth	<i>Salmonella</i>	michigan	Typical Growth
<i>Salmonella</i>	berta	Typical Growth	<i>Salmonella</i>	Minnesota	Typical Growth
<i>Salmonella</i>	blockley	Typical Growth	<i>Salmonella</i>	Montevideo	Typical Growth
<i>Salmonella</i>	bongori	Typical Growth	<i>Salmonella</i>	Muenchen	Typical Growth
<i>Salmonella</i>	Braenderup	Typical Growth	<i>Salmonella</i>	muenster	Typical Growth
<i>Salmonella</i>	Brandenburg	Typical Growth	<i>Salmonella</i>	Newington	No growth*
<i>Salmonella</i>	Bredeney	Typical Growth	<i>Salmonella</i>	Newport	Typical Growth
<i>Salmonella</i>	Breukelen	Typical Growth	<i>Salmonella</i>	Ngili	Typical Growth
<i>Salmonella</i>	california	Typical Growth	<i>Salmonella</i>	Ochsensoll	Typical Growth
<i>Salmonella</i>	Cerro	Typical Growth	<i>Salmonella</i>	Ohio	Typical Growth
<i>Salmonella</i>	chaco	Typical Growth	<i>Salmonella</i>	oranienberg	Typical Growth
<i>Salmonella</i>	Chailey	Typical Growth	<i>Salmonella</i>	panama	Typical Growth
<i>Salmonella</i>	Chester	Typical Growth	<i>Salmonella</i>	Paratyphi A	Typical Growth
<i>Salmonella</i>	Choleraesuis	Typical Growth	<i>Salmonella</i>	Paratyphi B var L(+) tartrate +	Typical Growth
<i>Salmonella</i>	Cubana	Typical Growth	<i>Salmonella</i>	paratyphi C	No growth*
<i>Salmonella</i>	Derby	Typical Growth	<i>Salmonella</i>	Parera	Typical Growth
<i>Salmonella</i>	Diarizonae Lactose (+) H2S (-)	Typical Growth	<i>Salmonella</i>	Phoenix	Typical Growth
<i>Salmonella</i>	Dublin	Typical Growth	<i>Salmonella</i>	Pomona	Typical Growth
<i>Salmonella</i>	Enteritidis	Typical Growth	<i>Salmonella</i>	Poona	Typical Growth
<i>Salmonella</i>	essen	Typical Growth	<i>Salmonella</i>	pullorum	Typical Growth
<i>Salmonella</i>	Etterbeek	Typical Growth	<i>Salmonella</i>	Putten	Typical Growth
<i>Salmonella</i>	Ferlac	Typical Growth	<i>Salmonella</i>	Reading	Typical Growth
<i>Salmonella</i>	Florida	Typical Growth	<i>Salmonella</i>	San Diego	Typical Growth
<i>Salmonella</i>	gallinarum	Typical Growth	<i>Salmonella</i>	Senftenberg	Typical Growth
<i>Salmonella</i>	give	Typical Growth	<i>Salmonella</i>	simsbury	Typical Growth
<i>Salmonella</i>	Hadar	Typical Growth	<i>Salmonella</i>	Sloterdijk	Typical Growth
<i>Salmonella</i>	Harmelen	Typical Growth	<i>Salmonella</i>	St. Paul	Typical Growth
<i>Salmonella</i>	Heerlen	Typical Growth	<i>Salmonella</i>	stanley	Typical Growth
<i>Salmonella</i>	Heidelberg	Typical Growth	<i>Salmonella</i>	tallahassee	Typical Growth
<i>Salmonella</i>	Hilversum	Typical Growth	<i>Salmonella</i>	Tennessee	Typical Growth
<i>Salmonella</i>	Hooggraven	Typical Growth	<i>Salmonella</i>	typhi	No growth*
<i>Salmonella</i>	houten	Typical Growth	<i>Salmonella</i>	Typhimurium	Typical Growth
<i>Salmonella</i>	Illinois	Typical Growth	<i>Salmonella</i>	Typhimurium Var Copenhagen	Typical Growth
<i>Salmonella</i>	indiana	Typical Growth	<i>Salmonella</i>	Typhimurium Var O:5-	Typical Growth
<i>Salmonella</i>	indica	Typical Growth	<i>Salmonella</i>	typhisuis	Typical Growth
<i>Salmonella</i>	Infantis	Typical Growth	<i>Salmonella</i>	Urbana	Typical Growth
<i>Salmonella</i>	Inverness	Typical Growth	<i>Salmonella</i>	Vellore	Typical Growth
<i>Salmonella</i>	Javiana	Typical Growth	<i>Salmonella</i>	Worthington	Typical Growth

* *S. Newington*, typhi (6439) and paratyphi C showed typical growth on Brilliance Salmonella Agar when plated directly from blood agar (no growth in ONE Broth directly). *S. abortusovis* appeared to grow in ONE Broth but did not grow on Brilliance Salmonella Agar when plated.

CONCLUSIONS

This method comparison evaluation of the Salmonella Rapid Culture Method (Oxoid ONE Broth-Salmonella/Brilliance Salmonella Agar) clearly demonstrated that this method is equivalent to the USDA/FSIS reference method for the detection and presumptive identification of *Salmonella* spp. at spiked levels ranging from low (1 colony forming unit (CFU)/25g) to high (1.1 CFU/g) in ground beef, ground chicken, lettuce, shrimp, and shell eggs. In addition, the Salmonella Rapid Culture method was found to be superior to the reference method in incubation time and ease of identifying typical colonies. The inclusivity data demonstrated that the Salmonella Rapid Culture method (Oxoid ONE Broth Salmonella/Brilliance Salmonella Agar) detected essentially all species and serovars of *Salmonella* tested. The exclusivity data confirmed that the Salmonella Rapid Culture method was able to discriminate *Salmonella* spp. from non-*Salmonella* micro-organisms. Lot-to-lot comparability and stability data along with the ruggedness data verified that the Salmonella Rapid Culture method was robust and can provide reproducible results over a range of culture conditions (time and temperature).

© 2009, Oxoid Ltd.; copyrights to photographs held separately; contact Oxoid Ltd. for details. Photographs may not be extracted or reproduced in any way. Folio no.: LT1331/PR/09/09