



# Welcome!


Learn more about HiMedia

**HIMEDIA**<sup>®</sup>

For life is precious



# Introduction



Production  
Building

# Introduction



R & D Building




# Introduction



Production Equipment

# Introduction



# Company presentation



**ISO 9001:2015 by IQNet, Quality Austria & Bureau of Indian Standards (BIS)**



**ISO 13485-2012 by IQNet, Quality Austria**



**WHO GMP Certification by Quality Austria**



**CE Marking**



**ISO/IEC 17025:2005 by NABL**





## Microbiology

### Dehydrated Culture Media

- Powdered Media
- Granulated Media
- Encapsulated Media
- 100% Animal Free Media
- Chemically Defined Media



## HiMedia Today

- Largest range of culture media in the world
- Export to over 150 countries
- OEM supplier to reputed international brands
- Joint label with VWR USA
- Revolutionizing the field of microbiology







# Microbiology

## Microbiology Product Range

Pharmaceutical Industry

Agriculture

Brewery & Fermentation

Environmental & Sanitary

Veterinary

Cosmetics

Clinical Microbiology

Water, Food, Dairy

The screenshot displays the Himedia website's 'Applications in Microbiology' section. The URL is [www.himedialabs.com/intl/en/applications/Clinical-Microbiology/4#](http://www.himedialabs.com/intl/en/applications/Clinical-Microbiology/4#). The navigation bar includes 'Products', 'Applications in Microbiology' (highlighted), 'About Us', 'Careers', and 'Support'. The 'Applications in Microbiology' dropdown menu lists various industries: Pharmaceutical Industry, Agriculture, Brewery and Fermentation, Environmental & Sanitary, Textile and Petroleum, Cosmetics Industry, Molecular Biology Growth Media, Clinical Microbiology, Water, Food, Dairy, Veterinary, and Vaccine Production. The 'Clinical Microbiology' section is currently selected, showing a detailed view of the product range. The main heading is 'Microbial Examination in Clinical Microbiology'. Below this, there is a list of products: Diagnostic Media for Bacteria, Diagnostic Media for Fungi, Diagnostic Media for Protozoa, Susceptibility Testing, Stains, Indicators and Reagents, Biochemical and Identification, Differential Discs, General Purpose Media, and Transport Media. The text describes clinical microbiology as a branch of medicine dealing with the study of microorganisms and their role in causing diseases in human beings. It mentions that clinical microbiology includes the study of microorganisms, disease pathology, and immunology. The text also states that clinical microbiology is the most widely studied and for a wide range of products that can be used in various branches of Clinical Microbiology.



# Microbiology

## Product categories

### Dehydrated Culture Media (DCM)



- Animal Origin
- Vegetable Origin
- Chemical Origin

### Supplements



- Growth and Selective Supplements
- Vitamin & Growth Factor Supplements
- Antibiotic Mixtures
- Egg yolk Supplements
- Serum & other derivatives



# Microbiology

## Granulated Media

- More than 150 granulated medias listed in catalogue
- Combines the high throughput technology of granulation and production of dehydrated culture media
- Has similar quality attributes with several added benefits in physical
- Safe for use and less dust formation
- Only two company has Granulated media (HiMedia and Merck)







## Granulated Media (HiEncap)

- Customization of many medium in the granular form is possible
- Premeasured granulated media for 250 ml, 500ml & 1000ml
- Just drop the capsule in water and autoclave





# Microbiology

## Media

### HiVeg™

- Free of TSE/BSE risks
- Substituted with renewable sources of carbon & nitrogen
- Comparable productivity
- Lesser carbon foot print
- Ecofriendly
- More than 1500 products

### HiCynth™

- Chemically Defined Media
- Free of TSE/BSE and GMO risks
- Precise source of Carbon & Nitrogen
- Increased consistency
- Launched around 80 products
- Nutritional requirement of bacterial species can be accurately determined.

### HiCrome™

- Largest range of more than 80 such media
- Easy identification & differentiation
- Saves time
- Designed for significant bacterial identification, yeast identification
- Differentiating within a group of organisms




# Microbiology

## Why HiVeg™ Culture Media?

### Why vegetable instead of animal-based?

- Substituted with renewable sources of carbon & nitrogen
- Lesser carbon foot print
- Ecofriendly








# Microbiology

## Why HiVeg Culture Media?


### European Pharmacopoeia reference for use non-animal origin products

“When manufacturers have a choice, the use of materials from 'non TSE-relevant animal species' or non-animal origin is preferred“

Source EP 8.0,2014 <Section 5.2.8> Minimising the risk of transmitting TSE via medicinal products



# HiFill™ Test HiVeg™



MV2018G

Microbial contamination is indicated by colour change from light yellow to maroon-red

- Beneficial to pharmaceutical sectors
- In this line the HiFill™ Test Medium with the addition of MFT indicator, helps to verify the microbiological growth in aseptic production process.
- MFT Indicator in the medium is utilized by all microorganisms and the microbial contamination is indicated by colour change
- Easier method for detection of contamination with less time consumption.
- *Developed against Biomerieux equivalent : Media Fill 3P® It also contains a unique color indicator, changing from yellow to maroon-red in case of contamination.*



# Soyabean HiVeg™ Medium



**Why contaminate your sterile line with Animal Products??**

When it can be swift and safe with HiFill™ HiVeg™ Media

## MFT Media Fill Trial

Media fill studies, simulates the filling process during production and helps in detecting contamination in the production line, if any. Generally the commercial media is prepared, autoclaved and after filtering through a 0.2 micron sterilizing filter is used to investigate presence or absence of contamination. To make the process faster, efficient and safer; HiMedia provides **gamma irradiated dehydrated culture media** which can be directly used. Soyabean HiVeg™ medium sterile powder, γ-irradiated from vegetable source can be used.

### 3 - General Considerations 3-1: Scientific Principles for Minimising risk

*"When manufacturers have a choice, the use of materials from 'non TSE-relevant animal species' or non-animal origin is preferred"*

Source EP 8.0.2014 -Section 5.2.8> Minimising the risk of transmitting TSE via medicinal products

### Animal Free Media **MV011G**

Soyabean HiVeg™ Medium, Sterile Powder

### Gamma Irradiated Medium for Media Fill Trial

Also available classical Tryptone Soya Broth as per USP **MH011G / GMH011G**

\* Granulated form

Switch now to

HiMediaLaboratories™  
himedialabs.com

**HIMEDIA**  
For Life is Precious

## Media Fill : Maximum Benefits & Minimizing Risks with HiVeg™ Gamma Irradiated TSB.



With the spurt in number of BSE symptoms across global bovine population & its exhibit CJD in humans concerns were raised about bovine origin products.

Elimination of BSE/TSE Risk can be achieved by use of raw material from right origin & right parts of the animal. Definition of Risk Categories by EU:

- Category A: High Infectivity (e.g. brain, spinal cord)
- Category B: Moderate infectivity (e.g. spleen, lung, liver)
- Category C: No infectivity found (e.g. milk, bile, skeletal muscle, heart, skin)

HiMedia only sources from risk category 'C' for its products. Moreover as per the Definition of Geographical BSE Risk by EU, raw material sourced from India has no listings. In spite of such a proven track record of quality, a step further to provide more secure process HiVeg™ culture media was launched. Both USP & EP preferred or recommend that alternative, non-animal source ingredients be substituted for animal-source ingredients whenever possible.

The risk of Mycoplasma is always lurking in the raw material. Moreover Mycoplasma can move through 0.2 mm filters & Reach high titers ( $10^7$  —  $10^8$  cfu/ml) without producing pH changes or media turbidity proving itself as invisible threat. In such cases a prudent step ahead to provide medium quality assurance is to provide γ-irradiated TSB.

γ-irradiation does not affect product performance, and results in a Contaminant-free material, this has been evaluated by comparative studies on growth performance of pharmacopoeia listed pathogens. Thus HiVeg™ γ-irradiated TSB is the choice of a prudent quality system.

Introduced gamma irradiated HiFill™ Test Medium recommended for the evaluation of sterility in manufacturing process for easy detection of contamination. The medium is designed with TSB containing an MFT indicator wherein the colour change is from yellow to pink red.

#### Reference:

- \* The USP Perspective to Minimize the Potential Risk of TSE- Infectivity in Bovine-derived Articles Used in the Manufacture of Medical Products; with Ian DeVaux and Roger Dabbeh. Pharmacopoeial Forum. 30(5):1911-1921.2004
- \* European Pharmacopoeia (Supplement 6.3), 2008, European Department, for the Quality of Medicines

HiMedia Laboratories Pvt. Ltd.  
www.himedialabs.com

#### CORPORATE OFFICE -

A-516, Swastik Disha Business Park, Via Vadhvani Indl Est, LBS Marg, Mumbai - 400 085, India.  
Tel.: 90-91-22-6147 1919 / 2500 3747 | Fax: 00-91-22-6147 1920 / 2500 5764 | Email: hdmexpo@himedialabs.com

HiMedia No.	Product Range for Media Fill trials
M011G-500G M011G-2.5KG M011G-5KG	Soyabean Casein Digest Medium, Sterile Powder γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
MV011G-500G MV011G-2.5KG MV011G-5KG	Soyabean <b>HiVeg</b> Medium, Sterile Powder γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
GMV011G-500G	Soyabean <b>HiVeg</b> Medium, Granulated, Sterile γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
MH011G-500G	Soyabean Casein Digest Medium, Sterile powder γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
GMH011G-500G	Soyabean Casein Digest Medium, Granulated, Sterile γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
M1655G-500G M1655G-2.5KG	Soyabean Casein Digest Medium w/ Mannitol, Sterile Powder γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process. It can also be used for cultivation of a wide variety of microorganisms.
M1655G-500G M1655G-2.5KG M1655G-5KG	Soyabean Casein Digest Medium w/ BCP, Sterile Powder γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process.
M010G-500G M010G-2.5KG M010G-5KG	Alternative Thioglycollate Medium, Sterile Powder γ-irradiated sterile powder recommended for evaluation of sterility in manufacturing process.
MV010G-500G MV010G-2.5KG MV010G-5KG	Alternative Thioglycollate <b>HiVeg</b> Medium, Sterile Powder γ-irradiated sterile powder recommended for evaluation of sterility in manufacturing process.
MU010G-500G MU010G-2.5KG MU010G-5KG	Alternative Thioglycollate Medium, Sterile Powder γ-irradiated sterile powder recommended for evaluation of sterility in manufacturing process in accordance with USP.
M2018G-500G	HiFill™ Test Medium γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process for easy detection of contamination by Media Fill Test.
MV2018G-500G	HiFill™ Test <b>HiVeg</b> Medium γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process for easy detection of contamination by Media Fill Test.
MCD2018G-500G	HiFill™ Test HiCynth™ Medium γ-irradiated sterile powder recommended for the evaluation of sterility in manufacturing process for easy detection of contamination by Media Fill Test.
RM565G-5KG RM565G-50KG	Lactose monohydrate, Sterile (γ irradiated sterile powder)
RM565GT-5KG	Lactose monohydrate, Sterile Powder (γ irradiated Triple Pack)
RM570G-5KG RM570G-50KG	D-Mannitol, A.R. sterile (γ irradiated)

Bulk packing are available on request for all above products.




**HIMEDIA**  
For Life is Precious

**HIMEDIA**  
For life is precious



# HiCrome™ - Chromogenic Media

- Largest range of more than 80 such media
- Easy identification & differentiation
- Saves time
- Designed for significant bacterial identification, yeast identification
- Differentiating within a group of organisms
- HiVeg™ HiCrome™ Media also available



M1991I suggest for detection of *Escherichia coli* and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014 .



**HIMEDIA®**

For Life is Precious

**HiCrome™**

**Chromogenic Coliform**

**Agar (CCA Agar) (M1991I)**

**in accordance with ISO 9308-1:2014**

**Use :**

Recommended for detection of *Escherichia coli* and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014

**Advantages :**

- ▶ Simultaneous detection of *Escherichia coli* and total coliforms.
- ▶ Excellent recovery from water with low count
- ▶ Easy interpretation due to colour differentiation
- ▶ Available in dehydrated and ready prepared form

**Principle :**

- ▶ Mixture of three chromogens, easily detects  $\beta$ -galactosidase and  $\beta$ -glucuronidase enzymes.
- ▶ IPTG is added to enhance colour detection.
- ▶ L-Tryptophan - improved indole reaction (helps in easy detection).
- ▶ Tergitol-7 for selectivity - Gram positive bacteria inhibited

**Interpretation :**

- ▶ *E.coli* - dark blue to violet
- ▶ Other coliforms - pink to red
- ▶ Other gram negative bacteria - colourless
- ▶ Organisms with weak glucuronidase activity but no beta galactosidase activity produce light blue to turquoise colonies.



*Escherichia coli*  
ATCC 25922

*Salmonella enteritidis*  
ATCC 13076

*Escherichia coli*  
ATCC 25922

*Salmonella enteritidis*  
ATCC 13076

*Escherichia coli*  
ATCC 25922

**M1991I CCA Agar**

\* Formerly known as Enterobacter aerogenes

# Folic Acid Casei Medium

Available Globally

HiMedia's  
Folic Acid Casei Medium, M543

For the microbiological assay of folic acid in blood serum

Using *Lactobacillus casei* ATCC 7469 as the test organism.

\* Only HiMedia has this Medium







# Folic Acid Casei Medium

## How does it work?

Cell density of *Lactobacillus casei* is dependent on the folic acid content in the growth medium. By measuring the turbidity of the growth medium containing the sample, the folic acid content of the sample can be determined.


Users of Folic Acid Casei Medium are found in the clinical, pharmaceutical and food industry.



The folic acid concentration of the sample is determined by means of a standard curve.

# LQ104 - Rappaport Vassiliadis Medium

- Sterile Medium in Ready prepared form
- Convenient Pack sizes of 5ML, 10ML, 20ML
- For direct inoculations of samples
- Recommended for Selective Enrichment Medium for Salmonella
- Sterilized by autoclaving at 115°C as per validated cycle
- Sterility assurance level assured by using BI strips
- Formula in compliance with as per Pharmacopoeias USP, EP, BP, JP
- For use in Pharmaceuticals, Clinical and Food testing



# LQ104 - Rappaport Vassiliadis Medium

## 2.6.13. Test for specified micro-organisms

EUROPEAN PHARMACOPOEIA 6.3

Table 2.6.13.1 – Growth promoting, inhibitory and indicative properties of media

	Medium	Property	Test strains
Test for bile-tolerant gram-negative bacteria	Enterobacteria enrichment broth-Mossel	Growth promoting	<i>E. coli</i> <i>P. aeruginosa</i> <i>S. aureus</i>
		Inhibitory	
	Violet red bile glucose agar	Growth promoting + indicative	<i>E. coli</i> <i>P. aeruginosa</i>
Test for <i>Escherichia coli</i>	MacConkey broth	Growth promoting	<i>E. coli</i>
		Inhibitory	<i>S. aureus</i>
	MacConkey agar	Growth promoting + indicative	<i>E. coli</i>
Test for <i>Salmonella</i>	Rappaport Vassiliadis <i>Salmonella</i> enrichment broth	Growth promoting	<i>Salmonella enterica</i> ssp. <i>enterica</i> serotype typhimurium or <i>Salmonella enterica</i> ssp. <i>enterica</i> serotype abony
		Inhibitory	<i>S. aureus</i>
	Xylose, lysine, deoxycholate agar	Growth promoting + indicative	<i>Salmonella enterica</i> ssp. <i>enterica</i> serotype typhimurium or <i>Salmonella enterica</i> ssp. <i>enterica</i> serotype abony
		Indicative	<i>E. coli</i>
Test for <i>Pseudomonas aeruginosa</i>	Cetrimide agar	Growth promoting	<i>P. aeruginosa</i>
		Inhibitory	<i>E. coli</i>
Test for <i>Staphylococcus aureus</i>	Mannitol salt agar	Growth promoting + indicative	<i>S. aureus</i>
		Inhibitory	<i>E. coli</i>
Test for clostridia	Reinforced medium for clostridia	Growth promoting	<i>Cl. sporogenes</i>
	Columbia agar	Growth promoting	<i>Cl. sporogenes</i>
Test for <i>Candida albicans</i>	Sabouraud dextrose broth	Growth promoting	<i>C. albicans</i>
	Sabouraud dextrose agar	Growth promoting + indicative	<i>C. albicans</i>

# Rappaport Vassiliadis Medium

## MH1491 & GMH1491

- Recommended for selective enrichment of Salmonella from pharmaceutical products
- In accordance with the Tests for specified organisms-Nonsterile products (formerly microbial limit testing) by harmonized methodology of USP/EP/BP/JP.
- Advised to sterilize by autoclaving at 115°C as per validated cycle as given in pharmacopoeia
- Available Pack sizes of 100G, 500G, 2.5KG , 5KG
- g/l of HiMedia MH1491 & GMH1491 is 27.11\*  
while g/l of Merck 1.07666.0500 is 42.5

\*Loss of water molecules of magnesium chloride, hexahydrate 29.0 g in medium is accounted in dehydrated medium



# Rappaport Vassiliadis Medium

## MH1491 & GMH1491

EUROPEAN PHARMACOPOEIA 6.3

2.6.13. Test for specified micro-organisms

Neutral red	30.0 mg	Heat to boiling for 1 min with shaking. Adjust the pH so that after sterilisation it is $7.2 \pm 0.2$ at $25^\circ\text{C}$ . Sterilise in an autoclave using a validated cycle.
Crystal violet	1 mg	
Purified water	1000 ml	
Adjust the pH so that after sterilisation it is $7.1 \pm 0.2$ at $25^\circ\text{C}$ . Boil for 1 min with constant shaking then sterilise in an autoclave using a validated cycle.		
<b>Rappaport Vassiliadis <i>Salmonella</i> enrichment broth</b>		
Soya peptone	4.5 g	
Magnesium chloride hexahydrate	29.0 g	
Sodium chloride	8.0 g	
Dipotassium phosphate	0.4 g	
Potassium dihydrogen phosphate	0.6 g	
Malachite green	0.036 g	
Purified water	1000 ml	
Dissolve, warming gently. Sterilise in an autoclave using a validated cycle, at a temperature not exceeding $115^\circ\text{C}$ . The pH is to be $5.2 \pm 0.2$ at $25^\circ\text{C}$ after heating and autoclaving.		
<b>Reinforced medium for clostridia</b>		
Beef extract	10.0 g	
Peptone	10.0 g	
Yeast extract	3.0 g	
Heat to boiling for 1 min with shaking. Adjust the pH so that after sterilisation it is $7.4 \pm 0.2$ at $25^\circ\text{C}$ . Sterilise in an autoclave using a validated cycle.		

2. Methods of analysis



## Technical Data

### Rappaport Vassiliadis *Salmonella* Enrichment Broth

MH1491

#### Intended use

Rappaport Vassiliadis *Salmonella* Enrichment Broth is recommended for selective enrichment of *Salmonella* species from pharmaceutical products in accordance with the microbial limit testing by harmonized methodology of USP/EP/BP/IP.

#### Composition\*\*

Ingredients	Gms / Litre
Soya peptone	4.500
Sodium chloride	8.000
Dipotassium phosphate	0.400
Potassium dihydrogen phosphate	0.600
Magnesium chloride, hexahydrate	29.000
Malachite green	0.036
pH after sterilization ( at $25^\circ\text{C}$ )	$5.2 \pm 0.2$

\*\*Formula adjusted, standardized to suit performance parameters

#### Directions

Suspend 27.11 grams of dehydrated medium(the equivalent weight of dehydrated medium per litre) in 1000 ml purified/ distilled water. Heat if necessary to dissolve the medium completely. Dispense as desired into tubes and sterilize by autoclaving at  $115^\circ\text{C}$  as per validated cycle

#### Principle And Interpretation

Rappaport Vassiliadis *Salmonella* Enrichment Medium is designed according to the revised formulation by Van Schothorst et al (1) and is recommended for the selective enrichment of *Salmonellae* from pharmaceutical products. This medium can

**HIMEDIA**

For life is precious



## **MEDIA FOR FOOD & DAIRY APPLICATIONS**



## Few Typical Media used for Food & Dairy products Testing

### For Total Bacterial Count

- M091 Plate Count Agar
- M1884 DEV Nutrient Agar

### For Yeasts & Moulds

- M640 Rose Bengal Chloramphenicol Agar

### For Meat & Poultry Antimicrobial Inhibitor Test Agar

- M1631 (pH 6.0), M1601 (pH 7.2), M1632 (pH 8.0)





## Few Typical Media used for Food & Dairy products Testing

### For Lactobacilli

- M1163 MRS Agar
- M641 Lactobacillus MRS Agar

### For Bifidobacteria

- M1734 Bifidobacterium Selective Count Agar Base

### For Osmophilic Organism

- M594 MY 40 Agar (Osmophilic Organisms)

### For Thermophilic Flat Sour Sporeformers



- M1104- M-Dextrose Tryptone Broth



# Microbiology

## Common Culture media for suspending, sample preparation and enrichment

Code	Product
GM525	Ringer salt solution
M028, M028I	Peptone Water
M1494, GM1494I	Buffered Peptone Water
M461	Phosphate Buffer
M618, M618I	Alkaline Peptone Water
M080	Lauryl Sulphate Broth
M149	Cooked Meat Medium

**Technical Data**

---

**Buffered Peptone Water, Granulated**GM1494I

Buffered Peptone Water, granulated is used as pre-enrichment medium for increasing the recovery of injured *Salmonella* species from foods prior to selective enrichment and isolation. The composition and performance criteria of this medium are as per the applications laid down in ISO 6579-2002.

**Composition\*\***

Ingredients	Gms / Litre
Enzymatic digest of casein	10.000
Sodium chloride	5.000
Disodium hydrogen phosphate, 12H <sub>2</sub> O	9.000
Potassium dihydrogen phosphate	1.500
Final pH ( at 25°C)	7.0±0.2

\*\*Formula adjusted, standardized to suit performance parameters

**Directions**

Suspend 20.07 grams (the equivalent weight of dehydrated medium per litre) in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

# Buffered Peptone Water

Universal Medium for use as diluent and preparation of samples in food industries

- g/l of M614, MV614, MCD614: 20 g/l
- g/l of M1494I, GM1494I (as per ISO), MV1494I, MCD1494I: 20.07 g/l

*9 g  $\text{Na}_2\text{HPO}_4 \cdot 12 \text{H}_2\text{O}$  on dehydration is equivalent to 3.57 g  $\text{Na}_2\text{HPO}_4$  anhydrous*

- g/l of Merck 1.07228: 25.5 g/l



# Buffered Peptone Water



## Technical Data Sheet

### GranuCult™

### Buffered Peptone Water

acc. ISO 6579, ISO 21528, ISO 22964, FDA-BAM  
and EP

Ordering number: 1.07228.0500 / 1.07228.5000

For the preliminary non-selective enrichment of bacteria, particularly pathogenic *Enterobacteriaceae* such as *Salmonella* and *Cronobacter*, from food and animal feed, water and other materials.

This culture medium complies with the specifications given by EN ISO 6579, EN ISO/IDF 6579-1, EN ISO 6785 I IDF 93, EN ISO 19250, EN ISO 21528-1, ISO/TS 22964 I IDF/DRM 210, FDA-BAM, APHA and EP.

#### Mode of Action

The broth is rich in nutrients and produces high resuscitation rates for sublethally injured bacteria and intense growth. The phosphate buffer system prevents bacterial damage caused by changes in the pH of the medium. Peptone including enzymatic digest of casein acts as a source of carbon, nitrogen, vitamins and minerals whilst sodium chloride maintains the osmotic balance.



The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

1 of 5



#### Typical Composition

Specified by ISO 6579, ISO/IDF 6579-1, ISO 19250, ISO 21528, ISO 22964		Specified by FDA-BAM M182		Specified by EP 2.8.31, ISO 6785 I IDF 93		GranuCult™ Buffered Peptone Water acc. ISO 6579, ISO 21528, ISO 22964, FDA-BAM and EP	
Enzymatic Digest of Casein*	10 g/l	Peptone	10 g/l	Peptone	10 g/l	Peptone (Includes Enzymatic Digest of Casein)	10 g/l
NaCl	5 g/l	NaCl	5 g/l	NaCl	5 g/l	NaCl	5 g/l
Na <sub>2</sub> HPO <sub>4</sub> x 12 H <sub>2</sub> O	9 g/l	Na <sub>2</sub> HPO <sub>4</sub> **	3.5 g/l	Na <sub>2</sub> HPO <sub>4</sub> x 12 H <sub>2</sub> O	9 g/l	Na <sub>2</sub> HPO <sub>4</sub> x 12 H <sub>2</sub> O	9 g/l
KH <sub>2</sub> PO <sub>4</sub>	1.5 g/l	KH <sub>2</sub> PO <sub>4</sub>	1.5 g/l	KH <sub>2</sub> PO <sub>4</sub>	1.5 g/l	KH <sub>2</sub> PO <sub>4</sub>	1.5 g/l
Water	1000 ml/l	Water	1000 ml/l	Water	1000 ml/l	Water	n/a
pH at 25 °C	7.0 ± 0.2	pH at 25 °C	7.0 ± 0.2	pH at 25 °C	7.0 ± 0.2	pH at 25 °C	7.0 ± 0.2

\* ISO/IDF 6579-1 specifies: Peptone - for example, enzymatic digest of casein.

\*\* 3.57 g Na<sub>2</sub>HPO<sub>4</sub> anhydrous is equivalent to 9 g Na<sub>2</sub>HPO<sub>4</sub> x 12 H<sub>2</sub>O

#### Preparation

**Dissolve 25.5 g in 1 l of purified water.** If desired dispense into smaller vessels and autoclave 15 min at 121 °C.

The prepared medium is clear and yellowish. The pH value at 25 °C is in the range of 6.8-7.2.

#### Experimental Procedure and Evaluation

Depend on the purpose for which the medium is used.

Incubate the inoculated broth under aerobic conditions, e.g. acc. to EN ISO 6579 36-38 °C for 16-20 h, acc. to EN ISO/IDF 6579-1 at 34-38 °C for 16-20 h.

Transfer material from the resulting culture to a selective enrichment culture medium following the method given by the appropriate standard.

According to EN ISO/IDF 6579-1, it is permissible to store the pre-enriched sample after incubation at +2 to +8 °C for a maximum of 72 h.

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Binding laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not release our customers of their own responsibility for observing the suitability of our products for the envisaged purpose.  
Merck, GranuCult, ReadyCult and ReadySorb are trademarks of Merck KGaA, Darmstadt, Germany. Anasorb, AnasorbCult, Redibest, ChromoCult, Cluspath, FluoroCult, ReadyCult, ReadySorb and ReadyTag are registered trademarks of Merck KGaA, Darmstadt, Germany.

2 of 5

**HIMEDIA**

For life is precious

# Buffered Peptone Water



## Technical Data

### Buffered Peptone Water

M1494I

#### Intended use

Buffered Peptone Water is used as pre-enrichment medium for increasing the recovery of injured *Salmonella* species from foods prior to selective enrichment and isolation. The composition and performance criteria of this medium are as per the applications laid down in ISO 6579-2017, ISO 6887 and ISO 21528-2017.

#### Composition\*\*

Ingredients	Gms / Litre
Tryptone #	10.000
Sodium chloride	5.000
Disodium hydrogen phosphate.12H <sub>2</sub> O	9.000
Potassium dihydrogen phosphate	1.500
Final pH ( at 25°C)	7.0±0.2

\*\*Formula adjusted, standardized to suit performance parameters

# Equivalent to Enzymatic digest of casein

#### Directions

Suspend 20.07 grams (equivalent weight of dehydrated medium) in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

#### Principle And Interpretation

Microorganisms that are subjected to environmental stresses may become structurally or metabolically damaged or injured. These microorganisms are unable to replicate in selective environments. Therefore these injured organisms must be resuscitated or permitted to repair the damage by incubation in an appropriate, non-selective environment (1). Edel and Kampelmacher (2) noted that sublethal injury to *Salmonellae* may occur in many food preservation processes. Enriching injured cells in Lactose



## Technical Data

### Buffered Peptone Water, Granulated

GM1494I

Buffered Peptone Water, granulated is used as pre-enrichment medium for increasing the recovery of injured *Salmonella* species from foods prior to selective enrichment and isolation. The composition and performance criteria of this medium are as per the applications laid down in ISO 6579-2002.

#### Composition\*\*

Ingredients	Gms / Litre
Enzymatic digest of casein	10.000
Sodium chloride	5.000
Disodium hydrogen phosphate, 12H <sub>2</sub> O	9.000
Potassium dihydrogen phosphate	1.500
Final pH ( at 25°C)	7.0±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Directions

Suspend 20.07 grams (the equivalent weight of dehydrated medium per litre) in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

#### Principle And Interpretation

Microorganisms that are subjected to environmental stresses may become structurally or metabolically damaged or injured. These microorganisms are unable to replicate in selective environments. Therefore these injured organisms must be resuscitated



For life is precious



## **MEDIA FOR WATER APPLICATIONS**





# Microbiology

## Media used for Water Testing

### For Enumeration

- M091 Plate Count Agar

### For Heterotrophic Plate Count

- M962 R2A Agar

### For Total Count & Gelatin Liquifying Organism

- M1609 Gelatin DEV Agar

### For Yeasts & Moulds

- M063 Sabouraud Dextrose Agar

### For Aerobic Bacteria

- M290 Soyabean Casein Digest Agar

**Granulated**

**HIMEDIA** REF GMH290-500G Net Content 500g

Storage :  
Hygroscopic; Keep tightly closed away from bright light. On receipt store at 30°C

© Equivalent to MU290/ME290/MM290/M290B

**Soybean-Casein Digest Agar (Casein Soyabean Digest Agar), Granulated**

Directions:  
Suspend 40.0 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes or as per validated cycle. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Use  
Recommended as a general purpose medium for cultivation of a wide variety of microorganisms from pharmaceutical products in accordance with the harmonized method of USP/EP/BP/JP/CP (Medium 2).

**\*\*Standard Formula**

Ingredients	Gms/litre
Tryptone \$	15.00
Soya peptone*	5.00
Sodium chloride	5.00
Agar	15.00

\*pH after sterilization 7.3±0.2  
\*pH can also be measured after sterilization at 25°C  
\*\*Formula adjusted, standardized to suit performance parameters  
\$ Pancreatic digest of casein  
\*Papain digest of soyabean meal

\*For More Information Refer Technical Data

COUNTRY OF ORIGIN-INDIA  
Company certified for ISO 9001:2008, ISO 13485:2003, WHO GMP

**HiMedia Laboratories Pvt. Ltd.**  
Reg.off. : 23, Vadhani Ind. Est., LBS Marg, Mumbai-400066, India.  
Works : B/4-6, M.I.D.C., Dindori, Nashik, India.  
Customer Care No.: 00-91-22-6116 9797  
Email: techhelp@himediailabs.com

FEB-2022 LOT 0000282057

8 902729 914632

Received \_\_\_\_\_ Opened \_\_\_\_\_  
Disposal : User must ensure safe disposal by autoclaving and/or by incineration used or unusable preparations of this product and derivatives thereof on completion of work to avoid contagion.

Expected performance during specified expiry period when material is duly maintained in the original powder form.



# Microbiology

## Culture media for detecting coliforms in water and waste water

### Coliforms

#### MPN Tests

- M080 Lauryl Sulphate Broth
- M083 MacConkey Broth Purple w/BCP
- M458 Violet Red Bile Broth

### Membrane Filter Technique

- M1106 M-Endo Agar LES
- M1111 M-FC Broth Base
- M1023 M-Lauryl Sulphate Broth
- M1066 M-Tergitol 7 Agar Base
- M1991I HiCrome™ Chromogenic Coliform Agar (CCA Agar) ISO 9308-1:2014.



# Microbiology

## Chemical Testing

### Multi Parameter Water Testing Kit - WT015

Comprehensive Lab-Free Qualitative and Quantitative Chemical Water Testing in a Single Kit

- Fluoride • Nitrate • Iron • Residual (Free) chlorine
- Chloride • Total hardness • Turbidity test • pH test



Bureau of Indian standards have set the requirements for essential and desirable characteristics to be tested for ascertaining the suitability of water in IS 10500-1991.

WT015 offered by HIMEDIA is a Multiparameter water testing kit determining levels of fluoride, nitrate, iron, residual (free) chlorine, chloride and total hardness besides measuring turbidity and pH.

#### Kit contents :

##### Type of test


Fluoride  
Nitrate  
Iron  
Residual (Free) chlorine  
Chloride  
Total hardness  
Turbidity test  
pH test

##### Range

0.0-5.0 mg/L (ppm) as Fluoride  
0.0-250 mg/L (ppm) as Nitrate ( $\text{NO}_3$ )  
0.0 - 2.0 mg/L (ppm) as Iron  
0.0-2.0 mg/L (ppm) as free chlorine  
10-200 mg/L (ppm) and 50-1000 mg/L (ppm) as Chloride  
25-600 mg/L (ppm) as  $\text{CaCO}_3$   
10-600 NTU, standards of 10 & 25 NTU  
pH test strips of range 6.5 to 9.0.

### Typical Tests Employed for Water Testing

Convenient and Easy Touch Plates



HiTouch *E. coli* Coliform Count Flexi Plate - FL022



Baird - Parker Agar - M043  
*Staphylococcus aureus* (ATCC 6538)



MacConkey Broth - M007

1. Uninoculated control
2. *Escherichia coli* (ATCC 25922)
3. *Enterobacter aerogenes* (ATCC 13048)
4. *Neisseria meningitidis* (ATCC 13883)
5. *Staphylococcus aureus* (ATCC 25923)
6. *Enterococcus faecalis* (ATCC 29212)

### Hi-Dip Slides for Lab-Free Testing



Rose Bengal Agar - HD008  
fungal growth observed




## CULTURE MEDIA FOR BREWERY & FERMENTATION








# Microbiology





# Microbiology



## Brewery & Fermentation

The Code Nos. Marked with Green Leaves are also available as HiVeg™ Media

Code	Product
<b>Controlling Fermentation Processes</b>	
<b>1. Agar Media</b>	
M1225	Acetate Agar
*M058	Actidione Agar Base w/o Actidione*
*M400	Actidione Agar w/ Actidione*
M642	Lysine Medium Base
FD123	50% Potassium Lactate
M253	Malt Agar
M1605	Modified NBB Agar (Pedi-Lacto Selective Agar base)
M1815	Pedi Lacto Selective Agar Base Modified
M1821	Pedi Lacto Selective Broth Base (Modified NBB Broth)
*M828	Sucrose Agar for Brewery Isolates
M1483	Universal Beer Agar, Modified
*M415	Universal Beer Agar (UB Agar)
*M115	WL Nutrient Medium
*M1060	WL Differential Agar
*M129	Wort Agar
<b>2. Liquid Media</b>	
M978	Anaerobic Fermentation Medium Base
*M885	Andrade Peptone Water
*M909	Andrade Peptone Water w/ HM Extract
M401	Blue Agar
*M676	Bromo Cresol Purple Broth Base (Yeast Fermentation Broth Base)
*M351	CHO Medium Base
M1662	Enteric Fermentation Base
M1589	Glucose Agar
*M860	Glucose Broth
M2000	Modified WL Nutrient Medium
*M028	Peptone Water
M1332	Universal Liquid Medium
*M1060	WL Differential Agar
*M410	WL Differential Broth
*M050	WL Nutrient Broth
*M333	Wort Broth
<b>Brewing</b>	
<b>1. Coliforms</b>	
M1905	Bromo Cresol Purple Agar w/ Lactose
*M026	Fluid Lactose Medium
*M081	MacConkey Agar w/ 0.15% Bile Salts, CV and NaCl

Code	Product
*M083	MacConkey Broth Purple w/ BCP
*M1316	Super Broth
<b>2. Lactobacilli / Yeasts</b>	
M1384	Lactic Acid Bacteria Selective Broth Base (Raka-Ray No. 3 Broth Base)
FD055	Lactic Supplement
FD155	Lactic Supplement Modified
*M1087	Lactic Bacteria Differential Agar
*M1086	Lactic Bacteria Differential Broth
*M1333	Lee's Multidifferential Agar
M609	Litmus SM Broth
M642	Lysine Medium Base
FD123	50% Potassium Lactate
M1846	MGYP Agar with Copper
M1504	M-BCG Yeast and Mould Agar
M1130	M-BCG Yeast and Mould Broth
M1741	M-BCG Yeast and Mould Broth, Modified
*M641	MRS Agar (Lactobacillus MRS Agar)
M1923	MRS Agar w/ pH 5.5
M639	Oxytetra Glucose Yeast Agar Base (OGYE Agar Base)
FD032	Oxytetra Selective Supplement
M1815	Pedi Lacto Selective Agar Base Modified
M1821	Pedi Lacto Selective Broth Base (Modified NBB Broth)
*M130	Rogosa SL Agar
*M640	Rose Bengal Chloramphenicol Agar
*M1331	SCHWARZ Differential Medium
M829	Tomato Juice Medium Base
FD098	Lactobacilli Supplement
*M415	Universal Beer Agar (UB Agar)
M1483	Universal Beer Agar, Modified
M1332	Universal Liquid Medium
*M1060	WL Differential Agar
*M410	WL Differential Broth
*M050	WL Nutrient Broth
*M115	WL Nutrient Medium
*M129	Wort Agar
*M424	Yeast Malt Agar (YM Agar) (ISP Medium No. 2)


Literature code: TL 25001, Culture media for brewery & fermentation (2022)



# Microbiology

## Ezy MIC™ Strip

- A single step Method For Antimicrobial Susceptibility Test
- According to EUCAST and CLSI
- Cellulose paper strip, printed with 30 MIC values on either side
- Coated with antibacterial agent in a concentration gradient manner.
- To be placed on seeded agar surface and incubate, shows accurate MIC value on incubation.
- Spare capacity for immediate ramp up of production
- Shows accurate MIC value on incubation





# Molecular Biology

HiMedia's has led to Progressions for

- Diagnostics
- Therapeutics
- Genetic screening

**HIMEDIA**<sup>®</sup>


For life is precious





# Molecular Biology

## Target Industry for Molecular Biology Products





# Molecular Biology

**HiMedia offers a Complete end to end Solution for Automation in Molecular Biology**





# Cross-Reference Search



**Thank You**



For life is precious



GenoChem world is our official distributor in Spain



For life is precious