



# Microbiology Focus

Volume 2.1, 2010



## "Fusariosaurus" - The Discovery of a New Organism?



*Fusariosaurus in Wheat*

Fusariosaurus.....	p2
Photo Competition.....	p2
<i>Fusarium</i> .....	p3
<i>Campylobacter</i> Overview .....	p5
More about <i>Campylobacter</i> .....	p6

# Fusariosaurus (cover)

By Jvo Siegrist, Product Manager Microbiology.... [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)

## Who says that scientists don't have a sense of humor?

The interesting and intriguing microscopic image on the front cover of this issue of *Microbiology Focus* is the winning photo from the Summer 2009 Fluka Microbiology Photo Competition. It captures a fungus of the species *Fusarium culmorum* with a notable example of a typical mycel-forming conidia with spores that have already germinated. Germ tubes, visible on image, develop later to abundant mycelium. Especially amusing was the description provided by Dr. Ela Suchowilska for that unique photo. She titled the image "Fusa-

riosaurus", based on the Pterosaurs which dominated the skies during the mesozoic era (late Triassic to the end of the Cretaceous period; 220 to 65.5 million years ago).

Dr. Ela Suchowilska is an expert on *Fusarium* and wrote an interesting article in this issue of *Microbiology Focus*. To learn more about this noteworthy topic, see "*Dangerous Fungi of the Genus Fusarium*".



Figure 1: Pterosaurs



Figure 2: "Fusariosaurus"

## Winning Entries of Photo Competition

Last summer we held our Fluka Microbiology Photo Competition and received some fascinating images from around the world. The photographic entry winning Best of Show is featured in the previous article. The entries winning second, third, and fourth place

prizes are showcased below. The aim of the competition was to encourage microbiologists to promote their work, with the condition that entries should illustrate any microorganisms (living or dead) or a microbiologist in action at work. Sigma-Aldrich would like to thank

all who entered our competition and also our independent jury members Dr. Antje Breitenstein (Scanbec GmbH), Prof. Dr. Corinne Gantenbein (ZHSW), and Prof. Dr. Mohammad Manafi (Medical University of Vienna).



### Second Place Winner:

This is the 1000X Microscopic view of two isolates from Indian Ocean nodules capable of metal oxidation from its sulphides and oxides.

#### Abhilash Pillai

National Metallurgical Laboratory, India



### Third Place Winner:

Typical colony of *Mycobacterium tuberculosis* seen under a microscope with 10x magnification.

#### Martin Anandi

Institute of Tropical Medicine, Antwerp, Belgium



### Fourth Place:

"Jogging with germinating *Fusarium culmorum* conidia"

#### Ela Suchowilska

University of Warmia and Mazury in Olsztyn, Poland

# Dangerous Fungi of the Genus *Fusarium*

By Dr. Ela Suchowilska, Department of Plant Breeding & Seed Production, University of Warmia and Mazury in Olsztyn-Poland....  
ela.suchowilska@uwm.edu.pl

Fungi are heterotrophic organisms that were previously included in the Kingdom Plantae, and today are classified as a separate kingdom, MYCOTA. The Kingdom Fungi comprises unicellular organisms (e.g. yeasts), molds (e.g. fungi of the genera *Aspergillus*, *Penicillium*, *Fusarium*) and the Basidiomycetes, including the well-known edible and poisonous cap mushrooms. Molds are composed of filamentous hyphae that form an interconnected network known as mycelium, on which species-specific spores of the imperfect stage (most often asexual conidiospores) are produced, sometimes accompanied by sexual spores. Their typical size is from several to several dozen micrometers, so they are referred to as “microscopic fungi”, since their structures can only be seen under the microscope.

Fungi of the genus *Fusarium* are of great economic significance due to their widespread occurrence and high pathogenicity to all crop species grown throughout the world. This is of particular importance in the case of cereals, as head and kernel infections may drastically reduce grain yield and quality. *Fusarium* head blight (FHB), caused by approximately 16 species of the genus *Fusarium* (Parry et al. 1995), can be devastating, with an overall decrease in yield reaching 70%. In addition, toxic secondary metabolites produced by *Fusarium* species can be found in FHB-affected grain. The most common among them are trichothecenes (mostly type B, deoxynivalenol and nivalenol) (Ueno 1977, Joffe 1986), zearalenone, fumonisins and other fusariotoxins (Nelson 1993) (Table 1). The above substances exhibit very strong phytotoxic and zootoxic effects. The development of FHB-resistant varieties is an impor-

tant consideration and the main goal of numerous breeding programs across the world. This disease has been increasing in incidence and severity in recent years, due to the implementation of simplified crop production technologies including intensive monoculture systems and simplified crop rotation (in particular with respect to wheat and maize), the lack of effective fungicides for *Fusarium* control and the absence of resistant varieties (Mesterhazy 2002). Therefore, the current direction of research studies is to develop new, *Fusarium*-resistant varieties that could be relied on as a safe source of food for humans and animals.

Mycotoxins	<i>Fusarium</i> species
Type A trichothecenes (e.g. T2-toxin, HT-2 toxin, diacetoxyscirpenol) and type B trichothecenes (e.g. deoxynivalenol and its acetyl derivatives, nivalenol)	<i>Fusarium sporotrichioides</i> , <i>F. poae</i> , <i>F. graminearum</i> , <i>F. culmorum</i> , <i>F. crookwellense</i> , <i>F. acuminatum</i>
Fumonisin	<i>F. moniliforme</i> , <i>F. proliferatum</i>
Zearalenone	<i>F. graminearum</i> , <i>F. culmorum</i> , <i>F. crookwellense</i>
Cyclodepsipeptides (enniatins, beauvericin)	<i>F. avenaceum</i> , <i>F. sporotrichioides</i>

**Table 1:** Most common mycotoxins produced by *Fusarium* species, found in infected grain (Nelson 1993).

The programs aimed to breed *Fusarium*-resistant varieties involve determining pathogen-host interactions, the types and mechanisms of resistance, as well as genes responsible for resistance to fungi of the genus *Fusarium*. Due to the low host specificity of these pathogens, FHB resistance is conferred by multiple genes (Parry et al., 1995). Mesterhazy (2002) reported five different types of resistance to FHB, i.e. resistance to initial infection (type I), resistance to the spread of the pathogen (type II), ability to degrade mycotoxins (type III), resistance to kernel infection (type IV), and tolerance (type V). FHB resistance may be conditioned phytochemically by the presence of phenolic compounds inhibiting mycelial growth in host tissues. An analysis of their concentrations could facilitate the selection of resistant varieties (Suchowilska 2008).

Recent research has greatly expanded our knowledge about the role of fusariotoxins in pathogenesis and their impact on humans and animals. The presence of the above metabolites in the human food chain is a serious worldwide concern. Trichothecenes are known to act by blocking protein biosynthesis in the cell. The phytotoxic effects of trichothecenes, mostly DON and NIV, have been demonstrated in a number of experiments, but their role in pathogenesis seems to be most significant. A relationship has been found between the FHB resistance of crops and the accumulation of fusariotoxins, particularly DON. Although DON does not necessarily cause typical FHB symptoms, the fungi that produce this toxin are usually highly virulent. The presence of DON also contributes to the spread of the fungus within the spike. NIV-producing chemotypes are less aggressive than DON-producing chemotypes (Miedaner and Reinbrecht, 2001). Many programs have been initiated to prevent and reduce mycotoxin contamination of wheat and other cereals. The co-occurrence of multiple mycotoxins can produce a combined synergistic effect. Resistant cereal lines inoculated with pathogens of the genus *Fusarium* contained lower concentrations of DON in the grain than susceptible cultivars (Miller et al., 1985, Mirocha et al., 1994). Harris et al. (1999) reported that the trichothecenes (especially DON) may act as virulence factors to enhance the spread of *F. graminearum* on maize. Cultivars, breeding strains and lines show a great variability in susceptibility to infection, which is reflected in various levels of yield reduction in different traits (Wiśniewska et al., 2004). Contemporary common wheat cultivars are all susceptible to *Fusarium* infection to a greater or lesser extent (Mesterhazy et al., 1999). The ingestion of food containing trichothecenes may induce numerous diseases termed mycotoxicoses whose symptoms resemble





**Figure 1:** FHB symptoms- field experiments with *Fusarium culmorum* infection of spring wheat



**Figure 2:** FHB symptoms, zoomed



**Figure 3:** *Fusarium culmorum* mycelium on agar plate

those of radiation sickness, including diarrhea, vomiting, nausea, leucocytosis (Alimentary Toxic Aleukia– ATA), bleeding in the gastrointestinal tract, and – in extreme cases – circulatory shock leading to death. Exposure to low doses of trichothecenes (encountered most frequently) results in weight loss, loss of appetite leading to malnutrition, neuroendocrine and immunological changes.

### References

- Harris LJ, Desjardins AE, Plattner RD, Nicholson P, Butler G, Young JC, Weston G, Proctor RH and Hohn TM. 1999. Possible role of trichothecene mycotoxins in virulence of *Fusarium graminearum* on maize. *Plant Dis* 83:954-960.
- Joffe A. 1986. *Fusarium* species. Their Biology and Toxicology. Red J. Willey & Sons Ltd,
- Mesterházy A, Bartok T, Mirocha CG and Komoroczy R. 1999. Nature of wheat resistance to *Fusarium* head blight and the role of deoxynivalenol for breeding. *Plant Breed* 118: 97–110.
- Mesterházy A. 2002. Role of deoxynivalenol in aggressiveness of *Fusarium graminearum* and *F. culmorum* and in resistance to *Fusarium* head blight. *Eur J Plant Pathol* 108: 675–684.
- Miller JD, Young JC and Sampson RD. 1985. Deoxynivalenol and *Fusarium* head blight resistance in spring cereals. *J Phytopathol* 113: 359–367.
- Mirocha CJ, Weiping Xie, Yichun Xu, Wilcoxson RD, Woodward RP, Etebarian RH and Behele G. 1994. Production of trichothecenes mycotoxins by *Fusarium graminearum* and *Fusarium culmorum* on barley and wheat. *Mycopathologia* 128: 19–23.
- Nelson P. E. 1993. Fumonisin, Mycotoxins Produced by *Fusarium* Species: Biology, Chemistry and Significance. *Phytopathology*, 31: 233-252
- Parry D.W., Jenkinson, P. and McLeod L. 1995. *Fusarium* ear blight (scab) in small grain cereals-a review. *Plant Path.* 44: 207-238.
- Suchowilska E., Wiwart M. and Grabowska G. 2008. Antifungal activity of methanol extracts from spikes of *Triticum spelta* and *Triticum aestivum* genotypes differing in their response to *Fusarium culmorum* inoculation. *Biologia* 63/4: 471-476.
- Ueno Y. 1997. Trichothecenes: overview adress. W: Mycotoxins in human and animal health. Ed. Rodricks J.V., Hesselstine C.W., Mehlman M., A., Patotox Publishers Inc. p.189-207.
- Wiśniewska H, Perkowski J and Kaczmarek Z. 2004. Scab Response and Deoxynivalenol Accumulation in Spring Wheat Kernels of Different Geographical Origins Following Inoculation with *Fusarium culmorum*. *J Phytopathol* 152: 613–621.

**SAVE  
55%**

## Discount on all Media designated for Fungi Detection and Cultivation

Media	Brand	Cat. No.	Organisms/remarks	Isolation	Cultiv.	Ident.	Selective	Recommended
Czapek Dox Agar	Fluka	70185	fungi and yeasts, defined C-source and inorg. N-source		x			APHA
Dichloran Rose bengal Agar (Base)	Fluka	17147	yeasts and moulds	x		x	x	
Peptone Yeast Extract Agar	Fluka	77196	yeasts and moulds	x			x	
Potato Glucose Agar	Fluka	70139	yeasts and moulds	x	x			APHA, US FDA
Potato Glucose Rose bengal Agar (Base)	Fluka	17204	For promoting ascospore production.	x	x		x	APHA, US FDA
Rose bengal Chloramphenicol Agar	Fluka	17211	yeasts and moulds	x		x	x	
Sabouraud 2% Glucose Agar	Fluka	84086	dermatophytes, fungi and yeasts	x	x			
Sabouraud 4% Glucose Agar	Fluka	84088	dermatophytes, fungi and yeasts	x	x			USP, EP
Sabouraud 4% Glucose Agar Plates	Fluka	40376	dermatophytes, fungi and yeasts	x	x			USP, EP
Selective Agar for pathogenic fungi	Fluka	84886	pathogenic fungi, in particular dermatophytes	x				
Universal Beer Agar	Fluka	17226	yeasts and organisms which are responsible for beer spoilage (wild yeast, bacteria, moulds)	x				
WL Nutrient Agar	Fluka	17222	yeasts and contaminants (yeasts, bacteria, moulds)	x				
WL Nutrient Broth	Fluka	W2261	yeasts and contaminants (yeasts, bacteria, moulds)	x				
Wort Agar	Fluka	70196	yeasts and moulds	x	x			
Yeast Malt Agar	Sigma	Y3127	yeasts, moulds and other aciduric microorganisms	x	x			
YGC Agar	Fluka	95765	yeasts and molds in dairy products.	x	x	x	x	ISO 6611:1992

When placing your order please quote **promotion code 981** or contact your local Sigma-Aldrich sales office. Offer valid until March 31st, 2010.

# Campylobacter

By Georgina Manning, Emerging Pathogen Group, School of Science and Technology, Nottingham Trent University, Nottingham, NG11 8NS, UK.... [georgina.manning@ntu.ac.uk](mailto:georgina.manning@ntu.ac.uk)

## Campylobacter is a major human pathogen worldwide.

The two species *Campylobacter jejuni* and *C. coli* are responsible for the majority of human food-borne infections in the UK. In 2008 there were 49,880 reported cases of Campylobacteriosis in England and Wales ([www.hpa.org.uk](http://www.hpa.org.uk)). Currently *C. jejuni* is known to account for approximately 90% of all reported cases. *C. coli*, accounting for 7-10% of the Campylobacteriosis cases, is still the fourth most common cause of gastrointestinal illness in England and Wales and therefore remains a significant health burden.

Human food poisoning is acquired through the consumption of contaminated food. Poultry is a major food source of infection and infection can occur through the consumption of undercooked poultry meat or through cross-contamination of other food products in the kitchen. In a recent report from the Food Standards Agency the prevalence of *Campylobacter* in chicken at retail was 65.2% ([www.food.gov.uk](http://www.food.gov.uk)). Other sources of infection include contaminated raw milk and water.

Symptoms of the disease can occur up to two weeks after consuming the contaminated food. They usually consist of a fever and general feeling of unwellness which is accompanied or followed by abdominal cramps and then watery,

sometimes bloody, diarrhoea (Blaser and Engberg 2008). The disease is usually self-limiting after 2-7 days; however, in some cases, particularly in the very young, elderly or immuno-compromised, antibiotic treatment may be required. In more severe cases, hospitalization may result.

The mechanisms utilised by *Campylobacter* to cause disease are not yet fully defined. However, it is thought that the organism first colonises the small intestine before travelling to the large intestine, where it causes a local inflammatory response (Black et al. 1988). The organism is motile by way of polar flagella and it can propel itself through the mucous layer on top of the epithelial cells. Adhesion to the surface of the epithelial cells followed by invasion into the cells is also thought to occur before the onset of diarrhoea. *Campylobacter* does produce a cytotoxic distending toxin, however the role of this toxin in disease is not fully characterized (Ketley 1997).

*C. jejuni* and *C. coli* are Gram negative, spiral, rod shaped cells (Figure 1). They are microaerophilic, requiring oxygen levels between 2 and 10% for growth. These two species are classed as thermophilic campylobacters as they have an optimum temperature for growth of 42 °C. *Campylobacter* can often be a difficult organism to isolate and grow in the laboratory. The current method used for detection of campylobacters from food-stuffs is the EN/ISO 10272-1 method. This method relies on homogenization of a known amount of material in Bolton Broth (BB) with selective supplement and 5% lysed horse blood either microaerobically or in a Schott bottle filled with Bolton enrichment broth (leaving a




Figure 1: Electron micrograph of *Campylobacter jejuni* (Gaynor et al. 2004)




Figure 2: *Campylobacter* colonies on mCCDA

gap at the top of no more than 2 cm) aerobically at 41.5 °C for 48 hours. However in a recent publication (Moran et al. 2009) a change to the ISO method was recommended, as it was found that aerobic incubation during enrichment was not very successful in isolation of the organism. Therefore, microaerobic enrichment was suggested as the preferred method, with incubation for 24 hours being sufficient. Following incubation, a sample of the enrichment broth is then plated onto modified charcoal cefoperazone deoxycholate agar (mCCDA, see Figure 2) with the addition of a selective supplement. Identification is usually done by colony morphology and speciation by PCR-based methods.

## References

1. Black, R.E., Levine, M.M., Clements, M.L., Hughes, T.P. and Blaser, M.J. (1988) Experimental *Campylobacter jejuni* infection in humans. *J. Infect. Dis.* 157, 472-479.
2. Blaser, M. J. and Engberg, J. (2008) Clinical aspects of *Campylobacter jejuni* and *C. coli* infections. In *Campylobacter* ed. Nachamkin, I., Szymanski, C.M. and Blaser, M.J. Washington DC: ASM Press. pp. 99-121.
3. Gaynor, E.C., Cawthraw, S., Manning, G., MacKichan, J.K., Falkow, S. and Newell, D.G. (2004) The genome-sequenced variant of *Campylobacter jejuni* NCTC 11168 and the original clonal clinical isolate differ markedly in colonization, gene expression, and virulence-associated phenotypes. *J. Bacteriol.* 186, 503-17.
4. Ketley, J.M. (1997) Pathogenesis of enteric infection by *Campylobacter*. *Microbiology* 143 ( Pt 1), 5-21.
5. Moran, L., Kelly, C. and Madden, R.H. (2009) Factors affecting the recovery of *Campylobacter* spp. from retail packs of raw, fresh chicken using ISO 10272-1:2006. *Lett. Appl. Microbiol.* 48, 628-632.



## Did you know...

what *Campylobacter* means?

*Campylobacter* means "twisted bacteria" because these bacteria are spiral-shaped. (see Figure 1)

# More to *Campylobacter* Species...

By Jvo Siegrist, Product Manager Microbiology... [ivo.siegrist@sial.com](mailto:ivo.siegrist@sial.com)

... some helpful information for detection and differentiation of *Campylobacter*.

## Media

*Campylobacter* are generally very fastidious microorganisms and grow only on complex media that has been amended with diverse essential amino acids and supplements, such as pyruvate,  $\alpha$ -ketoglutarate, hemin, formate and other essential ions. For selective isolation of *Campylobacter*, the growth media can be supplemented with antibiotics such as cefoperazone, vancomycin, trimethoprim, amphotericin, cycloheximide, rifampicin, cefsulodin and polymyxin B sulfate. Typical agars and broths used for the detection, identification, differentiation, enumeration and cultivation of *Campylobacter* are listed in Table 1.




Figure 1: *Campylobacter jejunii* colonies on CCDA agar

## Biochemical Differentiation

The metabolism of *Campylobacter* is chemoorganotrophic, with amino acids and intermediates of the citric acid cycle serving as energy sources; typical carbohydrates cannot be used. *Campylobacter* reduce nitrate to nitrite, obtaining oxygen for their metabolism by this pathway. These distinctive metabolic reactions can be used for the differentiation and identification of *Campylobacter* species (see Table 2).

Nonselective Media	Brand	Cat. No.
Columbia Agar	Fluka	27688
Tryptic Soy Agar	Fluka	22091
Tryptic Soy Agar (ready media in flasks)	Fluka	79872
Tryptic Soy Agar Plates (Diameter 55 mm)	Fluka	57994
Tryptic Soy Agar, Vegitone	Fluka	14432
Tryptic Soy Broth	Fluka	22092
Tryptic Soy Broth (ready media in flasks)	Fluka	43592
Tryptic Soy Broth No. 2	Fluka	51228
Tryptic Soy Broth, Vegitone	Fluka	41298
Nonselective Media for biochemical Differentiation	Brand	Cat. No.
Blood Agar Base No. 2	Fluka	B1676
Hippurate Broth	Fluka	53275
OF Test Nutrient Agar	Fluka	75315
Selective Media with Differential System	Brand	Cat. No.
Mac Conkey Agar No 1	Fluka	70143
<i>Campylobacter</i> Selective Media	Brand	Cat. No.
Blood Free <i>Campylobacter</i> Broth	Fluka	59751
Blood Free <i>Campylobacter</i> Selectivity Agar (mCCDA)	Fluka	B2426
Bolton Broth Base	Fluka	67454
Brucella Broth Base	Fluka	B3051
<i>Campylobacter</i> Selective Agar (Base)	Fluka	21378
Karmali <i>Campylobacter</i> Agar (Base)	Fluka	17152

Table 1: Typical *Campylobacter* media

When placing your order please quote promotion code 981 or contact your local Sigma-Aldrich sales office. Offer valid until March 31st, 2010.

**SAVE  
55%**

Species		<i>C. jejuni</i>		<i>C. coli</i>	<i>C. lari</i>	<i>C. fetus</i>	
		subspecies <i>jejuni</i>	subspecies <i>doylei</i>			subspecies <i>fetus</i>	subspecies <i>veneralis</i>
Growth at	25°C	no	No	no	no	yes	Yes
	37°C	yes	Yes	yes	yes	yes	yes
	42°C	yes	No	yes	yes	partial	no
	Oxidase	yes	Yes	yes	yes	yes	yes
	Catalase	yes	Yes	yes	yes	yes	yes
	Nitrate Reduction	yes	No	yes	yes	yes	yes
	Hippurate Hydrolysis	yes	Yes	no	no	no	no
	Indoxyl Acetate Hydrolysis	yes	Yes	yes	no	no	no
	H <sub>2</sub> S Production	no	no	yes	yes	no	
Antibiotic Sensitivity							
	Nalidixic acid	S	S	S	R	R	R
	Cephalotin	R	S	R	R	S	S
	Penicillin	R	R	R	R	R	S
	TTC*	R	R	R	S	S	S

Table 2: Table of differentiating characteristics of *Campylobacter* species and subspecies


(\* triphenyltetrazolium chloride)

<b>Campylobacter Test</b>	<b>Brand</b>	<b>Cat. No.</b>
Catalase Test	Fluka	88597
Gram staining kit	Fluka	77730
Hippurate Disks	Fluka	40405
Hydrogen Sulfide Test Strips	Fluka	06728
Indoxyl Strips	Fluka	04739
Nitrate Reagent A	Fluka	38497
Nitrate Reagent B	Fluka	39441
Nitrate Reagent Disks	Fluka	08086
Oxidase Test	Fluka	70439
Oxidase Strips	Fluka	40560
Oxidase Reagent acc. Gaby-Hadley A	Fluka	07345
Oxidase Reagent acc. Gaby-Hadley B	Fluka	07817
Oxidase Reagent acc. Gordon-McLeod	Fluka	18502

**Table 3:** Sigma-Aldrich tests for identification and differentiation of *Campylobacter*

### Detection of *Campylobacter* by rRNA

A further rapid and simple option to screen *Campylobacter* is the detection of specific rRNA. With oligonucleotide probes, specific sequences on the rRNA are detected and made visible with a sandwich hybridisation followed by a chromogenic reaction. It is a very easy method and does not need any PCR, is not sensitive to sample matrix and detects only living cells, as rRNA is destroyed within a few hours. No special expensive equipment is needed and the test is done within approximately 2 to 2.5 hours. The test format is based on a 96-well microplate with 12 strips of 8 wells so that small numbers of samples can be screened. A positive result is visible to the naked eye, however it is also possible by a standard microplate reader to quantify the number of cells at 450 nm. With a pre-enrichment step (44 - 48 hours in Bolton broth), a sensitivity of 1 CFU/25 g can be reached. The HybriScan®D *Campylobacter* Kit specifically detects the most relevant *Campylobacter* spp. (*C. jejuni*, *C. coli*, *C. lari*, *C. upsaliensis*). HybriScanD *Campylobacter* was compared and validated (acc. to EN ISO 16140:2003) with the cultivation based method according to § 64 LFGB (official method in Germany). Five different food categories were tested. The results of the validation lead to a relative accuracy of 95.2%, a relative specificity of 97.5%, and relative sensitivity of 93%, respectively.



**Figure 2:** Scheme of the test procedure



## Sigma-Aldrich Worldwide Offices

### Argentina

Free Tel: 0810 888 7446  
Tel: (+54) 11 4556 1472  
Fax: (+54) 11 4552 1698

### Australia

Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Austria

Tel: (+43) 1 605 81 10  
Fax: (+43) 1 605 81 20

### Belgium

Free Tel: 0800 14747  
Free Fax: 0800 14745  
Tel: (+32) 3 899 13 01  
Fax: (+32) 3 899 13 11

### Brazil

Free Tel: 0800 701 7425  
Tel: (+55) 11 3732 3100  
Fax: (+55) 11 5522 9895

### Canada

Free Tel: 1800 565 1400  
Free Fax: 1800 265 3858  
Tel: (+1) 905 829 9500  
Fax: (+1) 905 829 9292

### Chile

Tel: (+56) 2 495 7395  
Fax: (+56) 2 495 7396

### China

Free Tel: 800 819 3336  
Tel: (+86) 21 6141 5566  
Fax: (+86) 21 6141 5567

### Czech Republic

Tel: (+420) 246 003 200  
Fax: (+420) 246 003 291

### Denmark

Tel: (+45) 43 56 59 00  
Fax: (+45) 43 56 59 05

### Finland

Tel: (+358) 9 350 9250  
Fax: (+358) 9 350 92555

### France

Free Tel: 0800 211 408  
Free Fax: 0800 031 052  
Tel: (+33) 474 82 28 88  
Fax: (+33) 474 95 68 08

### Germany

Free Tel: 0800 51 55 000  
Free Fax: 0800 64 90 000  
Tel: (+49) 89 6513 0  
Fax: (+49) 89 6513 1160

### Hungary

Ingyenes telefonszám: 06 80 355 355  
Ingyenes fax szám: 06 80 344 344  
Tel: (+36) 1 235 9063  
Fax: (+36) 1 269 6470

### India

#### Telephone

Bangalore: (+91) 80 6621 9400  
New Delhi: (+91) 11 4358 8000  
Mumbai: (+91) 22 2570 2364  
Hyderabad: (+91) 40 4015 5488  
Kolkata: (+91) 33 4013 8003  
Bangalore: (+91) 80 6621 9650  
New Delhi: (+91) 11 4358 8001  
Mumbai: (+91) 22 4087 2364  
Hyderabad: (+91) 40 4015 5488  
Kolkata: (+91) 33 4013 8000

### Ireland

Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: (+353) 402 20370  
Fax: (+353) 402 20375

### Israel

Free Tel: 1 800 70 2222  
Tel: (+972) 8 948 4100  
Fax: (+972) 8 948 4200

### Italy

Tel: (+39) 02 3341 7310  
Fax: (+39) 02 3801 0737

### Japan

Tel: (+81) 3 5796 7300  
Fax: (+81) 3 5796 7315

### Korea

Free Tel: (+82) 80 023 7111  
Free Fax: (+82) 80 023 8111  
Tel: (+82) 31 329 9000  
Fax: (+82) 31 329 9090

### Malaysia

Tel: (+60) 3 5635 3321  
Fax: (+60) 3 5635 4116

### Mexico

Free Tel: 01 800 007 5300  
Free Fax: 01 800 712 9920  
Tel: (+52) 722 276 1600  
Fax: (+52) 722 276 1601

### The Netherlands

Free Tel: 0800 022 9088  
Free Fax: 0800 022 9089  
Tel: (+31) 78 620 5411  
Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Norway

Tel: (+47) 23 17 60 00  
Fax: (+47) 23 17 60 10

### Poland

Tel: (+48) 61 829 01 00  
Fax: (+48) 61 829 01 20

### Portugal

Free Tel: 800 202 180  
Free Fax: 800 202 178  
Tel: (+351) 21 924 2555  
Fax: (+351) 21 924 2610

### Russia

Tel: (+7) 495 621 5828  
Fax: (+7) 495 621 5923

### Singapore

Tel: (+65) 6779 1200  
Fax: (+65) 6779 1822

### Slovakia

Tel: (+421) 255 571 562  
Fax: (+421) 255 571 564

### South Africa

Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: (+27) 11 979 1188  
Fax: (+27) 11 979 1119

### Spain

Free Tel: 900 101 376  
Free Fax: 900 102 028  
Tel: (+34) 91 661 99 77  
Fax: (+34) 91 661 96 42

### Sweden

Tel: (+46) 8 742 4200  
Fax: (+46) 8 742 4243

### Switzerland

Free Tel: 0800 80 00 80  
Free Fax: 0800 80 00 81  
Tel: (+41) 81 755 2828  
Fax: (+41) 81 755 2815

### United Kingdom

Free Tel: 0800 717 181  
Free Fax: 0800 378 785  
Tel: (+44) 1747 833 000  
Fax: (+44) 1747 833 313

### United States

Toll-Free: 800 325 3010  
Toll-Free Fax: 800 325 5052  
Tel: (+1) 314 771 5765  
Fax: (+1) 314 771 5757

### Vietnam

Tel: (+84) 3516 2810  
Fax: (+84) 6258 4238

### Internet

[sigma-aldrich.com](http://sigma-aldrich.com)



**Mixed Sources**  
Product group from well managed  
forests, controlled sources and  
recycled wood or fiber  
[www.fsc.org](http://www.fsc.org) Cert no. SGS-COC-XXXXXX  
© 1996 Forest Stewardship Council



*Accelerating Customers'  
Success through Innovation and  
Leadership in Life Science,  
High Technology and Service*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052  
Technical Service (800) 325-5832 • [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
Development/Custom Manufacturing Inquiries **SAFC** (800) 244-1173  
Safety-related Information [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

**World Headquarters**  
3050 Spruce St.  
St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

©2010 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

MEQ  
T410026