

| Organisms (ATCC)                          | Growth | Colony Appearance  |
|---|--------|--------------------|
| <i>Bifidobacterium longum</i> (15707)     | +++    | red-brown (maroon) |
| <i>Bifidobacterium infantis</i> (15697)   | +++    | red-brown (maroon) |
| <i>Streptococcus thermophilus</i> (14486) | -      | -                  |
| <i>Lactobacillus acidophilus</i> (314)    | -      | -                  |
| <i>Lactobacillus bulgaricus</i> (11842)   | -      | -                  |

**Table 1.** *Bifidobacterium* sp. cultural characteristics on BSM Agar

| Name                            | Brand | Cat. No.     | Pack Size       |
|---------------------------------|-------|--------------|-----------------|
| BSM Agar                        | Fluka | <b>88517</b> | 500 g           |
| BSM Broth                       | Fluka | <b>90273</b> | 500 g           |
| BSM Supplement                  | Fluka | <b>83055</b> | 5 g             |
| Wilkins Chalgren Anaerobic Agar | Fluka | <b>W1761</b> | 500 g           |
| Mupirocin                       | Sigma | <b>M7694</b> | 50 mg<br>100 mg |
| Lithium Mupirocin               | Sigma | <b>07188</b> | 1 g<br>100 mg   |

**Table 2.** Selective medium and supplements for *Bifidobacterium*

*Bifidobacterium* grow very well on this medium, while *Lactobacillus* and *Streptococcus* strains are inhibited. *Bifidobacterium* colonies grow within 24-48 hours (occasionally up to three days because of the highly selective conditions). The *Bifidobacterium* colonies are purple-brown and therefore are easy to differentiate from other organisms.

In a Swiss governmental evaluation study for the enumeration of bifidobacteria in sour milk products, the traditional method was compared to Wilkins-Chalgren Agar with 100 mg/L mupirocine and BSM Agar. The traditional method produced statistically significant differences, while Wilkins-Chalgren Agar and BSM Agar showed similar results without any significant variances. The study concluded, "On the BSM Agar, the bifidobacteria forms purple-brown colonies which made the enumeration easy" [3].

### **Lactobacillus Species**

Lactobacilli are rod-shaped, Gram-positive, fermentative, facultative anaerobic or microaerophilic organotrophs. Normally they form straight rods, but under certain conditions spiral or coccobacillary forms have been observed. In most cases, they form chains of varying length. Lactobacilli belong to the lactic acid bacteria and comprise the major part of this group. As their name implies, they produce lactic acid and derive energy from the fermentation of lactose, glucose and other sugars to lactate via homofermentative metabolism. About 85-90% of the sugar utilized in the fermentative process is converted to lactic acid. However, there are some heterofermentative lactobacilli that produce alcohol in addition to lactic acid from sugars. This acid-producing mechanism inhibits growth of other organisms and favors the growth of lactobacilli that thrive in low pH environments. ATP is generated during the process by non-oxidative substrate-level phosphorylation.

Some strains of lactobacilli were shown to produce, like bifidobacteria, a bacteriocin-like substance and are able to inhibit a broad range of pathogens [8]. Lactobacilli also produce adhesins (proteins), which perform a vital role in recognizing specific host components (extracellular matrix) important for the bacterial adhesion and colonization at host surfaces, as well as in bacterial interaction with physiological and immunological processes [8].

In the last several years, several new *Lactobacillus* species have been introduced as probiotics, including *Lactobacillus rhamnosus*, *Lactobacillus casei*, and *Lactobacillus johnsonii*. As of this time, probiotics have not been used in the pharmaceutical industry due to the many open questions that remain to be answered [7].



**Figure 1.** Yogurt sample cultured on BSM Agar. Bifidobacteria appears as purple-brown colonies.

### **Medium for Lactobacilli**

Since lactobacilli prefer acidic conditions, natural extracts and juices from tomatoes and oranges, as well as other single metabolic acids (e.g. malic acid), are often used as media ingredients. Casein and yeast extract provide rich amino acid sources, and the maltose is used as a carbohydrate source for lactobacilli, which cannot utilize glucose as fermentable sugar. Fructose is the carbohydrate source of *Lactobacillus fructivorans*. Polysorbate, sorbitan mono-oleate and other related compounds act as a source of fatty acids and stimulate the lactic acid bacteria.

Today, it is standard practice to differentiate lactobacilli based on their phenotype using selective media. Classical phenotypic tests for identification of lactobacilli are based on physiological characteristics, like motility, growth temperature, respiratory type, and growth in sodium chloride, as well as on diverse biochemical characteristics, such as fermentation type, metabolism of carbohydrate substrates, production of lactic acid isomers, coagulation of milk, and presence of specific enzymes like arginine dihydrolase. In Bergey's Manual, *Lactobacillus* is described as a Gram-positive rod, non-spore forming, acid fast negative and catalase negative. The colony morphology on certain media is taken for the presumptive identification.



**Figure 2.** 3D rendering of isolated *Lactobacillus*



| Cat. No.  | Name  | Description   |
|---|---|---|
| <b>Media for Differentiation</b>                        |   |   |
| <b>17153</b>  | LS Differential Agar<br>Antibiotic free skim milk powder (Fluka 70166)<br>TTC Solution (Fluka 17779)    | For the maximum growth and differentiation of lactobacilli and streptococci on the basis of colonial morphology, T.T.C. reduction and casein reaction.  |
| <b>17158</b>  | Litmus Milk   | For maintenance of lactobacilli and for determining the action of bacteria on milk.   |
| <b>17222</b>  | WL Nutrient Agar  | For the examination of materials encountered in brewing and for industrial fermentations containing mixed flora of yeasts and bacteria.   |
| <b>17215</b>  | WL Differential Agar  | For selective isolation and enumeration of bacteria encountered in breweries and industrial fermentations.  |
| <b>Nonselective Media</b>                               |   |   |
| <b>Y3127</b>  | Yeast Malt Agar   | Used for the isolation and cultivation of yeasts, molds and other aciduric microorganisms.  |
| <b>17216</b>  | Tomato Juice Agar   | For the cultivation and enumeration of lactobacilli.  |
| <b>80957</b>  | Plate Count Skim Milk Agar  | For the enumeration of bacteria in milk and dairy products.   |
| <b>17218</b>  | Tomato Juice Broth  | For cultivation of yeasts and other aciduric microorganisms.  |
| <b>T2188</b>  | Tryptone Glucose Yeast Extract Agar   | Recommended for enumeration of bacteria in water, air, milk and dairy products.   |
| <b>17123</b>  | Elliker Broth   | For culturing streptococci and lactobacilli of importance in the dairy industry.  |
| <b>Selective Media (by low pH and selective agents)</b> |   |   |
| <b>R1148</b>  | Rogosa SL Agar  | Used as a selective medium for cultivation of lactobacilli.   |
| <b>85515</b>  | Sorbic Acid Agar (Base)   | For the isolation and differentiation of lactobacilli from food, faeces etc. according to Reuter.   |
| <b>17154</b>  | Lactobacillus Bulgaricus Agar (Base)  | Used with acetate buffer for isolation and identification of <i>Lactobacillus bulgaricus</i> .  |
| <b>69966</b>  | MRS Broth   | For the enrichment and isolation of all species of lactobacilli from all types of material.   |
| <b>41782</b>  | MRS Agar, Vegitone  | This MRS Agar is free of animal derived material. It is recommended for the isolation and cultivation of <i>Lactobacillus</i> species.  |
| <b>69964</b>  | MRS Agar  | For the enrichment, cultivation and isolation of all species of <i>Lactobacillus</i> from all types of material according to De Man, Rogosa and Sharpe.   |
| <b>38944</b>  | MRS Broth, modified, Vegitone   | This MRS Broth contains plant peptone instead of animal peptone. It is recommended for the isolation and cultivation of <i>Lactobacillus</i> species.   |
| <b>30912</b>  | MRS Agar, original acc.<br>DeMan-Rogosa-Sharpe  | For the enrichment, cultivation and isolation of all species of <i>Lactobacillus</i> from all types of materials. Recommended by the "Schweizerisches Lebensmittelbuch" 5th ed., chapter 56A.             |
| <b>75405</b>  | Orange-serum Agar   | For the isolation, cultivation and enumeration of acid-tolerant spoilage microorganisms in fruit juice and fruit juice concentrates, in particular from citrus fruit, according to Hays, Troy and Beisel. |
| <b>17226</b>  | Universal Beer Agar   | For culturing microorganisms of significance in the brewing industry.   |
| <b>W2261</b>  | WL Nutrient Broth   | Recommended for the cultivation of bacteria encountered in breweries and industrial fermentations.  |
| <b>02538</b>  | Raka Ray Agar<br>Lactic Acid Supplement, modified (Fluka 14121)<br>Cycloheximide solution (Fluka 18079) | A medium for selective isolation of lactic acid bacteria from beer and brewing processes.   |
| <b>83920*</b>   | Rogosa Agar   | Selective agar for the isolation and enumeration of lactobacilli belonging to the buccal and intestinal flora, as well as from food such as milk and meat.  |
| <b>64198</b>  | NBB Agar Base, modified   | Selective medium used for the detection of contaminating/spoilage microorganisms in brewery.  |
| <b>50725</b>  | NBB Broth Base, modified  | Selective medium used for the detection of contaminating/spoilage microorganisms in brewery.  |

**Table 3. Media for lactobacilli detection, differentiation, isolation, enumeration and cultivation**

\*Available in Europe only.



## Modern Methods to Detect Lactobacilli

As a modern alternative, molecular biology-based methods, like PCR, can be consulted. However, they are often quite expensive. Under our Fluka Brand, Sigma-Aldrich provides a revolutionary molecular biology method that is rapid, easy and cost effective. Based on the detection of rRNA, this method completely avoids

the need for PCR amplification. The sandwich hybridization test, called HybriScan, is performed on a microtiter plate. The range of lactobacilli detected by HybriScan tests is listed in **Table 4**.

More information about the test and the technical principles behind it are found on our web site [sigma-aldrich.com/hybriscan](http://sigma-aldrich.com/hybriscan)

| Cat. No.     | Brand | Description                                  | Specifity   | Reactions |
|--------------|-------|--|---|-----------|
| <b>62533</b> | Fluka | HybriScan®D Beer                             | Detects all beer spoiling bacteria of the genera <i>Lactobacillus</i> , <i>Pediococcus</i> , <i>Pectinatus</i> and <i>Megasphaera</i>   | 96        |
| <b>68301</b> | Fluka | HybriScan®D Drinks                           | Detects yeast of the genera <i>Saccharomyces</i> , <i>Zygosacchomyces</i> , <i>Brettanomyces</i> , <i>Torulasporea</i> , <i>Pichia</i> , <i>Candida</i> and bacteria of the genera <i>Lactobacillus</i> , <i>Acetobacteraceae</i> and <i>Alicyclobacillus</i> | 96        |
| <b>59744</b> | Fluka | HybriScan®D Lactobac                         | Detects bacteria of the genera <i>Lactobacillus</i> and <i>Pediococcus</i> in fruit juices and non-alcoholic beverage.  | 96        |
| <b>75724</b> | Fluka | HybriScan®I<br><i>Lactobacillus brevis</i>   | Identification of <i>Lactobacillus brevis</i>   | 48        |
| <b>80065</b> | Fluka | HybriScan®I<br><i>Lactobacillus buchneri</i> | Identification of <i>Lactobacillus buchneri</i>   | 48        |
| <b>86827</b> | Fluka | HybriScan®I<br><i>Lactobacillus lindneri</i> | Identification of <i>Lactobacillus lindneri</i>   | 48        |

**Table 4.** HybriScan products (HybriScan®D = detection kit; HybriScan®I = identification kit)

## References

1. F. Ibrahim, S. Ruvio, L. Granlund, S. Salminen, M. Viitanen, A.C. Ouwehand, Probiotics and immunosenescence: cheese as a carrier, *FEMS Immunology & Medical Microbiology*, 59, 1, p53-59 (2010)
2. Rada V. and Koc J.: The use of mupirocin for selective enumeration of bifidobacteria in fermented milk products. *Milchwissenschaft*. 55, 65-67 (2000)
3. Bundesamt für Gesundheit: Schweizerisches Lebensmittelbuch (SLMB): Kapitel 56, «Mikrobiologie» Neuauflage 2000, Stand 2004
4. J. Nowroozi, M. Mirzaii, M. Norouzi, Study of *Lactobacillus* as Probiotic Bacteria, *Iranian J. Publ. Health*, Vol. 33, No. 2, pp.1 -7 (2004)
5. R.B. Vaughan, The romantic rationalist: A study of Elie Metchnikoff, *Medical History* 9: 201-15 (1965)
6. H. Tissier, Recherchers sur la flora intestinale normale et pathologique du nourrisson, Thesis, University of Paris (1900)
7. G.W. Tannock, Probiotics: time for a dose of realism, *Curr. Issues Intest. Microbiol.*, 4(2), 33-42 (2003)
8. Å. Ljungh, T. Wadström, *Lactobacillus Molecular Biology: From Genomics to Probiotics*, Caister Academic Press (2009)
9. K. Lindfors, et al., Live probiotic *Bifidobacterium lactis* bacteria inhibit the toxic effects induced by wheat gliadin in epithelial cell culture, *Clin. Exp. Immunol.*, 152(3), 552-8 (2008)
10. J.F. Kim, et al., Genome sequence of the probiotic bacterium *Bifidobacterium animalis* subsp. *lactis* AD011, *J. Bacteriol.*, 191(2), 678-9 (2009)

## Genus *Yersinia*

Jvo Siegrist, Product Manager Microbiology  
[jvo.siegrist@sial.com](mailto:jvo.siegrist@sial.com)

### Detection, identification, differentiation and cultivation of *Yersinia* species

*Yersinia* is a rod shaped and facultative anaerobic Gram-negative bacteria belonging to the family of Enterobacteriaceae. It has a fermentative metabolism, is oxidase-negative, mannitol-positive, glucose-positive and lactose negative. On the Bismuth sulfite Agar (Fluka 95388) *Yersinia* can be differentiated from *Salmonella* because it is not able to produce hydrogen sulphide. It is a psychrophilic organism, surviving and proliferating at low temperatures of 0-4 °C (e.g., on food products in a refrigerator). Some *Yersinia* species are also relatively highly heat resistant. However, they can be quite easily inactivated by oxidizing agents such as hydrogen peroxide and potassium permanganate.

Pigs, rodents, rabbits, sheep, cattle, horses, dogs, and cats are the natural sources of *Yersinia*. At the moment, most human illness cases caused by *Yersinia* originate from *Y. enterocolitica*. This organism is the cause of yersiniosis, an infectious disease with symptoms like fever, abdominal pain, and diarrhea. Other clinically important species of this genus are *Y. pseudotuberculosis* (symptoms similar to *Y. enterocolitica* except in most cases no diarrhea is seen) and *Y. pestis* (organism responsible for the bubonic plague). Most infections are acquired through contaminated food, like raw or undercooked pork products, seafood, vegetables, unpasteurized milk or untreated water. However, infections may also occur after contact with infected animals or faeces, or through transmission by fleas.

The bacteria received its name from A. E. J. Yersin, a Swiss microbiologist, who discovered the *Yersinia pestis* bacterium in 1894 in Hong Kong.



| Species                      | Differentiation Properties   |
|------------------------------|--|
| <i>Y. enterocolitica</i>     | catalase-positive, most strains are ornithine-positive, motile at 22-26 °C, non-motile at 37 °C, urea-positive, sorbitol- and cellobiose- positive, most strains are sucrose-positive  |
| <i>Y. pseudotuberculosis</i> | motile at 22-26 °C, non-motile at 37 °C, urea-positive, rhamnose-positive, esculin-positive, ornithine decarboxylase negative, indole negative   |
| <i>Y. pestis</i>             | produces two antiphagocytic components (antiphagocytic slime), catalase positive non-motile at 22-26 °C and 37 °C, esculin-positive, indole negative, ornithine decarboxylase negative |

Table 1. Some differentiation of the important *Yersinia* species C

| Cat. No.   | Description  |
|--|--|
| <b>Nonselective Broths</b>   |  |
| 53286  | Brain Heart Broth                                      |
| 40893  | Peptone Water, phosphate-buffered, Vegetone            |
| 77187  | Peptone Water, phosphate-buffered                      |
| <b>Selective Enrichment Broths</b>   |  |
| 17156  | ITC Broth (Base)                                       |
| 69965  | Mossel Broth   |
| 17192  | Peptone Sorbitol Bile Broth                            |
| <b>Basal Media for Carbohydrates Utilization</b>   |  |
| A0715  | Andrade Peptone Water                                  |
| 28943  | Andrade Peptone Water, Vegetone                        |
| <b>Medium for Autoagglutination</b>  |  |
| 39484  | Methyl Red Voges Proskauer Broth                       |
| <b>Nonselective Agars for Differentiation &amp; Confirmation</b>                             |  |
| 27048*   | Christensen's Urea Agar                                |
| 60787  | Kligler Agar   |
| M1428  | Motility Indole Ornithine Medium                       |
| 75315  | OF Test Nutrient Agar                                  |
| O5386  | Ornithine Decarboxylase Broth                          |
| 51463  | Urea Broth   |
| <b>Selective Agars for Detection and Isolation</b>   |  |
| 95760*   | Yersinia Selective Agar (CIN Agar)                     |
| <b>Selective Agars with Differential System for Differentiation, Detection and Isolation</b> |  |
| 95388  | Bismuth sulfite Agar                                   |
| 70135  | DCLS Agar  |
| 90035  | DCLS Agar No. 2  |
| 17213  | Violet Red Bile Glucose Agar without Lactose           |
| 53605  | Violet Red Bile Glucose Agar without Lactose, Vegetone |
| 95273  | VRB MUG Agar   |

Table 2. Media for *Yersinia*

\*Available in Europe only.

## References

1. K.J. Ryan, C.G. Ray (editors), Sherris Medical Microbiology, 4th ed., McGraw Hill, 368–70 (2004)
2. F.M. Collins, Pasteurella, Yersinia, and Francisella. In: Barron's Medical Microbiology (S. Barron et al, eds.), 4th ed., University of Texas Medical Branch (1996)
3. F. Bichai, P. Payment, B. Barbeau, Protection of waterborne pathogens by higher organisms in drinking water: a review, Can. J. Microbiol. 54 (7), 509-524 (2008)
4. Bergey's Manual of Systematic Bacteriology, 1994, 9th ed., Williams and Wilkins, Baltimore

More details about the media or tests can be found on our website [sigma-aldrich.com/microbiology](http://sigma-aldrich.com/microbiology).




Figure 2. Bacteriophage ØR1-37 attached to *Y. enterocolitica* YeO3-R1 cell (length of cell = 1,7 µm; Source Mikael Skurnik, University of Turku, Finland)

| Cat. No. | Supplements                   |
|----------|-------------------------------|
| 38266    | Ferrioxamine E                |
| 75258    | Yersinia selective Supplement |
| 17778    | Ticarcillin Supplement        |
| 17777    | Potassium chlorate Supplement |

Table 3. Selective and growth supplement for *Yersinia*

| Cat. No. | Test for <i>Yersinia</i> Diagnostics |
|----------|--------------------------------------|
| 80507    | Bile Esculin Disks                   |
| 88597    | Catalase Test                        |
| 56481    | Cellobiose Disks                     |
| 63367    | Dextrose (Glucose) Disks             |
| 05686    | DMACA Indole Disks                   |
| 49825    | DMACA Reagent                        |
| 67309    | Kovac's Reagent for Indoles          |
| 78719    | Kovac's Reagent Strips               |
| 28816    | Lactose Disks                        |
| 94438    | Mannitol Disks                       |
| 08714    | Methyl Red Solution                  |
| 07345    | Oxidase Reagent acc. Gaby-Hadley A   |
| 07817    | Oxidase Reagent acc. Gaby-Hadley B   |
| 18502    | Oxidase Reagent acc. Gordon-McLeod   |
| 40560    | Oxidase Strips                       |
| 70439    | Oxidase Test                         |
| 93999    | Rhamnose Disks                       |
| 92971    | Salicin Disks                        |
| 93998    | Sorbitol Disks                       |
| 94309    | Sucrose Disks                        |

Table 4. Test for identification and differentiation of *Yersinia*