



*Clostridium perfringens* — an indicator

Surface water a source of *Clostridium perfringens* (source *C. perfringens* FISH image: BioVisible BV, The Netherlands) Clostridium perfringens — an indicator......p2

*Burkholderia cepacia* complex: multilingual bacteria......p6

SIGMA-ALDRICH®

# Clostridium perfringes — an Indicator

By Prof. Mohammad Manafi, Medical University of Vienna mohammad.manafi@meduniwien.ac.at

Jvo Siegrist, Sigma-Aldrich ivo.siegrist@sial.com

# *Clostridium perfringens* is an effective indicator of faecal contamination. Unlike *E. coli, C. perfringens* can withstand processing treatments and unfavorable growth conditions.

*Clostridium perfringens* is found in undercooked or improperly sterilized canned foods (germination of endospores) and in water (surface water). The natural contamination source is human and animal faeces transmitted into food products primarily by water. *C. perfringens* produces an extensive range of invasins and exotoxins. The enterotoxins cause the undesirable, predominantly meat-associated, food poisoning, as well as wound and surgical infections that lead to gas gangrene.

*C. perfringens* plays a subsidiary role in water examination [6]. *Clostridia* are spore builders and are resistant to heating, chlorination and other stress factors. In contrast to vegetative cells like coliforms (*E. coli*, enterococci), which are less resistant, *C. perfringens* has the advantage of a robust capability to survive [8]. Therefore, while faecal contamination is detected primarily by coliforms as an indicator, which could disappear after a processing step, *C. perfringens* remains present. The organism is not a hazard in water; rather, it is problematic when the water comes in contact with food.

In consideration of the aforementioned facts, it is obvious that detection and identification of *C. perfringens* is an important step toward the control and eradication of this potent pathogen. *C. perfringens* is an anaerobic, Gram-positive, spore-forming rod-shaped (see Figure 7, pg 4) bacterium [3, 5]. Some characteristic enzymes of *C. perfringens* are hemolysins ( $\beta$ -hemolysis), lecithinase, extracellular proteases, lipases (phospholipase-C), phosphatase, collagenase, hyaluronidase, saccharolytic enzymes and enzymes to reduce sulphite to sulphide. These enzymes are also used as detection and differentiation targets. It is also notable that *C. perfringens* is a non-motile bacterium, and it is the most important of the sulphite reducing clostridia [5].

Also, *C. perfringens* normally grows at 44 °C, whereas some other clostridia are inhibited at this temperature. This property is used in ISO methods to give the medium more selectivity [4].

Early detection of *Clostridium* in food is important to control outbreaks. To facilitate detection, Fluka has introduced a new chromogenic media, CP *ChromoSelect* Agar, for enumeration and differentiation of *Clostridium sp.*, in particular *Clostridium perfringens*, in aqueous samples (Figure 1).

For the standard detection of *C. perfringens*, mCP and TSC agar have been recommended [6, 7]. However, there are problems associated with each of these media. CP *ChromoSelect* Agar is more reliable and easier to handle than m-CP and TSC agars. The color does not diffuse in the agar and confirmation is not required since the green coloration is specific for *C. perfringens* [9]. CP *ChromoSelect* Agar avoids the disadvantages of m-CP agar, including the presence of ammonia, which prevents subculturing the *C. perfringens* colonies, the overly-selective nature of m-CP agar, and the evanescence of the red color of colonies after the addition of ammonia, which makes further confirmation impossible. In addition, the recovery of *C. perfringens* was rejected by ISO in favor of methods based on TSC agar [1, 8]. CP *ChromoSelect* Agar also eliminates the excessive and variable blackening of the peripheral colonies encountered with TSC agar, which makes colony counting at lower dilutions difficult and leads to false positives (**Figure 2**). It is also more reliable at high bacteria counts, where the TSC agar can produce false negatives because of an interference with the other enzymatic mechanisms from acid production and oxygen contact (**Figure 3**). TSC detects all

Figure 1: Drinking water sample cultured on CP ChromoSelect Agar. C. perfringens appears as distinct green colonies.





Figure 2: Comparison of m-CP, TSC and CP ChromoSelect Agar (from left to right)



Figure 3: *C. perfringens* ATCC 10873 on CP *ChromoSelect* Agar (above) and TSC agar (below). Note the false negatives on the TSC agar.



**Figure 4:** TSC agar + fluorogenic substrate = TSCF agar





Figure 5: C. perfringens on mCP Agar before (above) and after (below) ammonia hydroxide vapors





sulphite-reducing clostridia, however, and not only *C. perfringens*. TSC agar can be improved with the addition of MUP, a fluorogenic substrate which helps to confirm the presence of phosphatase. This modification is designated as TSCF (Figure 4).

In the present study, three media types (mCP [Figure 5], TSCF and CP *ChromoSelect* Agar) were evaluated for recovery of *C. perfringens* in various surface water samples. Using a membrane filtration technique on 139 water samples, 131 samples (94.2%) were found to be presumably positive for *C. perfringens* in at least one of the culture media. Green colored colonies on CP *ChromoSelect* agar (Figure 1) were counted as presumably *C. perfringens* isolates.

Out of 483 green colonies on CP *ChromoSelect* Agar, 96.9 % (465 strains, indole negative) were identified as *C. perfringens*, 15 strains (3.1%) were indole positive and were identified as *C. sordelli, C. bifermentans* or *C. tetani*. Only 3 strains (0.6 %) gave false positive results and were identified as *C. fallax, C. botulinum*, and *C. tertium* (Table 1). Variance analysis of the data obtained shows no statistically significant differences in the counts obtained between media employed in this work (Figure 6 and Table 2).

 Table 1: All green colonies isolated from CP ChromoSelect Agar and identified with API system (n=483)

strains	Indole reaction	n
C. perfringens	-	465 (96.3%)
C. tertium	-	1 ( 0.2%)
C. botulinum	-	1 ( 0.2%)
C. fallax	-	1 ( 0.2%)
C. bifermentans	+	2 ( 0.4%)
C. sordelli	+	12 ( 2.5%)
C. tetani	+	1 ( 0.2%)

**Figure 6:** Comparison between TSCF agar and the mCP and CPC (CP *ChromoSelect* Agar) media for enumerating strains of *C. perfringens* in water samples.




Table 2: Analysis of variance of C. perfringens spore counts obtained on all media

Effect	Sum of squares	df	Mean square	F ratio	p value
Media	0.392	2	0.196	0.190	0.827
Residual	235.578	228	1.033		



3

In general, the identification of typical and atypical colonies isolated from all media demonstrated that CP *ChromoSelect* Agar was the most useful medium for *C. perfringens* recovery in water samples [2].

In addition to its advantages over m-CP and TSC agars, CP *ChromoSelect* Agar is an ideal growth media. It contains only vegetable peptones and, together with yeast extract, it is an excellent source of nitrogen, carbon, amino acids and vitamin B complex. Sucrose acts as the fermentable carbohydrate and reducing agents lower the redox potential of the media. Diverse salts provide the required ions for enzymatic reactions. Buffering agents stabilize the pH within the ideal growth range. Inhibitors D-cycloserine and polymyxin B give the medium its selectivity, while further selectivity is achieved by incubation under anaerobic conditions at 44 °C. Various promoters and substrates protect injured cells to improve recovery rate and enhance growth. The chromogenic enzyme substrates in the CP *ChromoSelect* Agar provide the differentiation, particularly for *C. perfringens* (Table 3). A negative indol reaction (Kovac's Reagent) confirms the presence of *C. perfringens*.

Table 3: Clostridium sp. cultural characteristics in CP ChromoSelect Agar

Organisms (ATCC)	Growth	Colony Appearance
Clostridium perfringens (13124)	+++	Green
Clostridium bifermentans (638)	+++*	Dark blue with violet halo
Clostridium sporogenes (8534)	-	-
Clostridium sordelli (9714)	++	Dark green with halo (change to red with Kovac's Reagent)
Enterococcus faecalis (29212)	++	Violet
Escherichia coli (25922)	-	-
Pseudomonas aeruginosa (27853)		Colorless
Staphylococcus aureus (25923)	-	-
Bacillus subtilis (6051)	-	-
Salmonella typhimurium (DSM 554)	++	Violet

## Did you know...

*Clostridium perfringens* has one of the highest reproduction rates in the bacteria kingdom. Under optimal conditions, the cell doubles every 10 minutes!




Figure 7: Scanning electron micrograph of C. perfringens grown on a silicon wafer (source: S. Melville, Department of Biological Sciences, Virginia Tech University)

\* Growth at 40 °C, but no growth at 44 °C

#### Table 4: Typical Media for Clostridium perfringens

Nonselective Media	Cat. No.	Description
AEA Sporulation Broth (Base), modified	17170	For early sporulation of Clostridium perfringens from foods.
Supplements: 4 g/l raffinose (Fluka 83400) + 930 mg/l sodium		
carbonate (Fluka 71351) + 42.6 mg/l cobalt chloride		
(Sigma-Aldrich 60818) + 20 mg/l sodium ascorbate (Sigma 11140)		
Brewer thioglycollate medium	B2551	Used for testing the sterility of biological products.
Clostridial Nutrient Medium*	27546	For the cultivation and enumeration of clostridia and other anaerobes as well
		as facultative microorganisms in food, clinical and other material acc. to Hirsch & Grinsted.
Cooked Meat Broth	60865	For the primary cultivation of aerobic, microaerophilic and anaerobic bacteria
		from clinical specimens.
Reinforced Clostridial Agar	91365	Used for the cultivation and enumeration of Clostridia.
Nonselective Differential Media	Cat. No.	Description
Clostridial Differential Broth*	27544	For the MPN enumeration of all clostridia in food and other material acc. to Gibbs and Freame.
Gelatin Iron Agar	G0289	Used for detecting gelatin liquefaction and hydrogen sulphide production.
Lactose Gelatin Medium (Base)*	61348	For the detection of lactose and gelatine metabolizing microorganisms
Supplement: 120 g/l gelatine* 180 Bloom (Fluka 48722)		(Cl. perfringens).
Meat Liver Agar	46379	For the cultivation of anaerobic microorganisms.
Motility Nitrate Medium*	14305	Selective medium for motile nitrate-utilizing microorganisms (Ps. aeruginosa,
Supplement: 5 ml/l glycerol (Sigma 49767)		Cl. perfringens).
Nutrient Gelatin	70151	Nutrient Gelatin is recommended for the determination of gelatin-liquefying
		microorganisms as well as for the enumeration of proteolytic organisms in water by the plate count test.
Sulfite Iron Agar*	86128	For the detection and enumeration of <i>Clostridium</i> species in meat and
Supplement: 1.4 g/l ferrous sulfate (Fluka 44980)		meat products.



Clostridium perfringes

#### Table 4: Typical Media for Clostridium perfringens (Contd.)

Selective Differential Media	Cat. No.	Description
CP ChromoSelect Agar Supplements: 2 vials/I M-CP selective Supplement I (Fluka 51962) or 2 vial/I of Perfringens T.S.C. Supplement (P9352)	12398	A chromogenic Agar for enumeration detection and differentiation of <i>Clostridium sp.</i> , in particular <i>Clostridium perfringens</i> , in aqueous samples.
m-CP Agar Base Supplements: 2 vial/1 m-CP Selective Supplement I (Fluka 51962) + 2 vial/1 m-CP Selective Supplement II (Fluka 82265)	75605	Recommended by the Directive of the Council of the European Union 98/83/ EC for isolation and enumeratin of <i>Clostridium perfringens</i> from water sample using membrane filtration technique.
Perfringens Agar Base Supplements: 50 ml/l Egg Yolk Emulsion (Fluka 17148) + 2 vial/l of Perfringens S.F.P. Supplement (53436) or 2 vial/l of Perfringens T.S.C. Supplement (P9352)	39727	Used for the presumptive identification and enumeration of <i>Clostridium perfringens</i> from food.
SPS Agar*	85627	Medium for detection, isolation and enumeration of <i>Clostridium perfringens</i> and <i>Clostridium botulinum</i> in food acc. to Angelotti et al. (1962).
SPS Agar, modified	17231	For the selective isolation and enumeration of <i>Clostridium perfringens</i> from foods.
TSC Agar* Supplements: TSC Agar: 400 mg/l D-cycloserine (Sigma 30020) or 2 vial/l Perfringens T.S.C. Supplement (Fluka P9352) and 50 ml/l Egg Yolk Emulsion (Fluka 17148)	93745	For the isolation and enumeration of vegetative forms as well as spores from <i>Clostridium perfringens</i> in food, clinical specimens and other material acc.to Harmon et al. (1971).
TSC Agar egg yolk free: 400 mg/l D-cycloserine (Fluka 30020) or 2 vial/l Perfringens T.S.C. Supplement (Sigma/Fluka P9352)		
TSC Fluorescence Agar: 5.5 g/l rehydrated TSC Agar supplement (Fluka 80548) SFP Agar (Shahidi-Ferguson Perfringens Agar): 2 vial/l of Perfringens S.F.P. Supplement (Fluka P9477) and 50ml/l Egg Yolk Emulsion (Fluka 17148)		
TSN Agar*	93735	Highly selective medium for the detection and enumeration of <i>Clostridium perfringens</i> in food and other material.

\* not sold in USA

#### References

1. Sartory, D. P; Field, M.; Curbishley, S.M.; Pritchard, A.M.; Evaluation of two media for the membrane filtration enumeration of Clostridium perfringens from water. *Lett. Appl. Microbiol.* **1998**, 27, 323-327.

2. Goshu, G.; Evaluation of microbial faecal indicators and quantifying the respective level of pollution in ground and surface water of Bahidar and peri-urban areas. MSc Thesis ES 07.36 Ethiopa UNESCO-IHE Institute for Water Education, Delft, The Netherlands, **2007**.

- 3. Sherris Medical Microbiology, 4th ed., K.J. Ryan and C.G. Ray Editors, McGraw Hill, 2004.
- 4. Wells, C. L.; Wilkins, T. D., Clostridia: Sporeforming Anaerobic Bacilli, In Baron's Medical Microbiology, 4th ed., S. Baron Editor, Univ. of Texas Medical Branch, 1996.
- 5. Elmer Koneman, W. et al Color Atlas and Textbook of Diagnostic Microbiology, 5th ed., Lippincott Williams & Wilkins, 1997.
- 6. ISO 6461 Part 2. Water Quality detection and enumeration of C. perfringens. (Revision of ISO 6461-2: 1986).
- 7. Council Directive 98/83/EC, Relating to the Quality of Water Intended for Human Consumption, Official Journal of the European Communities, p.L330-32-L330/53 (1998).
- 8. National Standard Method, Enumeration of Clostridium perfringens by membrane filtration, Issue no 3.1, Reference no W5i3.1 (2005).
- 9. Orenga S.; James A.L.; Manafi M.; Perry J.D.; Pincus D.H.; Enzymatic substrates in microbiology, J. Microbiol. Meth., Vol. 79, Issue 2, p. 139-155, 2009





Burkholderia cepacia

# The *Burkholderia cepacia* Complex: A Group of Multilingual Bacteria

By Leo Eberl, Department of Microbiology, Institute of Plant Biology, University Zürich leberl@botinst.unizh.ch

# In the world of microorganisms, some members of the genus *Burkholderia* are multilingual, utilizing at least three different communication systems.


Over the past 20 years, it has become evident that bacteria exist not only as individual cells; in addition, they often coordinate their activities and act in a concerted manner similar to that of multicellular organisms. Interactions between cells are brought about by cell-to-cell communication systems that depend on the production and detection of small diffusible signal molecules. This phenomenon, which has been termed 'quorum sensing' (QS) [1], provides bacteria with a regulatory mechanism to express certain phenotypic traits in a population density-dependent manner. To date, various signal molecules belonging to different chemical classes have been described.

Some bacteria utilize more than just one communication system and therefore can be considered multilingual. One well-investigated example are the bacteria belonging to the *Burkholderia cepacia* complex (Bcc). The Bcc currently comprises a group of 17 formally-named bacterial species that, although closely related, are phenotypically diverse. Strains of the Bcc are ubiquitously distributed in nature and have been isolated from soil, water, the rhizosphere of plants, industrial settings, hospital environments, and from infected humans. Some Bcc strains have emerged as problematic opportunistic pathogens in patients with cystic fibrosis as well as in immunocompromised individuals [2].

All Bcc species investigated so far produce N-octanoylhomoserine lactone (C8-HSL), a member of the N-acyl homoserine lactone (AHL) signal molecule family, by the aid of the AHL synthase Cepl. As the cell density increases, C8-HSL accumulates in the growth medium until a critical threshold concentration is attained. At this point, C8-HSL binds to its cognate LuxR-type receptor protein CepR, which, in turn, leads to the induction or repression of target genes (Figure 1). Previous work has identified several QS-regulated functions in strains of the genus Burkholderia, including the production of extracellular proteases, chitinases, siderophores and the fungicide pyrollnitrin, swarming motility, biofilm formation and pathogenicity [3]. Moreover, some Bcc strains harbor additional AHL-based QS systems that utilize different AHL species to either fine-tune the CepIR QS system or regulate separate sets of genes. In addition to AHLs, Bcc strains were demonstrated to produce diffusible signal factors (DSFs). While most strains synthesize cis-2-dodecenoic acid, which was named the *Burkholderia* diffusible signal factor (BDSF), some Bcc species also synthesize other DSF family molecules. Subsequent studies have shown that BDSF plays a role in the regulation of bacterial virulence and biofilm formation [4]. BDSF is synthesized by the gene product of Bcam0581 and appears to be sensed by the histidine sensor kinase BCAM0227. Interestingly, the AHL- and BDSF-dependent QS systems appear to act in conjunction in the regulation of *B. cenocepacia* virulence, as a set of the AHL-controlled virulence genes is also positively regulated by BDSF [4, 5]. Furthermore, mutation of Bcam0581 results in substantially retarded energy production and exhibits impaired growth in minimal medium, highlighting the dual roles of the BDSF QS system in physiology and infection.

Finally, Bcc strains were also shown to produce a third class of signal molecules, namely 4-hydroxy-2-alkylquinolines (HAQs). Recent work has shown that *Burkholderia ambifaria* (a member of the Bcc) HAQs contain predominantly an unsaturated aliphatic side chain and are typically methylated (at the 3 position), hence their designation as 4-hydroxy-3-methyl-2-alkylquinolines (HMAQs) [6]. Interestingly, inactivation of HAQs' production of *B. ambifaria* has led to an increased AHL production and, consequently, affected phenotypes that are under the control of QS in *B. ambifaria*, namely antifungal activity, siderophore production, and proteolytic activity.


**Figure 1:** Bcc strains can utilize up to three different cell-to-cell communication systems that rely on different classes of signal molecules: N-acyl-homoserine lactones (particularly C8-HSL), cis-2-dodecenoic acid, the Burkholderia diffusible signal factor (BDSF), and 4-hydroxy-3-methyl-2-alkylquinolines (HMAQs).




At present, it is unclear why Bcc strains synthesize multiple signal molecules. However, it should be kept in mind that these bacteria inhabit various habitats and thus may encounter a large diversity of bacteria. It is therefore tempting to speculate that the utilization of different communication systems will enable the bacteria to interact with a wide range of indigenous bacteria through chemical cross-talk. These interspecies interactions may be important for the fitness of the bacteria in their natural environment.

#### References

- 1. C. Fuqua, S.C. Winans, E.P. Greenberg, "Census and consensus in bacterial ecosystems: the LuxR-LuxI family of quorum-sensing transcriptional regulators", Annu. Rev. Microbiol. 50, 727-51 (1996)
- 2. E. Mahenthiralingam, T.A. Urban, J. B. Goldberg, "The multifarious, multireplicon Burkholderia cepacia complex," Nat. Rev. Microbiol. 3, 144-56 (2005)
- 3. P.A. Sokol, R.J. Malott, K. Riedel, L. Eberl, "Communication systems in the genus Burkholderia: global regulators and targets for novel antipathogenic drugs", Future Microbiol. 2, 555-63 (2007)
- 4. Y. Deng, C. Boon, L. Eberl, L.H. Zhang, "Differential modulation of Burkholderia cenocepacia virulence and energy metabolism by the quorum-sensing signal BDSF and its synthase", J. Bacteriol. 191(23), 7270-8 (2009)
- Y. McCarthy, L. Yang L, Twomey KB, Sass A, Tolker-Nielsen T, Mahenthiralingam E, Dow JM, Ryan RP, "A sensor kinase recognizing the cell-cell signal BDSF (cis-2-dodecenoic acid) regulates virulence in Burkholderia cenocepacia", Mol. Microbiol. 77(5), 1220-36 (2010)
- L. Vial, F. Lépine, S. Milot, M.C. Groleau, V. Dekimpe, D.E. Woods, E. Déziel, "Burkholderia pseudomallei, B. thailandensis, and B. ambifaria produce 4-hydroxy-2-alkylquinoline analogues with a methyl group at the 3 position that is required for quorum-sensing regulation", J. Bacteriol. 190(15), 5339-52 (2008)



## Sigma-Aldrich – Number One for Quorum Sensing Compounds



General chemical structure of an N-acyl homoserine lactone

Since 1998, Sigma-Aldrich has produced many Quorum Sensing Compounds and performed studies on their cell signaling reaction. The beginning efforts focused on simple N-acyl homoserine lactones (AHL) and an application test with the induction of violacein expression (violet pigment) in a *Chromobacterium violaceum* mutant (not able to produce AHL). Today the signaling molecules are emerging as more complex, while at the same time inhibiting compounds are growing in importance for research.

For more information, please visit our homepage: *sigma-aldrich.com/qs* 

#### **Quorum Sensing Substances**

	Cianallin a		
Cat. No.	Signalling Compounds	Substance	
49619	BDSE	cis-2-Dodecenoic acid	
10942	C4-HCT	N-Butyryl-DL-homocysteine thiolactone	
09945	C4-HSI	N-Butyryl-DL-homoserine lactone	
09926	C6-HSL	N-Hexanoyl-DL-homoserine lactone	
10939	C7-HSI	N-Heptanoyl-DL-homoserine lactone	
10940	C8-HSL	N-Octanoyl-DL-homoserine lactone	
17248	C10-HSL	N-Decanoyl-DL-homoserine lactone	
17247	C12-HSI	N-Dodecanoyl-DL-homoserine lactone	
10937	C14-HSL	N-Tetradecanoyl-DL-homoserine lactone	
42052	DSF	cis-11-Methyl-2-dodecenoic acid	
53727	3-Hydroxy-C12-HSL	N-(3-Hydroxydodecanoyl)-DL-homoserine lactone	
51481	3-Hydroxy-C14-HSL	N-(3-Hydroxytetradecanoyl)-DL-homoserine lactone	
K3255	3-Oxo-C6-HSL	N-(3-Oxohexanoyl)-DL-homoserine lactone	
K3007	3-Oxo-C6-HSL	N-(3-Oxohexanoyl)-DL-homoserine lactone	
O1639	3-Oxo-C8-HSL	N-(3-Oxooctanoyl)-DL-homoserine lactone	
01764	3-Oxo-C8-HSL	N-(3-Oxooctanoyl)-L-homoserine lactone	
O9014	3-Oxo-C10-HSL	N-(3-Oxodecanoyl)-L-homoserine lactone	
O9139	3-Oxo-C12-HSL	N-(3-Oxododecanoyl)-L-homoserine lactone	
09264	3-Oxo-C14-HSL	N-(3-Oxotetradecanoyl)-L-homoserine lactone	
07077	pC-HSL	N-(p-Coumaroyl)-L-homoserine lactone	
94398	PQS (Pseudomonas	2-Heptyl-3-hydroxy-4(1H)-quinolone	
	quinolone signal)		
Cat. No.	Antagonist/Inhibito	ors	
53796	(Z-)-4-Bromo-5-(bromomethylene)-2(5H)-furanone (Furanone C-30)		

68749 3-Oxo-N-(2-oxocyclohexyl)dodecanamide



## Sigma-Aldrich<sup>®</sup> Worldwide Offices

#### Argentina

Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

#### Australia

Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

#### Austria

Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

#### Belgium

Free Tel: 0800 14747 Free Fax: 0800 14745 Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

#### Brazil

Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

#### Canada

Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

#### Chile

Tel: (+56) 2 495 7395 Fax: (+56) 2 495 7396

#### China

Free Tel: 800 819 3336 Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

**Czech Republic** Tel: (+420) 246 003 200 Fax: (+420) 246 003 291

#### Denmark

Tel: (+45) 43 56 59 00 Fax: (+45) 43 56 59 05

#### Finland

Tel: (+358) 9 350 9250 Fax: (+358) 9 350 92555

#### France

Free Tel: 0800 211 408 Free Fax: 0800 031 052 Tel: (+33) 474 82 28 88 Fax: (+33) 474 95 68 08

#### Germany

Free Tel: 0800 51 55 000 Free Fax: 0800 64 90 000 Tel: (+49) 89 6513 0 Fax: (+49) 89 6513 1160

#### Hungary

Ingyenes telefonszám: 06 80 355 355 Ingyenes fax szám: 06 80 344 344 Tel: (+36) 1 235 9063 Fax: (+36) 1 269 6470

#### India

**Telephone** Bangalore: (+91) 80 6621 9400 New Delhi: (+91) 11 4358 8000 Mumbai: (+91) 22 2570 2364 Hyderabad: (+91) 40 4015 5488 Kolkata: (+91) 33 4013 8003

#### Fax

Bangalore: (+91) 80 6621 9550 New Delhi: (+91) 11 4358 8001 Mumbai: (+91) 22 4087 2364 Hyderabad: (+91) 40 4015 5488 Kolkata: (+91) 33 4013 8000

#### Ireland

Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: (+353) 402 20370 Fax: (+ 353) 402 20375

#### Israel

Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4100 Fax: (+972) 8 948 4200

#### Italy

Free Tel: 800 827 018 Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

#### **Japan** Tel: (+81) 3 5796 7300

Fax: (+81) 3 5796 7315

#### Korea

Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111 Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090

**Malaysia** Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

#### Mexico Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: (+52) 722 276 1600 Fax: (+52) 722 276 1601

 The Netherlands

 Free Tel: 0800 022 9088

 Free Fax: 0800 022 9089

 Tel: (+31) 78 620 5411

 Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666 Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

#### **Norway** Tel: (+47) 23 17 60 00 Fax: (+47) 23 17 60 10

Poland Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

#### Portugal

Free Tel: 800 202 180 Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

### Russia

Tel: (+7) 495 621 5828 Fax: (+7) 495 621 6037

#### **Singapore** Tel: (+65) 6779 1200 Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562 Fax: (+421) 255 571 564

#### South Africa Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

#### Spain

Free Tel: 900 101 376 Free Fax: 900 102 028 Tel: (+34) 91 661 99 77 Fax: (+34) 91 661 96 42

#### Sweden

Tel: (+46) 8 742 4200 Fax: (+46) 8 742 4243

#### Switzerland Free Tel: 0800 80 00 80

Free Fax: 0800 80 00 81 Tel: (+41) 81 755 2828 Fax: (+41) 81 755 2815

#### **United Kingdom** Free Tel: 0800 717 181 Free Fax: 0800 378 785

Tel: (+44) 1747 833 000 Fax: (+44) 1747 833 313

#### United States

Toll-Free: 800 325 3010 Toll-Free Fax: 800 325 5052 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

**Vietnam** Tel: (+84) 3516 2810 Fax: (+84) 6258 4238

Internet sigma-aldrich.com

Enabling Science to Improve the Quality of Life Order/Customer Service (800) 325-3010 • Fax (800) 325-5052 Technical Service (800) 325-5832 • sigma-aldrich.com/techservice Development/Custom Manufacturing Inquiries **SAFC**<sup>•</sup> (800) 244-1173 Safety-related Information sigma-aldrich.com/safetycenter World Headquarters 3050 Spruce St. St. Louis, MO 63103 (314) 771-5765 sigma-aldrich.com

©2010 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, LP. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, publications. Pworka publications. Pworkaser must determine the suitability of the product(S) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.