



Microbiology Focus

Volume 4.2, 2012



Legionella pneumophila, an Unpredictable Pathogen



Steam from a cooling tower is one possible way that *Legionella* can be spread.

| | |
|---------------------------------------|---|
| Legionnaires' Disease..... | 2 |
| Media Ingredients | 5 |
| iPad/Photography Competition | 7 |



Legionnaires' Disease

By Jvo Siegrist, Product Manager Microbiology.... ivo.siegrist@sial.com

A Potentially Fatal Form of Pneumonia

Legionnaires' disease is a dangerous and infectious pneumonia which can affect anybody, but those at highest risk are the elderly, sick, smokers, or other people with weak immune systems.

The causative organism is the *Legionella pneumophila*, which easily forms aerosols (fine airborne particles) which can be inhaled, causing infection. The natural source of *Legionella pneumophila* is all types of water, but particularly water containing algae, rust, sludge, scale, and other organic compounds. It is known that *Legionella* are able to build biofilm inside water pipes, where they are protected against antimicrobial treatments. The infected water droplets can be produced in such places as whirlpool spas or air conditioners, but they can also occur in freshwater such as rivers and ponds. However, most people exposed to *Legionella* do not become ill, and it is not possible for Legionnaires' disease to be transmitted from person to person.

The name Legionnaires' disease originated with an outbreak in Philadelphia, Pennsylvania at an American Legion convention (1976). In 1977, the new pathogen was identified and named *Legionella* (Figure 1). Symptoms are flu-like, with chills, muscle pains, headache, cough and high fever. Death can result from severe pneumonia. The Philadelphia outbreak sickened hundreds and resulted in 34 deaths. Additionally, *Legionella* is also the causative organism of Pontiac fever, which is a milder respiratory illness resembling acute influenza.

A serious outbreak of *Legionella pneumophila* occurred in 2009/2010 in Neu Ulm, Germany. With 64 cases and 5 deaths, it was the worst outbreak of legionellosis in Germany. Based on the molecular biological profile, it was assumed that the source was a large, wet cooling tower. About 400 incidences of legionellosis are reported in all of Germany each year.

In November 2011, a new drinking water ordinance was put into place in Germany which regulates large industrial water heating systems. Once a year a sample has to be tested, and the level of *Legionella* should be below 100 cfu per 100 mL.

Legionella pneumophila is an obligate aerobic gram-negative, rod-shaped bacterium with monopolar flagella. Often the size of the bacteria in the culture is very variable (2-20 µm, Figure 2), and it is commonly found in aquatic environments. The organism can survive a wide range of conditions, including temperatures from 0 to 63 °C, and pH values from 5.0 to 8.5. Because *Legionella* has cysteine as a growth requirement, it does not grow on common blood agar media. *Legionella pneumophila* primarily uses amino acids as an energy source and also needs iron (III) ions. This is probably the reason that in nature, *Legionella pneumophila* is often found together with iron bacteria and amoeba (Figure 3). The optimal temperature range for growth of *Legionella* is given at 25 to 50 °C.

Did you know ...

that smokers are easier targets for bacteria?

According to a new study, it was observed that beneficial bacteria are inhibited in the mouths of smokers, resulting in greater susceptibility to other bacteria (Source: Science Daily). Information about the ratio of beneficial flora in the lungs is not yet available, but this study may provide insight into the reason why *Legionella pneumophila* thrives in smokers' lungs.

Figure 1: Smoking person



Figure 2: Colorized scanning electron micrograph (magnification 10000X) depicted a number of gram-negative *Legionella pneumophila* bacteria. Colony was grown on BCYE medium without antibiotic addition. (Source: CDC/ Margaret Williams, PhD; Claressa Lucas, PhD; Tatiana Travis, BS, 2009).




Table 1: Discs with certified reference test strains of *Legionella* (Restrictions may apply in some countries)

| Vitroids™ Test Strains | Origin | Strain No. | CFU | Cat. No. |
|-------------------------------|--------|------------|-------|----------|
| <i>Legionella bozemanii</i> | NCTC | 11368 | 50000 | RQC02908 |
| <i>Legionella pneumophila</i> | ATCC | 12821 | 50000 | RQC02008 |



Figure 3: Colorized electron micrograph depicts an amoeba, *Hartmannella vermiformis* (orange) as it entraps a *Legionella pneumophila* bacterium (green) with an extended pseudopod. After it is ingested, the *Legionella pneumophila* bacterium can survive as a symbiont within what then becomes its protozoan host. The amoeba then becomes what is referred to as a "Trojan horse", for by harboring the pathogenic bacteria, the amoeba can afford them protection, and in fact, in times of adverse environmental conditions, are able to metamorphose into a cystic-stage enabling it, and its symbiotic resident pathogens to withstand such environmental stresses. (Source CDC/ Dr. Barry S. Fields)

| | |
|-----------|-------------------------|
| Kingdom: | Bacteria |
| Division: | Proteobacteria |
| Class: | Gamma Proteobacteria |
| Order: | Legionellales |
| Family: | Legionellaceae |
| Genus: | <i>Legionella</i> |
| Species: | <i>pneumophila</i> |



The traditional method to detect *Legionella* is based on buffered charcoal yeast extract (BCYE) agar. For growth, the organism requires the presence of cysteine and therefore does not grow on common blood agar media. Common laboratory procedures recommend concentrating the sample by centrifugation and/or filtration through 0.2 µm filters before inoculation. For selective isolation, antibiotics like polymyxin B, anisomycin, vancomycin, natamycin, cycloheximide, cefamandole, trimethoprim, colistin sulphate, amphotericin B, bromo thymol blue and bromo cresol purple are added. Additionally, selectivity can be attained by applying heat or acid. Some typical media can be found in **Tables 2 and 3**.

Table 2: *Legionella* Media

| Medium | Cat. No. |
|---|----------|
| Legionella Agar, Base | 74303 |
| Feeley Gorman Agar | 05598 |
| Buffered Charcoal Yeast Extract (BCYE) Agar, Base | 86558 |

Table 3: BCYE Agar based media

| | Base Medium | Supplements |
|---------------------------------|----------------------|---|
| Legionella MWY Selective Agar: | BCYE Agar, Base | Legionella Supplement (Twin Pack; Fluka 89166), Legionella Selective Supplement IV (Fluka 94029) |
| Legionella GVPC Selective Agar: | BCYE Agar, Base | Legionella (GVPN) Selective Supplement (Fluka 43509), Legionella Supplement (Twin Pack; Fluka 89166) |
| Legionella GVPN Selective Agar: | BCYE Agar, Base | Legionella (GVPC) Selective Supplement (Fluka 61025), Legionella Supplement (Twin Pack; Fluka 89166) |
| Legionella BMPA Selective Agar | BCYE Agar, Base | Legionella Supplement (Twin Pack; Fluka 89166), BMPA Supplement: cefamandole 4 mg/L, anisomycin 80 mg/L, polymyxin B 80'000 IU/L |
| BCYE Agar without L-cysteine: | Legionella Agar Base | Legionella BCYE Supplement without L-Cysteine (Fluka 43753) |

Table 4: *Legionella* Supplements

| Supplement Content (per vial sufficient for 500 mL) | Legionella Supplement (Fluka 89166) | Selective Supplement IV (Fluka 94029) | GVPN (Fluka 43509) | GVPC (Fluka 61025) | BCYE Supplement w/out cysteine (Fluka 43753) | Legionella Selective Supplement (Fluka 18284) |
|---|---|---|-----------------------|-----------------------|--|---|
| ACES Buffer /KOH | | | | | 5 g | |
| L-Cysteine hydrochloride | 200 mg | | | | | |
| Ferric pyrophosphate, soluble | 125 mg | | | | 125 mg | |
| α-Ketoglutarate | | | | | 0.5 g | |
| Polymyxin B sulphate | | 25,000 Units | 40,000 IU | 39,600 IU | | |
| Glycine | | 1.5 g | 1.5 g | 1.5 g | | |
| Anisomycin | | 40 mg | | | | |
| Vancomycin | | 0.5 mg | 0.5 mg | 0.5 mg | | 2.5 mg |
| Natamycin | | | 20 mg | | | |
| Cycloheximide | | | | 40 mg | | |
| Colistin sulphate | | | | | | 7,500 Units |
| Trimethoprim | | | | | | 1.25 mg |
| Amphotericin B | | | | | | 1.25 mg |
| Bromo thymol blue | | 5 mg | | | | |
| Bromo cresol purple | | 5 mg | | | | |



Figure 4: Detection and Enumeration of *Legionella* in Water acc. ISO 11731 (1998), ISO 11731-2 (2008)




Table 5: HybriScan® Kits for Detection and Identification of *Legionella* Species

| Brand | Name | Description | Cat. No. |
|-------|--|---|----------|
| Fluka | HybriScan ^D <i>Legionella</i> | Detection of <i>Legionella</i> , including <i>L. pneumophila</i> in water supplies and air-conditioning systems | 16593 |
| Fluka | HybriScan ^D <i>Legionella pneumophila</i> | Detection and identification of <i>L. pneumophila</i> in water supplies and air-conditioning systems | 07190 |
| Fluka | HybriScan ^I <i>Legionella pneumophila</i> | Identification of <i>Legionella pneumophila</i> | 49417 |


New Technology for *Legionella* Detection

The current procedure used to detect and identify the *Legionella* species takes a lot of time and requires several steps; therefore, new methods for the detection of *Legionella* are of great interest. A modern detection target is rRNA, which is more numerous than DNA, and is only present in living cells. The test is performed on a microtiter plate and takes less than 2.5 hours. Cell count quantification is possible with this kit using photometric methods. Compared to PCR, our system does not count dead cells, is much easier to use (see work flow in **Figure 5**), is less expensive, and is not affected by the sample matrix. Confirming its utility, our new test system won several innovation awards, and will surely become a routine method to detect *Legionella*. The rapidity, sensitivity, reliability, robustness, adaptability to sample matrix, and time-savings meet today's analytical microbiology demands.

Specificity is achieved by targeting conserved or unique rRNA sequences. A biotin-labeled capture probe is used to immobilize the target sequence on a solid support plate (streptavidin-coated microtiter plate). A digoxigenin-labeled detection probe provides an enzyme-linked optical signal read out. Detection results from application of anti-DIG-horseradish peroxidase Fab fragments. The bound complex is visualized by horseradish peroxidase substrate TMB (3,3',5,5'-tetramethylbenzidine). Photometric data are measured at 450 nm and compared with standard solutions.

For more details of our new *Legionella* detection kit, please visit sigma-aldrich.com/hybriscan

Figure 5: Work Flow Process of the HybriScan *Legionella* Kit



Base Ingredients of Microbiology Media

Jvo Siegrist, Product Manager Microbiology...jvo.siegrist@sial.com

The quality of a media and its technical features depend on the base ingredients, ranging from simple sugars to peptones, salts, antibiotics, and more complex indicators.

A medium has one primary function -- to promote the growth of organisms. The components of a medium are often based upon the organism's natural habitat. For example, an organism growing on meat may require meat peptone, and an organism growing on nutrients with a high carbohydrate content may thrive on malt extract. In addition to this growth purpose, media may serve in a number of other applications, including the differentiation and identification of organisms, the selective isolation or enrichment of organisms, and the study of a certain reaction of the organism. A vast array of peptones, extracts and other additives is available to promote and sustain the growth of most organisms.

Figure 1: Diverse Ingredients Build the Base for All Kinds of Media



Proteins (Protein Hydrolysate, Amino Acids, etc.)

Although synthetic growth media are available, most media still use complex compounds, such as peptone or yeast extract, since synthetic media lack the complexity and richness of nutrients. Peptones and protein extracts are excellent natural sources of amino acids, peptides, proteins and many other growth factors. They are most often obtained by enzymatic digestion or acid hydrolysis of natural protein sources, such as animal tissues, milk, plants or microbial cultures. The range of available peptones is extensive and comprises a major role in the growth conditions of most organisms (Table 1).

Carbohydrates (Extracts, Sugars, etc.)

Carbohydrates are an important energy source. Mono-, di-, oligo- and polysaccharides, as well as natural extracts like rice or malt extracts, provide a versatile possibility of substrates for mold or bacteria cultures. They can also be used to make the media more selective or to identify fermentation profiles. Today, a broad range of media with chromogenic substrates is available (Table 2).

Biological Acids

Pyruvate, one type of biological acid, is known to promote growth and to improve the recovery rate. Other acids such as orange extracts, citric or acetic acid are also used for selective growth.

Buffering Agents

Potassium phosphates are the primary agents used for the buffering system.

Salts

Sodium chloride is used primarily for osmotic balance; however, it can also be used to make the medium more selective to halophilic and halotolerant bacteria. In addition, other salts such as lithium chloride or ammonium bismuth citrate are also used to make the medium more species specific.

Fatty Acids and Lipids

Fatty acids and lipids, such as lecithin, are necessary nutrients and a valuable source of proteins. Fluka offers egg powder and liquid sterile egg supplements, as well as pure lipids and fatty acids.

Vitamins and Trace Elements

Yeast extract, present in numerous complex media, is the most common source for Vitamin B12 (Table 3). Yeast extract also contains a large number of amino acids, additional vitamins, and trace elements. Some media also commonly incorporate the addition of pure vitamins and trace elements.

Selective Agents (Detergents, Bile Salts, Antibiotics, etc.)

Bile is often used as an inhibitory agent against most gram-positive bacteria. Cholates (Table 5), biological detergent-like compounds with anti-microbial activity, are major constituents of bile. Alternatively, SDS and other detergents are used for the same purpose. For the most part, however, selective agents are comprised of antibiotics that are often added as a mixture in supplemental vials (Figure 1).

Indicators and Dyes

These help to indicate biochemical properties or metabolic pathways and are vital for the identification and differentiation of organisms.

Agar

Agar is the solidifying agent in solid growth media, and its selection should be carefully considered based upon certain criteria and dependent upon the application. For example, when high transparency and brightness is needed, as in nutritional studies (Vitamin Assay Media) and sensitivity testing procedures, or when high purity and efficient diffusion of substances is essential, a highly purified agar (Fluka 05038) is recommended. For identification and differentiation, we recommend using a purified or even highly purified agar. However, when isolating a single colony, a standard quality will suffice in most cases. Typical solid media have an agar concentration of 1.0 - 1.5% to accommodate the requirements of different applications and the growth habits of target microorganisms (Table 4).



Table 1: Common Protein Sources for Media

| Cat. No. | Description |
|----------|---|
| A2427 | Amicase |
| B4888 | Beef extract |
| B3551 | Biopeptone |
| 53283 | Brain Heart Infusion |
| 75917 | Brain Heart Infusion, Porcine* |
| 93491 | Broadbean Peptone |
| C7970 | Casein acid hydrolysate vitamin free |
| 22090 | Casein Hydrolysate |
| 39396 | Casein Yeast Peptone |
| C4773 | Corn gluten meal |
| 55871 | Egg powder* |
| 93490 | Fish Peptone* |
| 49760 | Gluten Hydrolysate from maize |
| 92498 | Heart extract from bovine heart* |
| C0501 | Hy-Case® Amino |
| C9386 | Hy-Case® SF |
| 57462 | Infusion powder from bovine heart |
| 61300 | Lactalbumin Hydrolysate |
| 61302 | Lactalbumin Hydrolysate* |
| 03077 | Liver Hydrolysate |
| 70164 | Meat Extract |
| C0626 | N-Z-Amine® A |
| C1026 | N-Z-Case® |
| 18332 | Peptone (vegetable) |
| 51841 | Peptone (vegetable) acid hydrolysate |
| 19942 | Peptone (vegetable), no. 1 |
| 61854 | Peptone (vegetable), no. 2 |
| 77180 | Peptone from animal proteins* |
| 70173 | Peptone from casein and other animal proteins |
| 70171 | Peptone from casein, acid digest* |
| 82303 | Peptone from casein, enzymatic digest |
| 70169 | Peptone from casein, pancreatic digest |
| 70172 | Peptone from casein, tryptic digest |
| 70951 | Peptone from gelatin, enzymatic digest |
| 70176 | Peptone from gelatin, pancreatic digest |
| P0521 | Peptone from Glycine max (soybean) |
| 93733 | Peptone from meat and soybean meal* |
| 82962 | Peptone from meat, enzymatic digest |
| 70174 | Peptone from meat, peptic digest* |
| 96174 | Peptone from pea |
| 93926 | Peptone from porcine heart |
| 83059 | Peptone from potatoes |
| 70178 | Peptone from soybean meal, enzymatic digest |
| 90765 | Peptone from soybean, enzymatic digest |
| 87972 | Peptone from soybean, enzymatic digest |
| P6463 | Peptone Hy-Soy® T |
| P6713 | Peptone N-Z-Soy® BL 7 |
| P4838 | Peptone Primatone® HS |
| P4963 | Peptone Primatone® RL |
| P5088 | Peptone Primatone® RL T |
| 68971 | Peptone special |
| 92976 | Peptone special (vegetable) |
| 77199 | Peptone, mycological |
| P8388 | Primatone® |
| 82514 | Protein Hydrolysate Amicase |
| 82524 | Protein Hydrolysate N-Z-Amine® AS* |
| 29185 | Proteose Peptone (vegetable) |
| P0431 | Proteose Peptone Enzymatic hydrolysate |
| 82450 | Proteose-Peptone |
| 70166 | Skim Milk Powder |
| S1674 | Soy protein acid hydrolysate |
| T7293 | Tryptone |
| 16922 | Tryptone (vegetable) |
| 95039 | Tryptone enzymatic digest from casein |

| | |
|-------|-----------------------------------|
| 61044 | Tryptone Plus |
| 70937 | Tryptose |
| T2813 | Tryptose |
| 12331 | Tryptose (vegetable) |
| 05138 | Vegetable Extract |
| 04316 | Vegetable Extract no. 1 |
| 49869 | Vegetable Extract no. 2 |
| 07436 | Vegetable hydrolysate no. 2 |
| 67381 | Vegetable Infusion powder |
| 95757 | Vegetable Special Infusion powder |
| 93492 | Wheat Peptone |

Table 2: Carbohydrate Sources in Microbial Media Quality

| Cat. No. | Description |
|----------|----------------------------------|
| 10850 | D-(-)-Arabinose |
| 22150 | D-(+)-Cellobiose |
| 22160 | D-(+)-Cellobiose octaacetate |
| 31400 | Dextrin from potato starch |
| 31405 | Dextrin from potato starch |
| 47740 | D-(-)-Fructose |
| 48260 | D-(+)-Galactose |
| 49159 | D-(+)-Glucose monohydrate |
| 49200 | alpha-D-(+)-Glucose pentaacetate |
| 57570 | myo-Inositol |
| 70167 | Malt Extract |
| 63560 | D-Mannitol |
| 63580 | D-(+)-Mannose |
| 63582 | D-(+)-Mannose |
| 63620 | D-(+)-Melezitose Monohydrate |
| 63630 | D-(+)-Melibiose |
| 66940 | Methyl alpha-D-glucopyranoside |
| 67770 | Methyl alpha-D-mannopyranoside |
| 07915 | Potato Extract |
| 83400 | D-(+)-Raffinose pentahydrate |
| 83650 | L-(+)-Rhamnose Monohydrate |
| 84100 | Sucrose |
| 90210 | D-(+)-Trehalose Dihydrate |

Table 3: Yeast Extracts for Media and Fermentations

| Cat. No. | Description |
|----------|---|
| 73145 | Yeast Autolysate |
| 92144 | Yeast Extract |
| 70161 | Yeast Extract (premium quality) |
| 09182 | Yeast Extract for technological purpose |
| 89526 | Yeast Extract micro-agglomerated |

Table 4: Agars for Microbiology

| Cat. No. | Description |
|----------|---|
| 05040 | Agar standard |
| 05039 | Agar purified |
| 05038 | Agar highly purified |
| 50524 | Agar anti-swarming |
| 91411 | Agar for Baird Parker Agar |
| 05729 | Agar for chromogenic media |
| 42146 | Agar for membrane filtration, low gel strength |
| 56763 | Agar high purity, low ionic content, low gel strength |

Table 5: Bile Salts

| Cat.No. | Description |
|---------|-------------------------------|
| B8381 | Bile from bovine and ovine |
| 48305 | Bile Salts* |
| 70168 | Ox-bile, dehydrated, purified |

Win an iPad in the 2012 Fluka Microbiology Photography Competition



This photography competition is sponsored by Sigma-Aldrich with the aim of encouraging microbiologists to promote some aspect of their work or their field of research. The best photographic entries with the best description of the photograph's subject will win prizes such as an iPad, a Swiss army knife, a USB flashdrive, and a laser pointer.

The winning images will be published in Microbiology Focus and the best one will have the distinction of being featured on the cover.



Rules of the Competition and Conditions of Entry

1. The competition is open to all residents worldwide.
2. Entries should illustrate any microorganisms (living or dead) or a microbiologist in action at work.
3. Picture size should be at least 400 dpi and 90 x 120 mm (max 3 MB). The file format must be in jpg, tiff or pdf!
4. The entries will be judged on:
 - clarity of presentation
 - composition
 - illumination and contrast
 - congruency of subject matter and title of photograph
 - scientific interest and relevance
 - originality
5. Winning entries will be retained by Sigma-Aldrich, who will have sole rights of publication, reproduction and display.
6. Closing date for contest entries will be July 31, 2012.
7. Entries received after the closing date will not be considered. Entries received incomplete, illegible, mutilated, altered or not complying exactly with the instructions and theme may be disqualified.
8. Decisions of the judges in all matters affecting the competition will be final and legally binding.

The competition will be judged by:

Dr. Lars Fieseler

Zurich University of Applied Sciences - ZHAW
Supervisor, Department Microbiology

Prof. Mohammad Manafi

Medical University of Vienna
Head of Department for Food Hygiene

Jvo Siegrist, Sigma-Aldrich

Product Manager, Microbiology

Method of Entry

There is no entry fee, but an entry form must be completed for each entry (a maximum of two entries may be entered)

Entry forms are available from our website
sigma-aldrich.com/fluka-mibi-competition

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com



*Enabling Science to
Improve the Quality of Life*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

©2012 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

ONL
T412028

SIGMA-ALDRICH®