Microbiology Focus Volume 5.3, 2013



Fluorescence Assays in Microbiology

Fluorescence assays are interesting tools used for the detection of microorganisms. Fluorogenic substrates and fluorescence labelling of organisms are possible, as are sensitive methods of detection either by using a fluorescence reader, viewing the sample under the fluorescence microscope, or by simply viewing the sample directly with the naked eye.

Image from CDC: Mixed culture of microorganisms stained with fluorescent antibodies (1959)

FluoroSELECT[™] Assays2

Fluorogenic Media4

Quantitative Certified Control Strains6

SIGMA-ALDRICH®

FluoroSELECT[™] Assays

Jvo Siegrist, Product Manager Microbiology — ivo.siegrist@sial.com

FluoroSELECT is a rapid test based on an assay with fluorogenic substrates.



The FluoroSELECT Assay system is a presumptive screen intended for the rapid detection of E. coli and Gram-negative bacteria from a surface swab sample using the FluoroSELECT handheld fluorometer and the fluorescent assay. It can also be used as a confirmation test after using a classical or modern method. E. coli is commonly found in the lower intestine of warm-blooded organisms and is therefore an indicator for fecal contamination. Most strains are harmless, but some serotypes can cause serious food poisoning in humans, and are occasionally responsible for product recalls due to food contamination. Gram-negative bacteria are bacteria that are defined by their inability to retain crystal violet dye in the Gram staining protocol. This test was and still is vital for the differentiation and the classification of bacteria. Compared with Gram-positive bacteria, Gram-negative bacteria are more resistant to antibiotics because of their impenetrable cell wall.

FluoroSELECT utilizes specific fluorogenic substrates which are cleaved from characteristic enzymes of *E. coli* or Gram-negative bacteria. During the typical peptide hydrolysis, the specific enzymes hydrolyze the fluorogenic substrates and produce a fluorescence which is read by a portable, low cost fluorometer. The excitation wavelength is at 360 nm (UV), while the emission wavelength is at 460 nm (blue light), the intensity of which is then measured to decide if the sample is positive or negative. It is a rapid test since the assay takes only 30 minutes, although it may require some incubation time. Because the enzymes are specific, the assay is highly specific and robust. In addition, since it is a typical fluorescence assay, it is also highly sensitive and small amounts of the cleaved fluorogen will already have been detected. The detection performance for both tests is immediate detection of 100,000 cfu (colony forming units), while 100 cfu are detected in 8 hours, and 1 cfu is detected within 10 hours of incubation time.

Simplified testing procedures:

- 1. A cotton swab is used to wipe a particular surface of the detection area. The swab is then inserted into 0.5 mL of distilled water for an immediate test, or into 0.5 mL of incubation solution for an incubation period ranging from 3 to 10 hours at 38.5 °C, before measuring. Add an enzyme inducer after inserting the swab into the water or incubation solution.
- 2. Then add a lysing agent to the solution and let incubate for 5 minutes.
- 3. Next, a substrate is added to the solution and allowed to stand for 5 minutes.
- 4. Afterwards, the solution is transferred into a glass test tube and then inserted into the Fluorometer. The first measurement is taken immediately, and after 20 minutes, another measurement is taken. The findings can then be used to determine a positive or negative result.









Did you know...

Gram-negative bacteria are generally less sensitive to antibiotics?

The outer membrane of Gram-negative bacteria contributes to this intrinsic resistance by acting as an efficient permeability barrier, because the narrow porin channels limit the penetration of hydrophilic solutes and the low fluidity of the lipopolysaccharide leaflet slows down the inward diffusion of lipophilic solutes.

[Source: P. Plésiat, et al., Mol. Microbiol. 1992]



| Cat. No. | Description | Qty. |
|----------|---|----------|
| 53649 | FluoroSELECT <i>E. coli</i> Assay Kit | 50 Tests |
| 91333 | FluoroSELECT Gram-Negative Assay Kit | 50 Tests |
| Z805726 | FluoroSELECT single channel fluorometer (λex 360 nm; λem 460 nm) | 1 ea. |
| Z805823 | Glass vials for FluoroSELECT fluorometer (O.D. \times L 6 mm \times 25 mm, volume nominal capacity 200 μ L) | 100 ea. |

Table 1: Product list of microbiology FluoroSELECT kits and equipment

There are more enzymatic tests available for quantitative detection of acetate, ammonia, ascorbic acid, formaldehyde, glycerol and lactose with the FluoroSELECT system. Visit **sigma-aldrich.com/fluoroselect** or refer to **Table 2**.

| Cat. No. | Description | Qty. |
|----------|---------------------------------------|----------|
| 53659 | FluoroSELECT Ammonia Kit | 50 Tests |
| 89872 | FluoroSELECT Formaldehyde Kit | 50 Tests |
| 00254 | FluoroSELECT Glycerol Kit* | 50 Tests |
| 76691 | FluoroSELECT Ascorbic acid Kit* | 50 Tests |
| 76786 | FluoroSELECT Acetate Assay Kit | 50 Tests |
| 91218 | FluoroSELECT Lactose Detection Assay* | 50 Tests |

 Table 2: Non microbiology FluoroSELECT assays

 *Inquire on availability.



Microbiology Special Catalog

This is the premiere Sigma-Aldrich catalog for microbiology, and it contains over 900 specific microbiology products. The tests, media, reagents, microorganisms and base ingredients are sorted to help microbiologists find the best test, media and additives for the desired microbiological analysis.

Contents:

- Microorganisms Standards (Vitroids™)
- Proficiency Testing
- Molecular Biological, Immunological, Biochemical Tests
- Microscopical Kits and Reagents
- Media for Specific Organism, Special Application, or ISO Guideline
- Supplements
- Base Ingredients for Media



*only available as long as supplies last



Analytica





Fluorogenic Media

Jvo Siegrist, Product Manager Microbiology — ivo.siegrist@sial.com

Fluorogenic media combine traditional methods with new knowledge for beneficial results.

Traditional validated media versus improved media formulations containing chromogenic and fluorogenic substrates is an important topic in the field of microbiology. The focus behind such developments is the production of media that would make the detection and identification of microorganisms more rapid and more reliable. Diverse chromogenic and fluorogenic substrates such as ONPG, X-Gal, X-Glu MUP or MUG are available and are already combined in today's media formulations. Fluorogenic substrates have the disadvantage of requiring a UV-lamp, but the fluorogenic media are more sensitive and give an earlier detectable result. They can also be combined with a chromogenic system without interference. These reasons give such media additional value, and together with the specified selectivity of the medium and the substrate, new possibilities are available to solve problems in a smart way.



Figure 2: Fluorogenic media containing *E. coli* strains and *Pseudomonas putida* (right at the bottom) and *Lactobacillus* spp. (left).

Did you know...

about upconversion luminescence?

In normal cases, the color of emitted light requires a higher energy light source (in most cases UV light, invisible to the human eye). But there are also some exceptions. Now it is possible to have an upconversion luminescence where infrared is used as the excitation wavelength.

sigma-aldrich.com/upconversion



The principle behind fluorogenic media is simple - the target organisms are characterized by enzyme systems that metabolize the substrates to release the fluorogen. The fluorogen can then be visually detected by direct observation under UV lamp. Direct confirmation of the target organism without further testing is sometimes possible. Today, it is also possible to detect and differentiate more than one organism on the same plate in combination with chromogenic substrates and the help of selective media.



Figure 4: An example of a fluorogenic reaction. In the presence of a β -galactosidase positive organism, the 4-methylumbelliferyl galactoside (MUG) is split and results in the fluorophore and in free galactose.

Advantage of fluorogenic media:

- Faster results (compared to traditional methods)
- Reliable visual detection (often no further testing required)
- Additional testing possible directly from the media

In recent years, great strides have been made in the sector of fluorescence. There are dyes, labels and substrates which have a stronger fluorescence and are more stable. Additional advancements in the knowledge of enzyme and species specificity have also occurred within the past year. These recent gains in the development of selective agents and diverse fluorescence compounds have led to an impressive range of fluorescence tests.

| Cat. No. | Medium | Description | | | |
|----------|--|--|--|--|--|
| 51413 | Plate Count MUG Agar | For determination of plate count of microorganisms in milk and other dairy products and especially <i>E. coli</i> by fluorogenic method. | | | |
| 95273 | VRB MUG Agar | Selective medium for the detection and enumeration of coliform bacteria, in particular <i>E. coli</i> . Gram-positive accompanying flora are extensively inhibited by crystal violet and bile salts. A co change to red indicates lactose-positive colonies. | | | |
| 63014 | MacConkey MUG Agar | For the isolation of <i>Salmonella, Shigella</i> and coliform bacteria, in particular <i>E. coli</i> , from diverse material. Bile salts and crystal violet extensively inhibit the Gram-positive flora. The presence of lactose and neutral red indicate lactose-positive colonies from which <i>E. coli</i> can be identified by fluorescence in the UV. | | | |
| 78996 | HiFluoro™ Pseudomonas Agar Base | Used as a selective medium for the isolation of <i>Pseudomonas aeruginosa</i> from pus, sputum ar drains etc. <i>Pseudomonas aeruginosa</i> breaks the fluorogenic compound to release the fluorog which produces a visible fluorescence under long wave UV light. | | | |
| 62634 | LST-MUG Broth | Fluorescent method for the detection of <i>E. coli</i> . | | | |
| 16016 | BRILA MUG Broth* | Bile and brilliant green extensively inhibit the growth of accompanying flora, in particular Gram- positive microorganisms. The presence of <i>E. coli</i> results in fluorescence in the UV. | | | |
| 44782 | E. coli O157:H7 MUG Agar* | Selective agar for the isolation and differentiation of enterohemorrhagic (EHEC) <i>E. coli</i> O157:H7- strains from food and clinical material. | | | |
| 51489 | HiCrome [™] Rapid Coliform Broth Used for detection and conformation of <i>E. coli</i> and coliforms on the basis of enzyme reaction from water samples, using a combination of chromogenic and fluorogenic | | | | |
| 09142 | HiCrome ECD Agar with MUG* For the detection of <i>E. coli</i> in water and food samples by using a combination of chromo fluorogenic substrate. | | | | |
| 44657 | ECD MUG Agar* | The bile-salt mixture in this <i>E. coli</i> Direct Agar extensively inhibits the non-obligatory intestinal accompanying flora. Fluorescence in the UV demonstrates the presence of <i>E. coli</i> . | | | |
| M1678 | MUG EC Broth | Used for the detection of <i>E. coli</i> by a fluorogenic procedure. | | | |
| 17165 | MUG Tryptone Soya Agar | A Agar For cultivation of fastidious and nonfastidious microorganisms by fluorogenic method. | | | |

Table 3: Fluorogenic media (*not available in U.S.A.)



Analytical



Vitroids[™] – The Quantitative and Qualitative Control Strains

Jvo Siegrist, Product Manager Microbiology — ivo.siegrist@sial.com

Control of test results, or general testing performance, is currently an important issue in microbiology.

Organizations like ISO and UKAS either have regulations or recommendations to regularly verify the accuracy of various test methods. One really important requirement is to use certified microorganism standards, or if not available, at least a reference strain from a recognized culture collection like ATCC, DSMZ and NCTC. The test strains should show all relevant properties used for the applied application.

According to ISO 11133-1, it is possible to culture reference strains one time to produce reference stock cultures which are then controlled for purity and for use in biochemical tests. They should be stored in a freezer or in a freeze-dried form in small aliquots; however, defrosted cultures should not be refrozen. It is preferable for working strain cultures to be made out of stock cultures, and they should not be subcultured again.

In UKAS LAB 31 it is written that to conform with ISO/IEC 17025:2005, it is necessary to use, where possible, a traceable control strain which is certified according to ISO 17025. The strains should be from a recognized national culture collection or from a reference material producer accredited to ISO Guide 34.

For detailed performance testing, organizations such as ISO recommend special protocols, an example being media testing under ISO/TS 11133-2.

Reference Strain

- From an official strain collection (e.g. ATCC, NCTC, DSMZ)
- Defined, described and cataloged, at least down to the genus and species

Stock Culture

Prepared from a reference strain

Working Culture

Produced from a stock culture and used for control

Figure 5: Definition of ISO for control strains

Microorganisms in Discs

Sigma-Aldrich's innovation, small discs called Vitroids, is a solution which supports all of these requirements. The patented technology allows the production of stable and highly reliable certified microorganisms. The discs contain a defined certified number of colony forming units (CFU) from defined bacteria, molds and yeasts from a recognized culture collection like ATCC and NCTC. The discs are easy-to-use since they can be placed directly in water, diluent, broth or even on agar plates.



Figure 6: Classical plate count

The Vitroids contain highly viable bacteria, and when placed in contact with media, they dissolve rapidly and start to grow without a lag-phase. The viability of the CFU in a disc is stable for at least one year (for most organisms, more than two years) when kept at -20 °C. It is not a problem if the product is transported at ambient temperature or if it is kept for a short time in the refrigerator. The discs are produced under the ISO guideline 34 and are certified according ISO 17025. Each disc is packed in an individual tube with some desiccant and the tubes are then packed in mylar foil. Each package comes with a comprehensive certificate of analysis reporting the CFU and standard deviation.

All of the above mentioned features help microbiologist to have reliable results, save a lot of time (labor, documentation) and lower costs.

Reference microorganisms in certified and defined colony forming units (CFU)

- Standards in concentrations of 30–50,000 CFU per disc
- Produced acc. ISO Guide 34
- Certified acc. ISO 17025
- Delivered with detailed certificate of analysis
- Reference strains from ATCC, NCTC, etc.
- Minimum one year shelf life at -20 °C (usually two years)
- No lag-phase
- Amazingly low standard deviation (e.g. 100 CFU+/- 4%)

The Vitroids are packed in microvials and are lens shaped. They are colored and can easily be seen on the top of the filter plug. Below the filter plug is a desiccant which should be yellow. Clear (opaque) desiccant indicates the silica gel capacity has been depleted and is no longer effective. The tubes are packed in a resealable aluminium foil bag.



Figure 7: Vitroids in the vial

Preparation

Rehydrate the disc with a common phosphate buffer, or place the disc onto a solid or into a liquid medium. The rehydration process takes approximately 10 minutes. On solid media, the disc forms a droplet that can be spread with a loop. Liquid media may simply be shaken to dissolve the disc. The discs can be rehydrated in as little as 100 μ L of water (preferably a buffer) or added into larger volumes, e.g. 100 mL, for general water testing methods (MF, MTF, Quanti-Tray, etc). It is also possible to add the disc to the media for pour plate techniques.



Figure 8: A Vitroid just put into buffer (left) and after 10 minutes it is completely dissolved (right).



Figure 9: A single Vitroids disc on an agar plate. After about ten minutes on a plate, it is rehydrated automatically and forms a droplet (no water addition is needed). The drop can be spread with a loop.

| viciolus lest scialits | Ongin | Strain # | CFU | Cat. NO. |
|---|-------|----------|--------|----------|
| Aspergillus brasiliensis | ATCC | 16404™ | 80 | RQC15003 |
| Bacillus subtilis | ATCC | 6633™ | 10,000 | RQC02258 |
| Bacillus subtilis | ATCC | 6633 | 80 | RQC16003 |
| Candida albicans | ATCC | 10231™ | 80 | RQC14003 |
| Citrobacter freundii | ATCC | 8090 | 200 | RQC02105 |
| Clostridium perfringens | NCTC | 10240 | 30 | RQC02351 |
| Clostridium perfringens | NCTC | 10240 | 200 | RQC02355 |
| Clostridium perfringens | NCTC | 10240 | 500 | RQC20106 |
| Clostridium sporogenes | ATCC | 19404™ | 80 | RQC19003 |
| Enterobacter aerogenes | ATCC | 13048™ | 50 | RQC01652 |
| Enterobacter aerogenes | ATCC | 13048 | 200 | RQC01655 |
| Enterobacter aerogenes | ATCC | 13048 | 1,000 | RQC01657 |
| Enterobacter aerogenes | ATCC | 13048 | 10,000 | RQC01658 |
| Enterococcus cloacae | ATCC | 35030™ | 50 | RQC21102 |
| Enterococcus faecalis | ATCC | 19433™ | 50 | RQC01772 |
| Enterococcus faecalis | ATCC | 19433 | 200 | RQC01774 |
| Enterococcus faecalis | ATCC | 19433 | 500 | RQC01775 |
| Enterococcus faecalis | ATCC | 19433 | 1,000 | RQC01777 |
| Escherichia coli | ATCC | 11775™ | 50 | RQC01702 |
| Escherichia coli | ATCC | 11775 | 200 | RQC01705 |
| Escherichia coli | ATCC | 11775 | 1,000 | RQC01707 |
| Escherichia coli | ATCC | 11775 | 10,000 | RQC01708 |
| Escherichia coli | ATCC | 8739™ | 80 | RQC11003 |
| Heterotrophic Organisms | | | 100 | RQC02504 |
| Legionella bozemanii | NCTC | 11368 | 50,000 | RQC02908 |
| <i>Legionella pneumophila</i> (serogroup 1) | NCTC | 12821 | 50,000 | RQC02008 |
| Listeria monocytogenes | ATCC | 19115™ | 30 | RQC01901 |
| Pseudomonas aeruginosa | ATCC | 9027™ | 30 | RQC02202 |
| Pseudomonas aeruginosa | ATCC | 9027 | 100 | RQC02204 |
| Pseudomonas aeruginosa | ATCC | 9027 | 50 | RQC12002 |
| Pseudomonas aeruginosa | ATCC | 9027 | 200 | RQC12005 |
| Pseudomonas aeruginosa | ATCC | 9027 | 1,000 | RQC12007 |
| <i>Salmonella enterica</i> subsp. Enterica serovar Abony | NCTC | 6017 | 80 | RQC18003 |
| Salmonella enterica subsp. Enterica serovar Typhimurium | ATCC | 14028 | 50 | RQC17002 |
| Salmonella goldcoast | NCTC | 13175 | 30 | RQC02301 |
| Staphylococcus aureus susp. Aureus | ATCC | 6538 | 50 | RQC13002 |
| Staphylococcus aureus susp. Aureus | ATCC | 6538 | 200 | RQC13005 |
| Vitroids Blank | | | 0 | RQC0001 |
| | | | | |

Table 4: Vitroids portfolio

ATCC Licensed Derivative

ATCC Licensed Derivative

Many of RTC's Vitroids microbiological standards are offered as Quality Control microbiology products under the American Type Culture Collection (ATCC[®]) Licensed

Derivative® program. Look for the ATCC Licensed Derivative emblem for products derived from ATCC cultures. The ATCC Licensed Derivative emblem signifies ATCC-derived products are endorsed by ATCC. Products displaying the emblem are the only ones for which the quality of the ATCC ingredient can be properly assured. ATCC does not endorse products including ATCC ingredients from companies that are not members of the LDP program.

The ATCC Licensed Derivative emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. RTC is licensed to use these trademarks and to sell products derived from ATCC cultures.





Sigma-Aldrich[®] Worldwide Offices

Argentina

Free Tel: 0810 888 7446 Tel: (+54) 11 4556 1472 Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097 Free Fax: 1800 800 096 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10 Fax: (+43) 1 605 81 20

Belgium

Tel: (+32) 3 899 13 01 Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425 Tel: (+55) 11 3732 3100 Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400 Free Fax: 1800 265 3858 Tel: (+1) 905 829 9500 Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395 Fax: (+56) 2 495 7396

People's Republic of China Free Tel: 800 819 3336

Tel: (+86) 21 6141 5566 Fax: (+86) 21 6141 5567

Czech Republic Tel: (+420) 246 003 200 Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00 Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250 Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408 Free Fax: 0800 031 052 Tel: (+33) 474 82 28 88 Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000 Free Fax: 0800 64 90 000 Tel: (+49) 89 6513 0 Fax: (+49) 89 6513 1169

Hungary Tel: (+36) 1 235 9055

Fax: (+36) 1 235 9068

India Telephone

Bangalore: (+91) 80 6621 9400 New Delhi: (+91) 11 4358 8000 Mumbai: (+91) 22 4087 2364 Pune: (+91) 20 4146 4700 Hyderabad: (+91) 40 3067 7450 Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550 New Delhi: (+91) 11 4358 8001 Mumbai: (+91) 22 2579 7589 Pune: (+91) 20 4146 4777 Hyderabad: (+91) 40 3067 7451 Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888 Free Fax: 1800 600 222 Tel: +353 (0) 402 20370 Fax: + 353 (0) 402 20375

Israel

Free Tel: 1 800 70 2222 Tel: (+972) 8 948 4222 Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018 Tel: (+39) 02 3341 7310 Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300 Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111 Free Fax: (+82) 80 023 8111 Tel: (+82) 31 329 9000 Fax: (+82) 31 329 9090

Luxembourg Tel: (+32) 3 899 1301 Fax: (+32) 3 899 1311

Malaysia

Tel: (+60) 3 5635 3321 Fax: (+60) 3 5635 4116

Mexico Free Tel: 01 800 007 5300 Free Fax: 01 800 712 9920 Tel: (+52) 722 276 1600 Fax: (+52) 722 276 1601

The Netherlands Tel: (+31) 78 620 5411 Fax: (+31) 78 620 5421

New Zealand Free Tel: 0800 936 666 Free Fax: 0800 937 777 Tel: (+61) 2 9841 0555 Fax: (+61) 2 9841 0500

Norway Tel: (+47) 23 17 60 00 Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00 Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180 Free Fax: 800 202 178 Tel: (+351) 21 924 2555 Fax: (+351) 21 924 2610

Russia Tel: (+7) 495 621 5828 Fax: (+7) 495 621 6037

Singapore Tel: (+65) 6779 1200

Fax: (+65) 6779 1822

Slovakia Tel: (+421) 255 571 562 Fax: (+421) 255 571 564

South Africa Free Tel: 0800 1100 75 Free Fax: 0800 1100 79 Tel: (+27) 11 979 1188 Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376 Free Fax: 900 102 028 Tel: (+34) 91 661 99 77 Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200 Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80 Free Fax: 0800 80 00 81 Tel: (+41) 81 755 2511 Fax: (+41) 81 756 5449

Thailand

Tel: (+66) 2 126 8141 Fax: (+66) 2 126 8080

United Kingdom Free Tel: 0800 717 181 Free Fax: 0800 378 785 Tel: (+44) 01747 833 000 Fax: (+44) 01747 833 574

United States Toll-Free: 800 325 3010

Toll-Free Fax: 800 325 5052 Tel: (+1) 314 771 5765 Fax: (+1) 314 771 5757

Vietnam Tel: (+84) 8 3516 2810 Fax: (+84) 8 6258 4238

Internet sigma-aldrich.com

Enabling Science to Improve the Quality of Life Order/Customer Service (800) 325-3010 • Fax (800) 325-5052 Technical Service (800) 325-5832 • sigma-aldrich.com/techservice Development/Custom Manufacturing Inquiries **SAFC**^{*} (800) 244-1173 Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters 3050 Spruce St. St. Louis, MO 63103 (314) 771-5765 sigma-aldrich.com

©2013 Sigma-Aldrich Co. LLC. All rights reserved. SAFC and SIGMA-ALDRICH are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. FLUKA is a registered trademark of Sigma-Aldrich GmbH. FluoroSELECT is a trademark of Sigma-Aldrich Production GmbH. Vitroids is a trademark of Sigma-Aldrich Co. LLC. ATCC and Licensed-Derivative are registered trademarks of American Type Culture Collection. 16404, 6633, 10231, 19404, 13048, 35030, 19433, 11775, 8739, 1911 5 and 9027 are trademarks of American Type Culture Collection. HiCrome and HiFluoro are trademarks of HiMedia Laboratories Pvt. Ltd. RapidCheck is a trademark of RT Corp. Fluka brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

PPT 11746 / T413024 1073