



Helping all people  
live healthy lives

Same quality  
Same quality **required!**  
required!

**BD BBL™ IC-XT PACK**

From the company that brought you Difco™ and RODAC™

# A Partner in your Microbiological Control


As a pharmaceutical manufacturer ourselves, we understand the strict requirements for manufacturing, quality control and compliance. Our goal is your goal: eliminate the risk of microbial contamination for the health and safety of your patients.

## Ideal Quality Experience

Quality is critical, which is why BD has implemented the following elements to ensure the best media products are available for your use:

- **Vertically Integrated Production Process**

BD maintains control of the production process from peptone digestion to the finished product to consistently provide you with high quality media.



- **Robust Quality System**

BD ensures lot-to-lot consistency of the highest quality media through our ISO 9000 quality systems and cGMP manufacturing practices.

- **Automated Change Notification System**

BD is committed to your awareness of changes to your media products, so that you can quickly assess any impact to your procedures.

- **Certificates of Suitability**


BD saves you valuable time and resources by providing information about our dehydrated culture media used in these sterile plates to support regulatory requests for BSE risk evaluation.



















# BD BBL™ IC-XT Pack Isolator Cleanroom, eXtended Temperature

The BD BBL™ IC-XT Pack sterile plated media for EM applications provide you with:

- **Storage Flexibility:** A storage temperature range of 2-25°C provides you the highest flexibility and resource utilization.
- **Resource Efficiency:** Maximize batch sizes and reduce incoming QC testing with up to 6 months shelf-life.
- **Sample Security:** Locking lid RODAC™ plates minimize risks during sampling and handling.
- **Optimal Media Appearance:** A robust sleeve material with desiccant reduces environmental effects on the media.
- **Sterility Assurance:** An SAL claim of  $10^{-5}$  minimizes the risk of false positives and protects your environments from the risk of contaminants.
- **Ease of Transfer:** Triple-wrapped packaging enables easy transfer and storage while maintaining sterility.
- **VHP Resistance:** Validated VHP resistance of the sleeve material and hanging holes on each bag improves handling in the isolator.
- **Ease of Opening:** Even with isolator gloves, opening and removal from the packaging is simple.
- **Unique Traceability and Media Differentiation:** new colored coded barcode labels according to GS1 standards help you to identify easily and track both your medium and sample type.





Catalog numbers	Description	Type	Color of the 2D Barcode Label	Units per Pack
257632	BD BBL™ IC-XT Pack Trypticase™ Soy Agar	90 mm		100
257635	BD BBL™ IC-XT Pack Trypticase™ Soy Agar with Lecithin and Polysorbate 80	90 mm		100
257633	BD BBL™ IC-XT Pack Trypticase™ Soy Agar with Lecithin and Polysorbate 80 and Penicillinase	90 mm		100
257616	BD BBL™ IC-XT Pack Sanitizer Neutralizing Agar	90 mm		100
257662	BD BBL™ IC-XT Pack Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	90 mm		100
257619	BD BBL™ IC-XT Pack Sabouraud Dextrose Agar with Lecithin and Polysorbate 80 and Chloramphenicol	90 mm		100
257614	BD BBL™ IC-XT Pack CDC Anaerobe Agar and Penicillinase	90 mm		100
257661	BD BBL™ IC-XT Pack R2A Agar	90 mm		100
257637	BD BBL™ IC-XT Pack Trypticase™ Soy Agar with Lecithin and Polysorbate 80	RODAC™ LL		100
257634	BD BBL™ IC-XT Pack Trypticase™ Soy Agar with Lecithin and Polysorbate 80 and Penicillinase	RODAC™ LL		100
257617	BD BBL™ IC-XT Pack Sanitizer Neutralizing Agar	RODAC™ LL		100
257630	BD BBL™ IC-XT Pack Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	RODAC™ LL		100
257628	BD BBL™ IC-XT Pack Sabouraud Dextrose Agar with Lecithin and Polysorbate 80 and Chloramphenicol	RODAC™ LL		100
257615	BD BBL™ IC-XT Pack CDC Anaerobe Agar and Penicillinase	RODAC™ LL		100
257636	BD BBL™ IC-XT Pack Trypticase™ Soy Agar with Lecithin and Polysorbate 80	150 mm		30
257629	BD BBL™ IC-XT Pack Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	150 mm		30

LL = Locking Lid Plate

BD has brought you Difco™, RODAC™ and now BD BBL™ IC-XT Pack ... what's next ?



For additional information, please contact your BD sales representative or visit [www.bd.com](http://www.bd.com)



**BD**  
Tullastr 8-12  
69126 Heidelberg  
Germany  
Tel. +49 6221 3050  
Fax +49 6221 305304  
[www.bd.com](http://www.bd.com)