

DILUMATTM

Automated sample dilution for today and tomorrow

Dilumat

PIONEERING DIAGNOSTICS

NB DILUMAT 9307691 006-GB-C FAB.indd 1 28/02/2018 1

DILUMATTM



Standardize the sample preparation and increase your productivity

SPEED UP THE PREPARATION

25g 1:10 dilution in 10 sec*
375g 1:4 dilution in under 30 sec**
375g 1:10 dilution in under 80 sec**

* High speed mode ** DUAL pump mode

IMPROVE QUALITY

- Automated rotating arm, fully open access to the bag reduces greatly risks of contamination
- High accuracy meets BAM and ISO 7218
 requirements
 - Bidirectional LIMS communication
 - Automatic traceability, including communication with RFID tags
 - Customizable label







EVOLVE WITH YOUR NEEDS Only one DILUMAT™ for samples from 3g to 375g 3 dilution speeds adapted to all sample types Connect up to 6 different diluents (2 built-in pumps + 4 e-pumps) for maximum flexibility Upgrade to XL configuration for large samples (375g) Compatible system with DILUTUBE, DILUBAG™ and MASTERCLAVE® Heads-up control unit stand saves bench space (option)

EXPERIENCE TOTAL EASE OF USE

- Smart software and color changing LEDs inform user about dilution status
- Plug-and-play! No tedious installation needed
- No operator intervention with Autostart mode
- Cleanable mobile touchscreen control unit
- Small footprint enables customized workspace organization
- All surfaces are easy to clean





BLUE LINE™

Boost the productivity of your microbiology laboratory!

DILUBAG™

• 3 and 5 liter ready-to-use broth diluent bags



DILUTUBE™

• Sterile ready-to-use tubing



SMASHERTM

- High performance for up to 400ml and samples up to 75g (1:4 dilution)
- Shortest time lowest noise

VIDAS UP®

- VIDAS® UP provides next-day pathogens results
- Minimal hands-on-time and simple protocols

DILUMATTM



START EXPERT EVO w/ XL pack











DILUMAT START and DILUMAT EXPERT evo specifications depend on the model

MASTERCLAVE®

 Automated preparators for broth and agar culture media 9, 28 or 60 liter capacity





SMASHER™ XL

- Large capacity blender for two 2L or one 4L bags.
- Handles samples up to 375g (1:10 dilution)

TEMPOTM

- Automated & innovative enumeration of quality indicators
- Next-day results

9307691 006/GB/C bioMérieux - Document and/or pictures not legally binding. Modifications by bioMérieux can be made without prior notice • BIOMERIEUX, the blue logo, BLUE LINE, DILUBAG, DILUBAG, DILUBAG, DILUBAG, point and common and volve registered trademarks belonging to bioMérieux, or one of its subsidiaries or one of its subsidiaries or one of its companies. By ACS Lyon / Printed in France by Masson Communication 4388319127 RCS Saint-Malo. Photographs by AES CHEMUNEX

bioMérieux S.A. • 69280 Marcy l'Étoile • France • Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90 www.biomerieux.com • www.biomerieux-industry.com