



100 mL

Colilert*
Pseudalert*
Enterolert*
Enterolert*-DW

Presence/Absence Results Quick Reference Guide†



Colilert



Pseudalert



Enterolert Enterolert-DW



Coliforms



E. coli



-



+



P. aeruginosa



-



+

Enterococci



-



+

Enterococci

Enterolert

Enterolert-DW

 **250** mL

Colilert* 250
Pseudalert* 250
Enterolert* 250

Presence/Absence Results Quick Reference Guide†



†This guide is for reference only and is not a substitution for the product insert.
Please refer to the product insert for the complete procedure.

©2012 IDEXX Laboratories, Inc. All rights reserved. • 09-71598-00
*Colilert, Enterolert, Pseudalert are trademarks or registered trademarks of
IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.
The IDEXX Privacy Policy is available at idexx.com.

 **Colilert 250**  36°C ± 2.0°C  24



 **Pseudalert 250**  38°C ± 0.5°C  24



 **Enterolert 250**  41°C ± 1.0°C  24

