

# ECUANTE COSMETICS

ECUANTE  
Proficiency Tests

Below is detailed information regarding **Round I** (Circuit CO-M-01):

| Submission date | Analysis date | Round closing date | Report issue date |
|-----------------|---------------|--------------------|-------------------|
| 20/04/2026      | 21-24/04/2026 | 01/05/2026         | 01/06/2026        |

The **parameters** to be analysed will be as follows:

- Count of:
  - Total viable microorganisms
  - Moulds and yeasts
- Detection of:
  - Escherichia coli*
  - Candida albicans*
  - Burkholderia cepacia*
  - Estafilococos coagulasa* +
  - Pseudomonas aeruginosa*

The sample will consist of:

- Cosmetic matrix
- Vial containing the stabilised microbial suspension



For detailed information on pricing, please visit the website: <https://ecuate.labnovasl.com/ListaRondas/Index/>



Once the sample has been received, the following **procedure** must be followed:

①

Add the suspension of microorganisms from the vial to the cosmetic matrix provided.



②

Add the inactivating diluent, mix thoroughly and leave to act.



③

Homogenise the sample and carry out the analysis in accordance with each laboratory's usual methodology.



### Possible cosmetic matrices:



LOTIONS AND MILKS



BATH GELS



SHAMPOOS



CREAMS

## HOW TO REGISTER?

Go to the website (<https://labnovasl.com/programa-de-ensayos-de-aptitud-ecuanter/>) and click on "*Registration and enrolment for aptitude tests*".



- Labnova Distribuciones Agroalimentarias, S.L.
- C/ Vitoria, 274 - nave 138, Complejo Naves Taglosa, Pol. Ind. Gamonal - Villimar, 09007, BURGOS
- Phone: +34 947 04 06 63
- Email: [info@ecuanter.com](mailto:info@ecuanter.com)

