

GMO Extraction Kit

Ref. 4466336, 50 isolations

GMO Extraction Kit provides a fast and easy silica-based DNA purification method to obtain DNA from a big variety of food samples, like raw material (seeds, flours,...), processed food, and feed.

instituto de medicina genómica
imegen



Specific

Adapted to food samples for transgenic material analysis

Safe

No toxic reagents

Reliable


Inhibitors removal;
Isolation of high pure DNA
($A_{260}/A_{280} > 1.8$)

Easy to use

Minimal handling steps

Fast

Total assay time: 2 hours



Related Products

TaqMan® GMO Screening Kit, Ref. 4466334

Detection of the main transgenic events described in databases TaqMan® Roundup Ready Soya Quantification Kit,

Ref. 4466335

Quantification of Roundup Ready Soya vs. total soya present in one sample

To learn more about our food testing solution go to www.lifetechnologies.com

Email: food@lifetech.com to try the GMO Extraction kit today