

## TaqMan® Roundup Ready Soya Quantification Kit

Ref. 4466335, 2 x 48 reactions

**TaqMan® Roundup Ready Soya Quantification Kit** allows determining, by Real-Time PCR, the percentage of *Soya Roundup Ready* (Soya RR) with respect total soya in any food sample or feed.

The Kit is optimized for its use in all Applied Biosystems Real-Time PCR Instruments.



### Efficient


Use of **TaqMan® MGB probes**; detection of small DNA fragments, highest PCR efficiency

### Flexible

Adapted to processed DNAs

### Sensible

Quantification limit  
20 DNA copies  
Detection limit  
3 DNA copies




### Related Products

**TaqMan® GMO Screening Kit,**  
Ref. 4466334

Detection of the main transgenic events described in databases

**GMO Extraction Kit,** Ref. 4466336

Isolation of high pure DNA from any kind of food or feed



**Instituto de Medicina Genómica, S.L.**

Parque Científico-Universidad de Valencia  
Catedrático Agustín Escardino, 9  
46980 Paterna (Valencia)

[info@imegen.es](mailto:info@imegen.es); [www.imegen.es](http://www.imegen.es)

Tel: 963 212 340, Fax: 963 212 341